Features of Inspire

Able to select appropriate disease detection mode
User can take a single image or multiple images at once depended upon the mode they choose.

Image quality

Dental handpiece design Most familiar design for dentists.

Metal body case Metal body case increases product durability.

Slim head design makes it easier to approach the oral cavity with camera.

Retraction head design
Retract patient's tongue and cheek to get better image.

Dual capture buttons Capture buttons are placed on the top and bottom side of camera to make for an easy capture in any position

LED indicators LED indicators show the device status and mode.

Mouse function

Built-in mouse function for user convenience.

Zoom function
Ability to enlarge the image with zoom feature to support user convenience. Support twain driver
Support twain driver allows for use in almost any dental imaging software.

Software compatibility
Able to be used with almost any dental imaging software in the market.

Detachable USB cable Detachable USB cable makes it easy to move the camera to other dental chairs.

Dental Camera Curing Light Apex Locator -IC-WHT60 -IC-WHW61 -IC-DCAM70 -IC-DCAU71 -IC-EZAM80 -IC-EZAU81 -IC-WHMA90 -IC-WHCD100

Zoom function

Ability to enlarge the image with zoom feature to support user convenience.





* Need to install "Zoom function" program on the PC.

Software compatibilityAble to be used with almost any dental imaging software in the market.



1 Specification

| Model | IC-WHCD100 | |
|-----------|--------------|----------------------------|
| Camera | Power source | 5V (USB power) |
| | Size | Approx. 129 x 21.5 x 32 mn |
| | Weight | Approx. 61g (2.2oz) |
| | Sensor type | 5M CMOS sensor |
| | Lighting | LED |
| | Focusing | Macro to infinity |
| | Focus range | |
| USB cable | cable length | Approx. 3 M (9.8ft) |
| | | |

2 Recommended PC system

| | PC system: | 3.3GHz CPU, 4GB RAM |
|--|------------|---------------------|
| | OS: | Windows 7 |
| | USB port: | USB 2,0 |
| | | |



GoodDrs



Inspire

Autofocus Dental Camera with Mouse



E-mail. info@gooddrs.com



Good Doctors Co., Ltd. #208, B-dong, 283 Bupyeong-daero, Bupyeong-gu (Woolim Lions Valley, Cheongchoen-dong), Incheon, 21315, Republic of Korea Tel: +82 32 424 6325 Fax: +82 32 424 6326

(€



Patent pending



Camera is able to focus from macro mode (close) to infinity that includes a full arch and face shot.









FHD image support

The device shows FHD (1080p) in live image supported by a 5M sensor.





Original image



Full HD

Adjust the color impression on camera to fit different monitors or screen.









Same color chart image with xrite





- This camera mode provides the ability to examine the inside of enamel by use of a high bio-transparency light source.
- Useful for detecting proximal caries and micro cracks,
- Advanced technology that allows for an image showing enamel penetration without x-ray radiation.











- Disease detection technology using natural fluorescence that occurs under specific illumination.
- The FL mode provides simple ways to detect disease: User will see green fluorescence on normal enamel and red fluorescence where there is bio-film such as plaque or calculus,
- Florescence loss where there is early caries (white spot).
- The FL mode provides the ability to detect early caries, tooth crack, and secondary caries.
- The FL mode is able to detect residual composite resin, · The FL Mode supports 2 mode types: F mode and S mode,









Limits on each detection mode

Each disease detection method has advantages as well as has limits.



Hard to detect early stage caries,

TR mode

No quantitative measurement is produced

Subjective interpretation of whether a dark area on the gray-scale images is caries or a dental abnomaly leads to many false positives

FL mode Fluorescence mode Bacterial porphyrins, stains, tartar, food debris, and prophylaxis paste all fluoresce under the wavelengths used in these device, whether or not caries is presents, they can lead to false positive readings and unnecessary treatment. In addition,

Streptococcus mutans and lactobacilli, the key bacterial initiators of caries, do not have porphyrins that fluoresce when exposed to the light emitted by the device,

Inspire takes advantages of each disease detection mode (MI, TR, FL), and to overcome limits of each detection mode it provide 3 in1 shooting for comparative analysis.



3in1 3in1 shooting example.

