



[ptfesutures.com](http://ptfesutures.com)

# Golnit® PTFE Suture

Periodontal, Bone Graft  
and Implant Surgery

The ideal suture for  
dental surgery



## KEY ADVANTAGES

- Completely non-absorbable, 100% pure PTFE
- Monofilament – prevents bacterial wicking
- Biologically inert and chemically nonreactive
- Decreased inflammation and faster wound repair
- Soft – comfortable for patients
- Great handling characteristics and knot security

## APPLICATIONS

- Dental bone grafting
- Soft tissue grafting
- Implant surgery
- Periodontal surgery
- Ridge augmentation
- Guided tissue regeneration (GTR)

## WHY GOLNIT® PTFE SUTURES?

- The largest variety of PTFE sutures ever offered in the market
- Custom-made sutures for dental and aesthetic surgeries
- Needles made of highest grade, 302 series stainless steel wire
- Excellent needle penetration and smooth surface

We carry many more specifications for USP 6/0, 5/0, 4/0, 3/0, 2/0. Please contact us for additional sizes and lengths.

| NEEDLE DESCRIPTION             | NEEDLE LENGTH | AVAILABLE USP SIZES |
|--------------------------------|---------------|---------------------|
| Reverse Cutting 3/8 Circle ▼   | 24mm          | 2/0; 3/0            |
| Reverse Cutting 3/8 Circle ▼   | 21mm          | 2/0; 3/0; 4/0       |
| Reverse Cutting 3/8 Circle ▼   | 20mm          | 2/0; 3/0; 4/0       |
| Reverse Cutting 3/8 Circle ▼   | 19mm          | 2/0; 3/0; 4/0; 5/0  |
| Reverse Cutting 3/8 Circle ▼   | 18mm          | 2/0; 3/0; 4/0       |
| Reverse Cutting 3/8 Circle ▼   | 16mm          | 3/0; 4/0; 5/0; 6/0  |
| Reverse Cutting 3/8 Circle ▼   | 13mm          | 4/0; 5/0; 6/0       |
| Taper Point PERIO 1/2 Circle ☐ | 13mm          | 4/0; 5/0; 6/0       |
| Reverse Cutting 3/8 Circle ▼   | 12mm          | 4/0; 5/0; 6/0       |
| Reverse Cutting 3/8 Circle ▼   | 9mm           | 5/0; 6/0            |
| Taper Point 3/8 Circle ☐       | 9mm           | 5/0; 6/0            |