


# GeoMEDI

Implant Solutions  
PREMIUM CAD/CAM PRODUCT



Implant Library System  
(3Shape, exocad, Dentalwings)

LL/LS [Multibase]

CMFit [Premill, Blank]

3D-Analog

SC / GMS / GTS - Scanbody

Sirona Multibase

CCM Casting Abutment

Geo Screw Driver

# We are exploring excellent dental implant solutions

Since 2009, we have been developing patient-customize abutment, and all of our dentists and dental technicians are also trying to make the fittest implant prosthesis for each individual person.

## For Ideal Implant Restorations

We have a commitment to provide the best implant solution at a reasonable price, easy and fast way for everyone.

Co-Cr Abutment Wax up & Casting

CAD/CAM Abutment CAD design & CAM milling

## For Digital Implant Solution

It is possible to make implant prosthesis by using a scanner, dental CAD, milling machine and 3D printer as developing our own digital solution for dental clinic and laboratory.

# Contents

■ Company History / Certification -----	04
■ Dental CAD/CAM System Workflow -----	05
■ Design concept -----	06
GMS Scanbody /SC Scanbody -----	06
TiN screw / Geo LS Multibase -----	08
Geo LL Multibase -----	09
Geo CMFit / Reverse Jig -----	10
Geo Multi-Holder for CMFit-----	11
3D-Analog-----	12
Geo cad/cam prosthetic system-----	13
with digital library	
■ Compatible to Implant Systems-----	14
Osstem Implant TS -----	14
Osstem Implant US -----	15
Osstem Implant SS -----	16
Dentium Super-Line / Implantium -----	17
Dentium NR-Line -----	18
Megagen AnyRidge -----	19
Megagen AnyOne -----	20
Megagen EZ-plus -----	21
Megagen Ex-Feel [Internal] -----	22
Megagen MiNi -----	22
DIO SM -----	23
DIO UF / UFII -----	23
NeoBiotech IS -----	24
Dentis S-Clean / OneQ -----	25
Straumann Bone Level -----	26
Straumann Tissue Level -----	27
Straumann BLX -----	27
Nobel Biocare Active -----	28
Nobel Biocare Branemark -----	29
Nobel Biocare Replace -----	30
Dentsply Xive / Frialit2 -----	31
Dentsply Ankylos -----	32
Biomet 3i Certain -----	33
Biomet 3i Osseotite -----	34
Implant Direct 3.0 -----	35
Zimmer Tapered Screw-vent -----	36
Zimmer SPLINE -----	37
Astra Tech Osseospeed TX -----	38
Astra Tech Osseospeed EV -----	39
CAMLOG -----	40
Conelog -----	40
MIS SEVEN -----	41
MIS C1 -----	42
SIC ace / max -----	43
Thommen SPI -----	44
BEGO -----	45
MYTIS Arrow Implant -----	46
Kyocera POI -----	46
Biohorizons Laser-Lok 3.0 -----	47
Anthogyr AXIOM -----	48
WARANTEC Oneplant -----	49
Tekka In-Kone -----	50

Keystone PrimaConnex / Max-TL -----	50
Neudent CM -----	51
Neudent GM -----	51
BICON Implant -----	52
Medentika IPS -----	53
GC Aadvia Implant -----	53
BIOTECH KONTACT -----	54
Medentis medical ICX -----	55
B&B Dental -----	56
Bredent Medical Sky -----	56
NEOSS ProActive -----	57
Champions (R)evolution -----	57

■ MultiUnit(Multi-base) -----	57
■ GTS Transfer Scanbody -----	60
■ Sirona Multibase -----	62
■ CCM Casting Abutment -----	66
■ Geo Screw Driver -----	71
■ Instructions of CAD/CAM Abutment -----	74
■ Instructions of CCM Casting Abutment -----	75
■ Geo Library -----	77
■ Selection of Library -----	78
■ Geo Bridge-Single Library -----	79
■ Products List -----	81




# Company History & Certification



- 2009 • 01 Released Geo CCM Abutment
- 2010 • 03 Established GeoMedi. Co., Ltd.
- 2011 • 12 Released Geo Link, CMFit, Library, Geo Driver
- 2012 • 10 The 1st Geo Symposium at Seoul National University
- 2013 • 08 Obtained KGMP  
08 Obtained KFDA - implant, endosseous, laboratory supplies  
10 Obtained KFDA - Geo Link Abutment, Geo CCM Abutment  
12 Obtained KFDA - Geo Cmfit Abutment
- 2014 • 07 Obtained KFDA - Geo Tcp  
09 The 2nd Geo Symposium at Seoul Trade Exhibition & Convention
- 2015 • 02 Obtained ISO13485, International Organization for Standardization' Medical devices Certification  
02 Obtained CE, the European quality certification  
04 Released Geo CAD/CAM system, Geo Mill  
06 Certified as a venture business  
07 Expansion & Relocation of corporation  
09 Established research institute affiliated with a corporation
- 2016 • 03 Obtained KFDA - Geo Scanbody, Geo 3D Analog, Geo screw  
06 Obtained KFDA - Geo Multibase Abutment
- 2017 • 07 Obtained CE, the European quality certification - Geo CCM Abutment, Geo CMFit Abutment  
07 Obtained EN ISO13485 : 2012  
11 Obtained certification of INNOBIZ - Ministry of SMEs and Startups
- 2018 • 05 Obtained CE, the European quality certification - Geo Multibase Abutment  
11 Obtained EN ISO13485 : 2016 (31720001)  
12 Achieved 1 Million dollars in export
- 2019 • 09 Obtained certification of MainBIZ - Ministry of SMEs and Startups
- 2020 • 02 Obtained KFDA - GTS Scanbody, Multi Screw Abutment
- 2021 • 04 Being selected as a Small-Ginat enterprises by Ministry of Employment and Labor  
05 Obtained KFDA - OMEGA System  
05 Obtained CE - GTS Scanbody & Geo Multiscrew Abutment  
12 Certified as a venture business - Innovative Growth

# Dental CAD/CAM System Workflow



# Design concept

## GMS Scanbody (Intra-oral)

- Biocompatible Titanium (Ti-6Al-4V ELI) from CARPENTER Co, guarantees durability which can be autoclavable and reusable.
- special coating allows to scan without any scan spray.
- provides with two different sizes. – Short(7.5mm) / Long(11.5mm)
- consistent performance and visible in x-rays.
- Various implant system available.
- Highly precise Geo library is provided for CAD s/w, (3shape, exocad, Dental Wings)
- Customized Abutment Design with CNC or Premilled abutment is available.
- Zirconia prothesis design with Geo Link.



GMS Short (7.5mm)

GMS Long (11.5mm)




# Design concept

## GMS Scanbody (Intra-oral)



GMS Scanbody / SC Scanbody

## SC Scanbody(Model)




Dental Lab Scan Body.  
Scan spray not necessary.  
Semipermanent scanbody which assures accurate scan data.

# Design concept

## TiN Screw


Surface treatment of TiN-coating  
- increased of abrasion resistance  
- improved fixation of screw



Screw made of Titanium Grade ELI  
Two Screw - TiN coating screw & Titanium screw  
TiN coating screw - Final Setting

## LS Multibase

Screw channel path can be changed for different angle  
(ball type angled driver necessary)



Geo Multibase & Screw made of Titanium Grade ELI  
Two Screw - TiN coating screw & Titanium screw



# Design concept

## LL Multibase


LL Multibase height is editable on Geo Library.



Tin-coating - High esthetic.  
LL Multibase for Posterior.

NEW Item

Gingival height - 2.0mm  
height - 6.0mm



Grooves - Good adhesive strength.  
No need for sandblasting.

Ideal design for making emergence profile.



Gingival height-1.3mm --> height-5.5mm

Gingival height-2.0mm --> height-6.0mm

Two different gingival height is available.

Zirconia - Screw retained restorations  
Customized Zirconia abutment  
Geo Multibase & Screw made of Titanium Grade ELI  
Two Screw - TiN coating screw & Titanium screw



Geo LL Multibase(Link)

# Design concept

## Geo CMFit [ Premill, Blank ]



Diameter of  $\varnothing 10$  &  $\varnothing 14$  size available

Laser marked Geo & Implant System code

Provide with two different screw (TiN coating & Titanium)

Qualified material - Ti-6Al-4V ELI made by CARPENTER co.

CMFit [Premill, Blank] & Reverse Jig

## Reverse Jig



Occlusal milling is possible depend on machine flexibility


Support pin is not necessary

Fix pre-milled abutment with ordinary screw

Addtional Geo holder(jig) is necessary


# Design concept

## Geo Multi-Holder for CMFit



Multi-holder and Reverse holder

Geo Multi-Holder for CMFit




Customized Geo Multi holder compatible with leading Milling Machines.  
Available for all type of implant system.



# Design concept

## Geo 3D Analog

3D Analog



Pioneer of providing 3D digital Analog

Can be used both stone model & 3D printed model

Reasonable price with excellent durability

Provided all type of implant system

Stainless Steel



Cross section  
3D Analog / LL -Link / Zirconia



# Design concept

Geo cad/cam prosthetic system with digital library

Dental CAD s/w

3shape ▶

exocad

dental wings


Dental CAM s/w

FOLLOW-ME!  
TECHNOLOGY GROUP  
hyperDENT

sum3D  
MILL BOX  
One step beyond

vhf

Geo Scanbody




Scan body that provide the precise library matching according to the user's environment.

CMFit Abutment




Geamedi Library provides prosthesis definite sitting of a custom abutment





Dedicated libraries for the various gingival height of Geo Ti-base abutment is well known for their emphasis on aesthetic and functional design.

3D Analog



Retention groove on 3D analog and extra space to fill up wax provides stable fix

Multifunctional analog for ordinary stone model and 3D printed model.



Precise direct interface milling from Geo geometry in cooperation with dental CAM s/w companies.

Geo ca/cam prosthetic system with digital library

# Compatible to Implant Systems

## Osstem Implant TS

### Osstem Implant TS

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm)  (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.2Hex
				20Ncm Ø3.0 / 3.5
				30Ncm Ø4.0 / 4.5 / 5.0

Description-Library code	0101 - SUM	0102 - SURO
Product Connection Platform Ø	TS Mini  Ø3.0 / 3.5	TS Standard  Ø4.0/4.5/5.0
SC-Scanbody	SC-SUM	SC-SURO
GMS-Scanbody(Short)	GMS-SUMS	GMS-SUROS
GMS-Scanbody(Long)	GMS-SUML	GMS-SUROL
LL-Multibase (Single)	GH : 1.3 LL2-SUM13-H2 GH : 2.0 LL2-SUM20-H2 GH : 3.0 -	LL2-SURO13-H2 LL2-SURO20-H2 LL2-SURO30-H2
LL-Multibase (Bridge)	GH : 1.3 LL2-SUM13-N2 GH : 2.0 LL2-SUM20-N2 GH : 3.0 -	LL2-SURO13-N2 LL2-SURO20-N2 LL2-SURO30-N2
LS-Multibase / Single	LS2-SUM-H2	LS2-SURO-H2
LS-Multibase / Bridge	LS2-SUM-N2	LS2-SURO-N2
CMFit Abutment(S) / Ø10	CMFit10-SUM-H2	CMFit10-SURO-H2
CMFit Abutment(B) / Ø10	CMFit10-SUM-N2	CMFit10-SURO-N2
CMFit Abutment(S) / Ø14	CMFit14-SUM-H2	CMFit14-SURO-H2
CMFit Abutment(B) / Ø14	-	CMFit14-SURO-N2
3D-Analog	3D-SUM	3D-SURO
Screw(5ea)	AC-SUM5	AC-SURO5

# Compatible to Implant Systems

## Osstem Implant US

\*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.2Hex 30Ncm

Description-Library code	0103 - EXM	0104 - EXR	0105 - EXW
Product \ Connection Platform Ø	US Mini	US Regular	US Wide
SC-Scanbody	SC-EXM	SC-EXR	SC-EXW
GMS-Scanbody(Short)	GMS-EXMS	GMS-EXRS	GMS-EXWS
GMS-Scanbody(Long)	GMS-EXML	GMS-EXRL	GMS-EXWL
LL-Multibase / Single	LL2-EXM-H2	LL2-EXR-H2	LL2-EXW-H2
LL-Multibase / Bridge	LL2-EXM-N2	LL2-EXR-N2	LL2-EXW-N2
LS-Multibase / Single	LS2-EXM-H2	LS2-EXR-H2	LS2-EXW-H2
LS-Multibase / Bridge	-	-	-
CMFit Abutment(S) / Ø10	CMFit10-EXM-H2	CMFit10-EXR-H2	CMFit10-EXW-H2
CMFit Abutment(B) / Ø10	-	CMFit10-EXR-N2	-
CMFit Abutment(S) / Ø14	CMFit14-EXM-H2	CMFit14-EXR-H2	CMFit14-EXW-H2
3D-Analog	3D-EXM	3D-EXR	3D-EXW
Screw(5ea)	AC-EXM5	AC-EXR5	AC-EXW5

# Compatible to Implant Systems

Osstem Implant SS

## Osstem Implant SS

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(10mm)  (7.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.2Hex 30Nm

Description-Library code	0106 - ISR	0107 - ISW
Connection Platform Ø	SS Regular	SS Wide  Ø6.0
Product	Ø4.8	
SC-Scanbody	SC-ISRW	SC-ISRW
GMS-Scanbody(Short)	GMS-ISRS	GMS-ISWS
GMS-Scanbody(Long)	GMS-ISRL	GMS-ISWL
LL-Multibase / Single	LL2-ISR-02	LL2-ISW-02
LL-Multibase / Bridge	LL2-ISR-N2	LL2-ISW-N2
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	CMFit10-ISR-02	CMFit10-ISW-02
CMFit Abutment / Ø14	CMFit14-ISR-02	CMFit14-ISW-02
3D-Analog	3D-ISR	3D-ISW
Screw(5ea)	AC-EXR5	AC-EXR5

# Compatible to Implant Systems

## Dentium Super-Line / Implantium

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm)  (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) 1.27Hex
				25~30Nm

Description-Library code	0201 - DER	
Product \ Connection Platform Ø	Universal	Ø3.7/4.0/4.5
SC-Scanbody	SC-DER	
GMS-Scanbody(Short)	GMS-DERS	
GMS-Scanbody(Long)	GMS-DERL	
LL-Multibase (Single)	GH : 1.3	LL2-DER13-H2
	GH : 2.0	LL2-DER20-H2
	GH : 3.0	LL2-DER30-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-DER13-N2
	GH : 2.0	LL2-DER20-N2
	GH : 3.0	LL2-DER30-N2
LS-Multibase / Single	LS2-DER-H2	
LS-Multibase / Bridge	LS2-DER-N2	
CMFit Abutment(S) / Ø10	CMFit10-DER-H2	
CMFit Abutment(B) / Ø10	CMFit10-DER-N2	
CMFit Abutment(S) / Ø14	CMFit14-DER-H2	
CMFit Abutment(B) / Ø14	CMFit14-DER-N2	
3D-Analog	3D-DER	
Screw(5ea)	AC-DER5	

# Compatible to Implant Systems

## Dentium NR-Line

\*set (Scanbody/Multibase/CMFit Abutment + screw)

OS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase	
	-	-	-	
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-	-		1.0square 20Ncm

Description-Library code	0202 - DNR
Product	Connection Platform Ø Standard
OS-Scanbody	OS-DNR
LL-Multibase / Single	-
LL-Multibase / Bridge	-
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-DNR-0
CMFit Abutment / Ø14	-
3D-Analog	-
Screw(Sea)	AC-DNR5

# Compatible to Implant Systems

## Megagen AnyRidge

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm)  (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.2Hex
				35Ncm

Description-Library code	0301 - MARR	
Product	Connection Platform Ø	Universal
		Ø3.5/4.0/4.5
SC-Scanbody	SC-MARR	
GMS-Scanbody(Short)	GMS-MARRS	
GMS-Scanbody(Long)	GMS-MARRL	
LL-Multibase (Single)	GH : 1.3	LL2-MARR-H2
	GH : 2.0	LL2-MARR20-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-MARR-N2
	GH : 2.0	LL2-MARR20-N2
LS-Multibase / Single	LS2-MARR-H2	
CMFit Abutment / Ø10	CMFit10-MARR-H2	
CMFit Abutment / Ø14	CMFit14-MARR-H2	
3D-Analog	3D-MARR	
Screw(5ea)	AC-MARR5	

# Compatible to Implant Systems

Megagen AnyOne

## Megagen AnyOne

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral) *Using 1.27Hex Driver*	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm)  (7.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				25Ncm

Description-Library code	0302 - DERA	
Product	Connection Platform Ø	Universal  Ø 3.5/4.0/4.5/5.0
SC-Scanbody	SC-DER	
GMS-Scanbody(Short)	GMS-DERS	
GMS-Scanbody(Long)	GMS-DERL	
LL-Multibase (Single)	GH : 1.3	LL2-DERA13-H2
	GH : 2.0	LL2-DERA20-H2
	GH : 3.0	LL2-DERA30-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-DERA13-N2
	GH : 2.0	LL2-DERA20-N2
	GH : 3.0	LL2-DERA30-N2
LS-Multibase / Single	-	
LS-Multibase / Bridge	-	
CMFit Abutment(S) / Ø10	CMFit10-DERA-H2	
CMFit Abutment(B) / Ø10	CMFit10-DERA-N2	
CMFit Abutment(S) / Ø14	CMFit14-DERA-H2	
CMFit Abutment(B) / Ø14	CMFit14-DERA-N2	
3D-Analog	3D-DER	
Screw(5ea)	AC-MA05	

# Compatible to Implant Systems

## Megagen EZ-plus

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-	-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) 1.2Hex 35Nm
	-	-		

Description-Library code	0303 - MEPM	0304 - MEPR
Connection Platform Ø	Ø 3.3	Ø 4.0/4.5/5.0
Product		
SC-Scanbody	SC-MEPM	SC-MEPR
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Multibase / Single	-	-
LL-Multibase / Bridge	-	-
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	CMFit10-MEPM-E	CMFit10-MEPR-E
CMFit Abutment / Ø14	-	-
3D-Analog	-	-
Screw(5ea)	AC-MEPM5	AC-MEPR5

Megagen EZ - plus

# Compatible to Implant Systems

## Megagen ExFeel Internal

\*set (Scanbody/Multibase/CMFit Abutment + screw)

Megagen ExFeel Internal

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.2Hex
				30Nm

Description-Library code	0305 - ITRK
Connection Platform Ø	Regular  Ø4.8
Product	
SC-Scanbody	SC-ITRW
GMS-Scanbody(Short)	-
GMS-Scanbody(Long)	-
LL-Multibase / Single	LL2-ITRK-02
LL-Multibase / Bridge	LL2-ITRK-N2
LS-Multibase / Single	-
CMFit Abutment(S) / Ø10	CMFit10-ITRK-02
CMFit Abutment(B) / Ø10	CMFit10-ITRK-N2
CMFit Abutment(S) / Ø14	CMFit14-ITRK-02
3D-Analog	3D-ITR
Screw(5ea)	AC-EXR5

## Megagen MiNi

\*set (Scanbody/Multibase/CMFit Abutment + screw)

Megagen MiNi

SC-Scanbody(Model)	Single LL-Multibase	Ø10 CMFit Abutment	Screw	Driver (Nm)
	-			1.2Hex
				15Nm Mini 3.0/3.25

Description-Library code	0306 - MMN
Connection Platform Ø	MINI  Ø3.0/3.25
Product	
SC-Scanbody	SC-MMN
LL-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-MMN-H
Screw(5ea)	AC-MMNS

# Compatible to Implant Systems

## DIO SM

SC-Scanbody(Model)	$\varnothing 10$ CMFit Abutment	Single LS-Link	Screw	Driver (Ncm)
 		-		1.7 Torx 30Ncm

Description-Library code	0401 - DSMM	0402 - DSMR
Connection Platform Ø	Mini	Regular
Product	$\varnothing 3.8 / 4.1$	$\varnothing 4.5 / 5.0$
Scanbody	SC-DSMM	SC-DSMR
LS-Link / Single	-	-
CMFit Abutment $\varnothing 10$	CMFit10-DSMM-T	CMFit10-DSMR-T
Screw(5ea)	AC-DSMMS5	AC-DSMR5

## DIO UF / UFII

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Link	Bridge LL-Link	Single LS-Link
	(11.5mm)  (7.5mm)			-
$\varnothing 10$ CMFit Abutment	$\varnothing 14$ CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
			1.2Hex	20Ncm $\varnothing 3.0 / 3.3$ 30Ncm $\varnothing 3.8 / 4.0 / 4.5 / 5.0$

Description-Library code	0403 - DN30	0404 - SURD
Connection Platform Ø	UFII Narrow $\varnothing 3.0 / 3.3$	Standard $\varnothing 3.8 / 4.0 / 4.5 / 5.0$
Product		
SC-Scanbody	SC-MMN	SC-SURO
GMS-Scanbody(Short)	-	GMS-SUROS
GMS-Scanbody(Long)	-	GMS-SUROL
LL-Link / Single	GH : 1.3 GH : 2.0	- - LL2-SURD-H2 LL2-SURD20-H2
LL-Link / Bridge	GH : 1.3 GH : 2.0	- - LL2-SURD-N2 LL2-SURD20-N2
CMFit Abutment(S) / $\varnothing 10$	CMFit10-DN30-H	CMFit10-SURD-H2
CMFit Abutment(B) / $\varnothing 10$	-	CMFit10-SURD-N2
CMFit Abutment(S) / $\varnothing 14$	-	CMFit14-SURD-H2
CMFit Abutment(B) / $\varnothing 14$	-	CMFit14-SURD-N2
3D-Analog		3D-SURO
Screw(5ea)	AC-DN305	AC-SURD5

DIO SM

DIO UF / UFII

# Compatible to Implant Systems

NeoBiotech IS

## NeoBiotech IS

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm) DERN (7.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.2Hex
				30Nm Ø 3.5/4.0/4.5/5.0

Description-Library code		0501 - DERN	
Product		Connection Platform Ø	Standard Ø 3.5/4.0/4.5/5.0
SC-Scanbody			SC-SURO
GMS-Scanbody(Short)			GMS-DERNS
GMS-Scanbody(Long)			GMS-DERNL
LL-Multibase (Single)	GH : 1.3		LL2-DERN-H2
	GH : 2.0		LL2-DERN20-H2
LL-Multibase (Single)	GH : 1.3		LL2-DERN-N2
	GH : 2.0		LL2-DERN20-N2
LS-Multibase / Single		-	
CMFit Abutment(S) / Ø10		CMFit10-DERN-H2	
CMFit Abutment(B) / Ø10		CMFit10-DERN-N2	
CMFit Abutment(S) / Ø14		CMFit14-DERN-H2	
CMFit Abutment(B) / Ø14		CMFit14-DERN-N2	
3D-Analog		3D-DER	
Screw(5ea)		AC-DERN5	

# Compatible to Implant Systems

## Dentis S-Clean / One Q

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Link
	(11.5mm)  (7.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.25Hex 30Nm

Description-Library code		0601- DERT	
Product	Connection Platform Ø	Regular	Ø 3.7/4.2/4.7/5.2/6.0
SC-Scanbody		SC-DER	
GMS-Scanbody(Short)		GMS-DERS	
GMS-Scanbody(Long)		GMS-DERL	
LL-Multibase (Single)	GH : 1.3	LL2-DERT13-H2	
	GH : 2.0	-	
LL-Multibase (Bridge)	GH : 1.3	LL2-DERT13-N2	
	GH : 2.0	-	
LS-Multibase / Single		-	
CMFit Abutment(S) / Ø10		CMFit10-DERT-H2	
CMFit Abutment(S) / Ø14		CMFit14-DERT-H2	
3D-Analog		3D-DER	
Screw(5ea)		AC-DERT5	

# Compatible to Implant Systems

## Straumann Bone Level

### Straumann Bone Level

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				Torx
				35Ncm

Description-Library code		0700 - BLS	0701 - BLM	0702 - BLR
Product	Connection Platform Ø	SC  Ø2.9	NC  Ø3.3	RC  Ø4.1/4.8
SC-Scanbody		-	SC-BLM	SC-BLR
GMS-Scanbody(Short)		-	GMS-BLMS	GMS-BLRS
GMS-Scanbody(Long)		GMS-BLSL	GMS-BLML	GMS-BLRL
LL-Multibase (Single)	GH : 1.3	LL2-BLS-R	LL2-BLM-R2	LL2-BLR-R2
	GH : 2.0	-	LL2-BLM20-R2	LL2-BLR20-R2
LL-Multibase (Bridge)	GH : 1.3	-	LL2-BLM-N2	LL2-BLR-N2
	GH : 2.0	-	LL2-BLM20-N2	LL2-BLR20-N2
LS-Multibase / Single		-	LS2-BLM-R2	LS2-BLR-R2
CMFit Abutment / Ø10		CMFit10-BLS-R	CMFit10-BLM-R2	CMFit10-BLR-R2
CMFit Abutment / Ø14		-	CMFit14-BLM-R2	CMFit14-BLR-R2
3D-Analog		3D-BLS	3D-BLM	3D-BLR
Screw(5ea)		AC-BLS5	AC-BLM5	AC-BLR5

# Compatible to Implant Systems

## Straumann Tissue Level

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(10mm)  (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) Torx 35Nm

Description-Library code	0703 - IT35	0704 - ITR	0705 - ITW
Product Connection Platform Ø	NN  Ø3.5	RN  Ø4.8	WN  Ø 6.5
SC-Scanbody	SC-IT35	SC-ITRW	SC-ITRW
GMS-Scanbody(Short)	-	GMS-ITRS	GMS-ITWS
GMS-Scanbody(Long)	-	GMS-ITRL	GMS-ITWL
LL-Multibase (Single)	GH : 1.3 GH : 2.0	LL2-ITR-02 LL2-ITR20-02	LL2-ITW-02 LL2-ITW20-02
LL-Multibase (Bridge)	GH : 1.3 GH : 2.0	LL2-ITR-N2 LL2-ITR20-N2	LL2-ITW-N2 LL2-ITW20-N2
LS-Multibase / Single	LS2-IT35-0	LS2-ITR-02	LS2-ITW-02
CMFit Abutment / Ø10	CMFit10-IT35-0	CMFit10-ITR-02	CMFit10-ITW-02
CMFit Abutment / Ø14	-	CMFit14-ITR-02	CMFit14-ITW-02
3D-Analog	-	3D-ITR	3D-ITW
Screw(5ea)	AC-IT355	AC-IT5	AC-IT5

## Straumann BLX

\*set (Scanbody/Multibase/CMFit Abutment + screw)

GMS-Scanbody (Intra-Oral)	Single LL-Multibase	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) Torx 35Nm
(10mm)					

Description-Library code	0710 - BLX
Product Connection Platform Ø	Ø 3.5 - 6.5
GMS-Scanbody(Long)	GMS-BLXL
LL-Multibase / Single	LL2-BLX-T
CMFit Abutment / Ø14	CMFit14-BLX-T
3D-Analog	3D-BLX
Screw(5ea)	AX-BLX5

# Compatible to Implant Systems

## Nobel Biocare Active

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm)  (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				 unigrip
				15Nm Ø3.0
				35Nm Ø3.5/4.3/5.0/6.0

Description-Library code		0801- AC30	0802- ACN	0803- ACR	0804- ACW
Product	Connection Platform Ø		NP	RP	WP
Scanbody		OS-AC30	SC-ACN	SC-ACR	SC-ACW
GMS-Scanbody(Short)		-	GMS-ACNS	GMS-ACRS	GMS-ACWS
GMS-Scanbody(Long)		-	GMS-ACNL	GMS-ACRL	GMS-ACWL
LL-Multibase (Single)	GH : 1.3	LL2-AC30-T2	LL2-ACN-H2	LL2-ACR-H2	LL2-ACW-H2
	GH : 2.0	-	LL2-ACN20-H2	LL2-ACR20-H2	-
LL-Multibase (Bridge)	GH : 1.3	-	LL2-ACN-N2	LL2-ACR-N2	LL2-ACW-N2
	GH : 2.0	-	LL2-ACN20-N2	LL2-ACR20-N2	-
LS-Multibase / Single		LS2-AC30-T2	LS2-ACN-H2	LS2-ACR-H2	-
CMFit Abutment / Ø10		CMFit10-AC30-T2	CMFit10-ACN-H2	CMFit10-ACR-H2	CMFit10-ACW-H2
CMFit Abutment / Ø14		-	CMFit14-ACN-H2	CMFit14-ACR-H2	CMFit14-ACW-H2
3D-Analog		3D-AC30	3D-ACN	3D-ACR	3D-ACW
Screw(5ea)		AC-AC305	AC-ACN5	AC-ACR5	AC-ACW5

# Compatible to Implant Systems

## Nobel Biocare Branemark

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)  35Nm

Description-Library code	0805 - EXMB	0806 - EXRB	0807 - EXWB
Product \ Connection Platform Ø	NP  Ø3.5	RP  Ø4.1	WP  Ø5.1
SC-Scanbody	SC-EXMB	SC-EXRB	SC-EXWB
GMS-Scanbody(Short)	-	-	-
GMS-Scanbody(Long)	-	-	-
LL-Multibase / Single	LL2-EXMB-H2	LL2-EXRB-H2	LL2-EXWB-H2
LL-Multibase / Bridge	LL2-EXMB-N2	LL2-EXRB-N2	LL2-EXWB-N2
LS-Multibase / Single	LS2-EXMB-H2	LS2-EXRB-H2	LS2-EXWB-H2
LS-Multibase / Bridge	-	-	-
CMFit Abutment(S) / Ø10	CMFit10-EXMB-H2	CMFit10-EXRB-H2	CMFit10-EXWB-H2
CMFit Abutment(B) / Ø10	-	CMFit10-EXRB-N2	-
CMFit Abutment(S) / Ø14	CMFit14-EXMB-H2	CMFit14-EXRB-H2	CMFit14-EXWB-H2
3D-Analog	3D-EXM	3D-EXR	3D-EXW
Screw(5ea)	AC-EXMB5	AC-EXRB5	AC-EXWB5

# Compatible to Implant Systems

## Nobel Biocare Replace

### Nobel Biocare Replace

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm)  (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				unigrip 35Nm

Product \ Connection Platform Ø	0808- RPM	0809- RPR	0810- RPW	0811- RP60
NP				
Scanbody	SC-RPM	SC-RPR	SC-RPW	SC-RP60
GMS-Scanbody(Short)	GMS-RPMS	GMS-RPRS	GMS-RPWS	GMS-RP60S
GMS-Scanbody(Long)	GMS-RPML	GMS-RPRL	GMS-RPWL	GMS-RP60L
LL-Multibase / Single	LL2-RPM-E2	LL2-RPR-E2	LL2-RPW-E2	LL2-RP60-E2
LL-Multibase / Bridge	LL2-RPM-N2	LL2-RPR-N2	LL2-RPW-N2	LL2-RP60-N2
LS-Multibase / Single	LS2-RPM-E2	LS2-RPR-E2	LS2-RPW-E2	LS2-RP60-E2
CMFit Abutment / Ø10	CMFit10-RPM-E2	CMFit10-RPR-E2	CMFit10-RPW-E2	CMFit10-RP60-E2
CMFit Abutment / Ø14	CMFit14-RPM-E2	CMFit14-RPR-E2	CMFit14-RPW-E2	CMFit14-RP60-E2
3D-Analog	3D-RPM	3D-RPR	3D-RPW	3D-RP60
Screw(5ea)	AC-RPM5	AC-RPR5	AC-RPW5	AC-RP60

# Compatible to Implant Systems

## Dentsply Xive/Frialit2

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) 1.22Hex 24Nm

Description-Library code	0901 - XV3	0902 - XVM	0903 - XVR	0904 - XVW	0905 - XV55
Product \ Connection Platform Ø	Ø3.0	Ø3.4	Ø3.8	Ø4.5	Ø5.5
SC-Scanbody	SC-XV3	SC-XVM	SC-XVR	SC-XWV	SC-XWV
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-XV3-H2	LL2-XVM-H2	LL2-XVR-H2	LL2-XWV-H2	LL2-XW55-H2
LL-Multibase / Bridge	LL2-XV3-N2	LL2-XVM-N2	LL2-XVR-N2	LL2-XWV-N2	LL2-XW55-N2
LS-Multibase / Single	LS2-XV3-H2	LS2-XVM-H2	LS2-XVR-H2	LS2-XWV-H2	LS2-XW55-H2
CMFit Abutment / Ø10	CMFit10-XV3-H2	CMFit10-XVM-H2	CMFit10-XVR-H2	CMFit10-XWV-H2	CMFit10-XW55-H2
CMFit Abutment / Ø14	-	CMFit14-XVM-H2	CMFit14-XVR-H2	CMFit14-XWV-H2	CMFit14-XW55-H2
3D-Analog	3D-XV3	3D-XVM	3D-XVR	3D-XWV	3D-XW55
Screw(5ea)	AC-XV35	AC-XV5	AC-XV5	AC-XV5	AC-XV5

# Compatible to Implant Systems

Dentsply Ankylos

## Dentsply Ankylos

\*set (Scanbody/Multibase/CMFit Abutment + screw)

OS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase	
			-	
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			 1.0Hex 15Nm

Description-Library code	0906 - AKLL	
Connection Platform Ø	universal  Ø 3.4/4.5/5.5	
Product	OS-Scanbody	
LL-Multibase (Single)	GH : 1.3	LL2-AKLL-T
	GH : 2.0	LL2-AKLL20-T
LL-Multibase (Bridge)	GH : 1.3	LL2-AKLL-N
	GH : 2.0	LL2-AKLL20-N
LS-Multibase / Single	-	
CMFit Abutment(S) / Ø10	CMFit10-AKLL-T	
CMFit Abutment(B) / Ø10	CMFit10-AKLL-N	
3D-Analog	3D-AKLL	
Screw	-	

# Compatible to Implant Systems

## Biomet 3i Certain

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.2Hex 20Ncm

Description-Library code	1001 - 3ICM	1002 - 3ICR	1003 - 3ICW	1004 - 3IC60
Product	Connection Platform Ø	Ø3.4	Ø4.1	Ø5.0
SC-Scanbody	SC-3ICM	SC-3ICR	SC-3ICR	SC-3ICR
GMS-Scanbody(Short)	GMS-3ICMS	GMS-3ICRS	GMS-3ICRS	GMS-3ICRS
GMS-Scanbody(Long)	GMS-3ICML	GMS-3ICRL	GMS-3ICRL	GMS-3ICRL
LL-Multibase / Single	LL2-3ICM-H2	LL2-3ICR-H2	LL2-3ICW-H2	LL2-3IC60-H2
LL-Multibase / Bridge	LL2-3ICM-N2	LL2-3ICR-N2	LL2-3ICW-N2	LL2-3IC60-N2
LS-Multibase / Single	LS2-3ICM-H2	LS2-3ICR-H2	LS2-3ICW-H2	LS2-3IC60-H2
CMFit Abutment / Ø10	CMFit10-3ICM-H2	CMFit10-3ICR-H2	CMFit10-3ICW-H2	CMFit10-3IC60-H2
CMFit Abutment / Ø14	-	CMFit14-3ICR-H2	CMFit14-3ICW-H2	CMFit14-3IC60-H2
3D-Analog	3D-3ICM	3D-3ICR	3D-3ICW	3D-3IC60
Screw(Sea)	AC-3IC5	AC-3IC5	AC-3IC5	AC-3IC5

# Compatible to Implant Systems

## Biomet 3i Osseotite

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm)  (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)  1.2Hex  20Ncm

Description-Library code	1005- EX3M	1006- EXR	1007- EXT	1008- EXTW
Connection Platform Ø				
Product	Ø3.4	Ø4.1	Ø5.0	Ø6.0
Scanbody	SC-EX3M	SC-EXR	SC-EXT	SC-EXTW
GMS-Scanbody(Short)	-	GMS-EXRS	GMS-EXRS	GMS-EXRS
GMS-Scanbody(Long)	-	GMS-EXRL	GMS-EXRL	GMS-EXRL
LL-Multibase / Single	LL2-EX3M-H2	LL2-EXR-H2	LL2-EXT-H2	LL2-EXTW-H2
LL-Multibase / Bridge	LL2-EX3M-N2	LL2-EXR-N2	LL2-EXT-N2	LL2-EXTW-N2
LS-Multibase / Single	LS2-EX3M-H2	LS2-EXR-H2	LS2-EXT-H2	LS2-EXTW-H2
LS-Multibase / Bridge	-	-	-	-
CMFit Abutment(S) / Ø10	CMFit10-EX3M-H2	CMFit10-EXR-H2	CMFit10-EXT-H2	CMFit10-EXTW-H2
CMFit Abutment(B) / Ø10	-	CMFit10-EXR-N2	-	-
CMFit Abutment(S) / Ø14	-	CMFit14-EXR-H2	CMFit14-EXT-H2	CMFit14-EXTW-H2
3D-Analog	-	3D-EXR	-	-
Screw(5ea)	AC-EXR5	AC-EXR5	AC-EXR5	AC-EXR5

# Compatible to Implant Systems

## Implant Direct 3.0

\*set (Scanbody/Multibase/CMFit Abutment + screw)

OS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase	
		-	-	
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-		 	 1.25Hex 30Nm

Description-Library code	1101 - TSVM
Connection Platform Ø	Implant Direct  Ø3.0
Product	
OS-Scanbody	OS-TSVM
LL-Multibase / Single	LL2-TSVM-H2
LL-Multibase / Bridge	-
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-TSVM-H2
CMFit Abutment / Ø14	-
3D-Analog	3D-TSVM
Screw(5ea)	AC-TSVM5

# Compatible to Implant Systems

## Zimmer Tapered Screw-vent

### Zimmer Tapered Screw-vent

\*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Link	Bridge LL-Link	Single LS-Link
	(11.5mm)  (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.25Hex 30Ncm

Description-Library code		1201- TSV3	1202- TSV4	1203- TSV5
Product	Connection Platform Ø			
Scanbody		SC-TSV3	SC-TSV4	SC-TSV5
GMS-Scanbody(Short)		GMS-TSV3S	GMS-TSV4S	GMS-TSV5S
GMS-Scanbody(Long)		GMS-TSV3L	GMS-TSV4L	GMS-TSV5L
LL-Link / Single	GH : 1.0	LL2-TSV3-H2	LL2-TSV4-H2	LL2-TSV5-H2
	GH : 2.0	LL2-TSV320-H2	LL2-TSV420-H2	-
LL-Link / Bridge	GH : 1.0	LL2-TSV3-N2	LL2-TSV4-N2	LL2-TSV5-N2
	GH : 2.0	LL2-TSV320-N2	LL2-TSV420-N2	-
LS-Link / Single		LS2-TSV3-H2	LS2-TSV4-H2	LS2-TSV5-H2
CMFit Abutment / Ø10		CMFit10-TSV3-H2	CMFit10-TSV4-H2	CMFit10-TSV5-H2
CMFit Abutment / Ø14		CMFit14-TSV3-H2	CMFit14-TSV4-H2	CMFit14-TSV5-H2
3D-Analog		3D-TSV3	3D-TSV4	3D-TSV5
Screw(5ea)		AC-TSV5	AC-TSV5	AC-TSV5

# Compatible to Implant Systems

## Zimmer Spline

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
-		-		 1.25Hex 30Nm

Description-Library code	1204 - SP40	1205- SP50
Product \ Connection Platform Ø	Ø 3.7 / 4.0	Ø 5.0
SC-Scanbody	SC-SP40	SC-SP40
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Multibase / Single	LL2-SP40-H	LL2-SP50-H
LL-Multibase / Bridge	-	-
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	-	-
CMFit Abutment / Ø14	CMFit14-SP40-H	CMFit14-SP50-H
3D-Analog	-	-
Screw(5ea)	AC-TSV5	AC-TSV5

Zimmer Spline

# Compatible to Implant Systems

## Astra Tech Osseospeed TX

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm)  (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.27Hex
				15Ncm Ø3.0
				20Ncm Ø 3.5/4.0
				25Ncm Ø 4.5/5.0

Description-Library code		1301- AS3	1302- ASR	1303- ASW
Product	Connection Platform Ø	Yellow  Ø3.0	Aqua  Ø3.5/4.0	Lilac  Ø4.5/5.0
Scanbody		SC-AS3	SC-ASR	SC-ASW
GMS-Scanbody(Short)		-	GMS-ASRS	GMS-ASWS
GMS-Scanbody(Long)		-	GMS-ASRL	GMS-ASWL
LL-Multibase (Single)	GH : 1.3	LL2-AS3-T2	LL2-ASR-H2	LL2-ASW-H2
	GH : 2.0	-	LL2-ASR20-H2	LL2-ASW20-H2
LL-Multibase (Bridge)	GH : 1.3	-	LL2-ASR-N2	LL2-ASW-N2
	GH : 2.0	-	LL2-ASR20-N2	LL2-ASW20-N2
LS-Multibase / Single		-	LS2-ASR-H2	LS2-ASW-H2
CMFit Abutment / Ø10		CMFit10-AS3-T2	CMFit10-ASR-H2	CMFit10-ASW-H2
CMFit Abutment / Ø14		-	CMFit14-ASR-H2	CMFit14-ASW-H2
3D-Analog		-	3D-ASR	3D-ASW
Screw(5ea)		AC-AS35	AC-ASR5	AC-ASW5

# Compatible to Implant Systems

## Astra Tech Osseospeed EV

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.27Hex  25Ncm

Description-Library code	1304 - AEV30	1305 - AEV36	1306 - AEV42	1307 - AEV48	1308 - AEV54
Connection Platform Ø	Green	Purple	Yellow	Blue	Brown
SC-Scanbody	SC-AEV30	SC-AEV36	SC-AEV42	SC-AEV48	SC-AEV54
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-AEV30-T2	LL2-AEV36-T2	LL2-AEV42-T2	LL2-AEV48-T2	LL2-AEV54-T2
LL-Multibase / Bridge	LL2-AEV30-N2	LL2-AEV36-N2	LL2-AEV42-N2	LL2-AEV48-N2	LL2-AEV54-N2
LS-Multibase / Single	-	-	-	-	-
CMFit Abutment / Ø10	CMFit10-AEV30-T2	CMFit10-AEV36-T2	CMFit10-AEV42-T2	CMFit10-AEV48-T2	CMFit10-AEV54-T2
CMFit Abutment / Ø14	-	CMFit14-AEV36-T2	CMFit14-AEV42-T2	CMFit14-AEV48-T2	CMFit14-AEV54-T2
3D-Analog	3D-AEV30	3D-AEV36	3D-AEV42	3D-AEV48	3D-AEV54
Screw(5ea)	AC-AEV305	AC-AEV365	AC-AEV425	AC-AEV48545	AC-AEV48545

# Compatible to Implant Systems

Camlog

## Camlog

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.28Hex
				20Nm

Description-Library code	1401- CAM33	1402- CAM38	1403- CAM43	1404- CAM50	1405- CAM60
Product Connection Platform Ø					
Scanbody	SC-CAM33	SC-CAM38	SC-CAM43	SC-CAM5060	SC-CAM5060
LL-Multibase / Single	LL2-CAM33-T2	LL2-CAM38-T2	LL2-CAM43-T2	LL2-CAM50-T2	LL2-CAM60-T2
LL-Multibase / Bridge	LL2-CAM33-N2	LL2-CAM38-N2	LL2-CAM43-N2	LL2-CAM50-N2	LL2-CAM60-N2
LS-Multibase / Single	-	-	-	-	-
CMFit Abutment / Ø10	CMFit10-CAM33-T2	CMFit10-CAM38-T2	CMFit10-CAM43-T2	CMFit10-CAM50-T2	CMFit10-CAM60-T2
CMFit Abutment / Ø14	-	CMFit14-CAM38-T2	CMFit14-CAM43-T2	CMFit14-CAM50-T2	CMFit14-CAM60-T2
3D-Analog	3D-CAM33	3D-CAM38	3D-CAM43	3D-CAM50	-
Screw(5ea)	AC-CAM5	AC-CAM5	AC-CAM5	AC-CAM50605	AC-CAM50605

Conelog

## Conelog

\*set (Scanbody/Multibase/CMFit Abutment + screw)

GMS-Scanbody (Intra-Oral)	Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
					1.28Hex
					20Nm

Description-Library code	1406- CN33	1407- CN3843	1408- CN50
Product Connection Platform Ø			
GMS-Scanbody(Long)	GMS-CN33L	GMS-CN3843L	GMS-CN50L
CMFit Abutment / Ø10	CMFit10-CN33-E	-	-
CMFit Abutment / Ø14	-	CMFit14-CN3843-E	CMFit14-CN50-E
3D-Analog	3D-CN33	3D-CN3843	3D-CN50
Screw(5ea)	AC-CN5	AC-CN5	AC-CN505

# Compatible to Implant Systems

## MIS SEVEN

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm)  (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.25Hex
				30Nm

Description-Library code	1501 - MS33	1502 - MSR	1503 - MSW
Connection Platform Ø	Narrow  Ø3.3	Standard  Ø 3.75 / 4.2	Wide  Ø 5.0 / 6.0
Product			
SC-Scanbody	SC-MS33	SC-TSV3	SC-TSV4
GMS-Scanbody(Short)	-	GMS-TSV3S	GMS-TSV4S
GMS-Scanbody(Long)	-	GMS-TSV3L	GMS-TSV4L
LL-Multibase / Single	LL2-MS33-H	LL2-MSR-H	LL2-MSW-H
LL-Multibase / Bridge	LL2-MS33-N	LL2-MSR-N	LL2-MSW-N
LS-Multibase / Single	LS2-MS33-H	LS2-MSR-H	LS2-MSW-H
CMFit Abutment / Ø10	CMFit10-MS33-H	CMFit10-MSR-H	CMFit10-MSW-H
CMFit Abutment / Ø14	-	CMFit14-MSR-H	CMFit14-MSW-H
3D-Analog	3D-MS33	3D-TSV3	3D-TSV4
Screw(5ea)	AC-MS33	AC-MS5	AC-MS5

# Compatible to Implant Systems

MIS C1

## MIS C1

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.25Hex 30Nm

Description-Library code	1504- C1M	1505- C1R	1506- C1W
Product Connection Platform Ø	C1  Ø 3.3	C1  Ø 3.75/4.2	C1  Ø 5.0
Scanbody	SC-C1M	SC-C1R	SC-C1W
GMS-Scanbody(Short)	GMS-C1MS	GMS-C1RS	GMS-C1WS
GMS-Scanbody(Long)	GMS-C1ML	GMS-C1RL	GMS-C1WL
LL-Multibase / Single	LL2-C1M-T	LL2-C1R-T	LL2-C1W-T
LL-Multibase / Bridge	LL2-C1M-N	LL2-C1R-N	LL2-C1W-N
LS-Multibase / Single	-	-	-
CMFit Abutment / Ø10	CMFit10-C1M-T	CMFit10-C1R-T	CMFit10-C1W-T
CMFit Abutment / Ø14	CMFit14-C1M-T	CMFit14-C1R-T	CMFit14-C1W-T
3D-Analog	3D-C1M	3D-C1R	3D-C1W
Screw(5ea)	AC-C1M5	AC-C15	AC-C15

# Compatible to Implant Systems

## SIC ace / max

\*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Link	Bridge LL-Link	Single LS-Link
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			1.2Hex 20Ncm

Description-Library code	1601- SICM	1602- SICR
Connection Platform Ø	Ø 3.3	Ø 4.2
Product		
Scanbody	SC-SIC	SC-SIC
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Link / Single	LL2-SICM-H	LL2-SICR-H
LL-Link / Bridge	LL2-SICM-N	LL2-SICR-N
LS-Link / Single	-	-
CMFit Abutment / Ø10	CMFit10-SICM-H	CMFit10-SICR-H
CMFit Abutment / Ø14	CMFit14-SICM-H	CMFit14-SICR-H
3D-Analog	3D-SICM	3D-SICR
Screw(5ea)	AC-SIC5	AC-SIC5

SIC ace / max

# Compatible to Implant Systems

Thommen SPI

## Thommen SPI

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				Four lobe 15Ncm Ø3.5 25Ncm Ø 4.0/4.5/5.0/6.0

Description-Library code	1701- SPI35	1702- SPI40	1703- SPI45	1704- SPI50	1705- SPI60
Product \ Connection Platform Ø	Ø3.5	Ø4.0	Ø4.5	Ø5.0	Ø6.0
Scanbody	SC-SPI35	SC-SPI40	SC-SPI45	SC-SPI45	SC-SPI45
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-SPI35-H	LL2-SPI40-H	LL2-SPI45-H	LL2-SPI50-H	LL2-SPI60-H
LL-Multibase / Bridge	LL2-SPI35-N	LL2-SPI40-N	LL2-SPI45-N	LL2-SPI50-N	LL2-SPI60-N
LS-Multibase / Single	LS2-SPI35-H	LS2-SPI40-H	LS2-SPI45-H	LS2-SPI50-H	LS2-SPI60-H
LS-Multibase / Bridge	-	-	-	-	-
CMFit Abutment / Ø10	CMFit10-SPI35-H	CMFit10-SPI40-H	CMFit10-SPI45-H	CMFit10-SPI50-H	CMFit10-SPI60-H
CMFit Abutment / Ø14	-	CMFit14-SPI40-H	CMFit14-SPI45-H	CMFit14-SPI50-H	CMFit14-SPI60-H
3D-Analog	3D-SPI35	3D-SPI40	3D-SPI45	3D-SPI50	-
Screw(5ea)	AC-SPI35	AC-SPI40	AC-SPI45	AC-SPI50	AC-SPI55

# Compatible to Implant Systems

## BEGO

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.25Hex  30Ncm

Description-Library code	1801 - BEGO37	1802 - BEGO41	1803 - BEGO45	1804 - BEGO55
Product \ Connection Platform Ø				
SC-Scanbody	SC-BEGO	SC-BEGO	SC-BEGO	SC-BEGO
GMS-Scanbody(Short)	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-
LL-Multibase / Single	LL2-BEG037-H	LL2-BEG041-H	LL2-BEG045-H	LL2-BEG055-H
LL-Multibase / Bridge	LL2-BEG037-N	LL2-BEG041-N	-	-
LS-Multibase / Single	LS2-BEG037-H	LS2-BEG041-H	-	-
CMFit Abutment / Ø10	CMFit10-BEG037-H	CMFit10-BEG041-H	CMFit10-BEG045-H	CMFit10-BEG055-H
CMFit Abutment / Ø14	-	-	CMFit14-BEG045-H	CMFit14-BEG055-H
3D-Analog	3D-BEG037	-	-	-
Screw(5ea)	AC-BEG05	AC-BEG05	AC-BEG05	AC-BEG05

# Compatible to Implant Systems

MYTIS Arrow Implant

## MYTIS Arrow Implant

\*set (Scanbody/Link/CMFit Abutment + screw)

OS-Scanbody(Intra Oral)	Ø10 CMFit Abutment	Screw	Driver (Nm)
			 1.3Hex
			30Nm

Description-Library code	1901 - AR3740	1902 - AR4760
Product	Connection Platform Ø	Connection Platform Ø
OS-Scanbody	 Ø 3.7/4.0	 Ø 4.7/6.0
	OS-AR3740	OS-AR4760
	CMFit10-AR3740-H	CMFit10-AR4760-H
Screw(5ea)	AC-AR5	AC-AR5

Kyocera POI

## Kyocera POI

\*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	Single LL-Link	Bridge LL-Link	Ø10 CMFit Abutment	Screw	Driver (Nm)
		-			 1.1Hex
					20Nm

Description-Library code	2001 - POI37	2002 - POI42
Product	Connection Platform Ø	Connection Platform Ø
SC-Scanbody	 Ø 3.7	 Ø 4.2
LL-Link / Single	SC-POI37	SC-POI42
LL-Link / Bridge	LL2-POI37-H	LL2-POI42-H
CMFit Abutment / Ø10	CMFit10-POI37-H	CMFit10-POI42-H
Screw(5ea)	AC-POI375	AC-POI425

# Compatible to Implant Systems

## Biohorizons Laser Lok 3.0

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			 1.25Hex 30Ncm

Description-Library code	2101 - BH30
Connection Platform Ø	 Ø 3.0
Product	
OS-Scanbody	OS-BH30
LL-Multibase / Single	LL2-BH30-H
LL-Multibase / Bridge	-
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-BH30-H
CMFit Abutment / Ø14	-
3D-Analog	3D-BH30
Screw(5ea)	AC-BH305

# Compatible to Implant Systems

Anthogyr AXIOM

## Anthogyr AXIOM

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			 1.18Hex  25Nm

Description-Library code	2201 - ANG
Product	Connection Platform Ø universal  Ø 3.4/4.0/4.6/5.2
SC-Scanbody	SC-ANG
GMS-Scanbody(Short)	-
GMS-Scanbody(Long)	-
LL-Multibase / Single	LL2-ANG-E
LL-Multibase / Bridge	LL2-ANG-N
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-ANG-E
CMFit Abutment / Ø14	-
3D-Analog	3D-ANG
Screw(5ea)	AC-ANG5

# Compatible to Implant Systems

## Warantec Oneplant

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	$\varnothing 10$ CMFit Abutment	Single LS-Multibase	Screw	Driver (Nm)
		-		 1.2Hex  30Nm

Description-Library code	2301 - WOP		
Product	Connection Platform Ø	Standard	 Ø 4.3/4.8/5.3
Scanbody	SC-WOP		
LS-Multibase / Single	-		
CMFit Abutment $\varnothing 10$	CMFit10-WOP-0		
Screw(5ea)	AC-WOPS		

## Tekka In-Kone

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
		-	-	
		-	-	
$\varnothing 10$ CMFit Abutment	$\varnothing 14$ CMFit Abutment	3D-Analog	Screw	 1.2Hex  15Nm $\varnothing 3.5/4.0/4.5/5.5$
				

Description-Library code	2401 - TEK		
Product	Connection Platform Ø	Standard	 Ø 3.5/4.0/4.5/5.5
SC-Scanbody	SC-TEK		
GMS-Scanbody(Short)	-		
GMS-Scanbody(Long)	-		
LS-Multibase / Single	LS2-TEK-H		
CMFit Abutment / $\varnothing 10$	CMFit10-TEK-H		
CMFit Abutment / $\varnothing 14$	-		
3D-Analog	3D-TEK		
Screw(5ea)	AC-TEK5		

Warantec Oneplant

Tekka In-Kone

# Compatible to Implant Systems

Keystone PrimaConnex / Max-TL

## Keystone PrimaConnex / Max-TL

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.25 Square 30Ncm

Description-Library code	2501 - KEP35	2502 - KEP41	2503 - KEP50	2504 - KEX50	2505 - KEX60
Connection Platform Ø	Ø 3.5	Ø 4.1	Ø 5.0 / 6.0	Ø 7.0	Ø 8.0
Product	SC-KEP35	SC-KEP41	SC-KEP50	SC-KEP50	SC-KEP50
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-KEP35-T	LL2-KEP41-T	LL2-KEP50-T	LL2-KEX50-T	LL2-KEX65-T
LL-Multibase / Bridge	LL2-KEP35-N	LL2-KEP41-N	LL2-KEP50-N	LL2-KEX50-N	LL2-KEX65-N
LS-Multibase / Single	LS2-KEP35-T	LS2-KEP41-T	LS2-KEP50-T	-	-
CMFit Abutment / Ø10	CMFit10-KEP35-T	CMFit10-KEP41-T	CMFit10-KEP50-T	-	-
CMFit Abutment / Ø14	CMFit14-KEP35-T	CMFit14-KEP41-T	CMFit14-KEP50-T	CMFit14-KEX50-T	CMFit14-KEX65-T
3D-Analog	3D-KEP35	3D-KEP41	3D-KEP50	3D-KEX50	3D-KEX65
Screw(5ea)	AC-KEP5	AC-KEP5	AC-KEP5	AC-KEP5	AC-KEP5

# Compatible to Implant Systems

## Neodent CM

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-	-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				0.9Hex 15Ncm

Description-Library code	2601 - NCM	
Product	Connection Platform Ø	Standard  Ø 3.5/4.3/5.0
SC-Scanbody	SC-NCM	
LS-Multibase / Single	-	
CMFit Abutment / Ø10	CMFit10-NCM-H	
CMFit Abutment / Ø14	-	
3D-Analog	3D-NCM	
Screw	-	

## Neodent GM

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				unigrip (not compatible with Nobel) 20Ncm

Description-Library code	2602 - NCG	
Product	Connection Platform Ø	Standard  Ø 3.5/3.75/4.0/4.3/5.0
SC-Scanbody	SC-NCG	
CMFit Abutment / Ø10	CMFit10-NCG-H	
CMFit Abutment / Ø14	CMFit14-NCG-H	
LL-Multibase / Single	LL2-NCG-H	
LL-Multibase / Bridge	-	
3D-Analog	3D-NCG	
Screw(5ea)	AC-NCG5	

Neodent CM

Neodent GM

# Compatible to Implant Systems

Bicon Implant

## Bicon Implant

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	Single & Bridge LL-Multibase	Single & Bridge LS-Multibase	$\varnothing 10$ CMFit Abutment	$\varnothing 14$ CMFit Abutment
		-		-
Description-Library code	2701- BC20	2702- BC25	2703- BC30	
Product Connection Platform Ø	<input type="radio"/> Ø 2.0	<input type="radio"/> Ø 2.5	<input type="radio"/> Ø 3.0	
SC-Scanbody	SC-BC20	SC-BC25	SC-BC30	
LL-Multibase / <small>Single &amp; Bridge</small>	LL2-BC20-N	LL2-BC25-N	LL2-BC30-N	
LS-Multibase / <small>Single &amp; Bridge</small>	-	-	-	
CMFit Abutment / $\varnothing 10$	CMFit10-BC20-N	CMFit10-BC25-N	CMFit10-BC30-N	
CMFit Abutment / $\varnothing 14$	-	-	-	

# Compatible to Implant Systems

## Medentika IPS

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	$\varnothing 10$ CMFit Abutment	$\varnothing 14$ CMFit Abutment	Screw	Driver (Nm)
		-		 1.3 Hex
				25Nm

Description-Library code	2801- MTR
Product	Connection Platform Ø  Ø3.5/4.3/5.0
SC-Scanbody	SC-MTR
CMFit Abutment / Ø10	CMFit10-MTR-R
Screw(5ea)	AC-MTR5

Medentika IPS

## GC Aadvia Implant

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	$\varnothing 10$ CMFit Abutment	$\varnothing 14$ CMFit Abutment	Screw	Driver (Nm)
		-		 Torx
				20Nm

Description-Library code	2901- GAD
Product	Connection Platform Ø  Regular / Wide Ø4.0 / 5.0
SC-Scanbody	SC-GAD
CMFit Abutment / Ø10	CMFit10-GAD-H
Screw(5ea)	AC-GAD5

GC Aadvia Implant

# Compatible to Implant Systems

## BIOTECH KONTACT

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-		-	
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			 1.2 Hex 20Ncm

Description-Library code	3001- BK30	3002- BK40
Product	Connection Platform Ø	
SC-Scanbody	 Ø 3.0	 Ø 3.6 / 4.2 / 4.8 / 5.4
LL-Multibase / Single	LL2-BK30-T	LL2-BK40-T
LL-Multibase / Bridge	-	-
LS-Multibase / Single	LS2-BK30-T	LS2-BK40-T
CMFit Abutment / Ø10	CMFit10-BK30-T	CMFit10-BK40-T
3D-Analog	3D-BK30	3D-BK40
Screw(5ea)	AC-BK305	AC-BK405

# Compatible to Implant Systems

## Medentis medical ICX-System

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.42Hex 30Nm

Description-Library code	3101 - ICX
Product	Connection Platform Ø
SC-Scanbody	universal  Ø3.75/4.1/4.8
LL-Multibase / Single	LL2-ICX-H
LL-Multibase / Bridge	LL2-ICX-N
LS-Multibase / Single	LS2-ICX-H
CMFit Abutment / Ø10	CMFit10-ICX-H
CMFit Abutment / Ø14	CMFit14-ICX-H
3D-Analog	3D-ICX
Screw(5ea)	AC-ICX5

# Compatible to Implant Systems

## B&B Dental

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase			
-			-	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)			
	-			 1.27 Hex			
Description-Library code				25Nm			
Connection Platform Ø		 Ø 4.0/4.5/5.0					
GMS-Scanbody(long)	GMS-BNBL						
LL-Multibase / Single	LL2-BNB-H						
LL-Multibase / Bridge	-						
CMFit Abutment / Ø10	CMFit10-BNB-H						
3D-Analog	3D-BNB						
Screw(5ea)	AC-BNB5						

## Bredent Medical Sky


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase			
-			-	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)			
				 Torx			
Description-Library code				25Nm			
Connection Platform Ø		 Ø 3.5(NP) / 4.0-6.0					
GMS-Scanbody(long)	GMS-SKRL						
LL-Multibase / Single	LL2-SKR-T						
LL-Multibase / Bridge	-						
CMFit Abutment / Ø10	CMFit10-SKR-T						
CMFit Abutment / Ø14	CMFit14-SKR-T						
3D-Analog	3D-SKR						
Screw(5ea)	AC-SKR5						

\*set (Scanbody/Multibase/CMFit Abutment + screw)

# Compatible to Implant Systems

## NEOSS ProActive

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-		-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)   (not compatible with Nobel)
-				32Nm
Description-Library code		3401- NEP		
Product Connection Platform Ø	 Ø3.5/4.0/4.5/5.0/5.5			
	GMS-Scanbody(long)	GMS-NEPL		
LL-Multibase / Single		-		
LL-Multibase / Bridge		-		
CMFit Abutment / Ø10		-		
CMFit Abutment / Ø14		CMFit14-NEP-T		
3D-Analog		3D-NEP		
Screw(5ea)		AC-NEP5		

## Champions (R)evolution

\*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-		-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)   1.2 Hex
-				30Nm
Description-Library code		3501- CHR		
Product Connection Platform Ø	 Ø 3.5/4.5/5.5			
	GMS-Scanbody(long)	GMS-CHRL		
LL-Multibase / Single		-		
LL-Multibase / Bridge		-		
CMFit Abutment / Ø10		-		
CMFit Abutment / Ø14		CMFit14-CHR-H		
3D-Analog		3D-CHR		
Screw(5ea)		AC-CHR5		

# Compatible to Implant Systems

## MultiUnit(Multibase)

### Dentium Screw Abutment

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Nm)
		-		1.27Hex 10Nm

Description-Library code	0203- MUDE45	0204- MUDE55
Connection Platform Ø	Ø 4.5	Ø 5.5
Product		
SC-Scanbody	SC-MUDE	SC-MUDE
LB-Multibase / Bridge	LB-MUDE45	LB-MUDE55
3D Analog	-	-
Screw(Sea)	AC-MUDE5	AC-MUDE5

### Straumann Tissue Level SynOcta coping

\*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Nm)
		-		Torx 35Nm

Description-Library code	0706- MUIT
Connection Platform Ø	Ø 4.8
Product	
SC-Scanbody	SC-MUIT
LB-Multibase / Bridge	LB-MUIT
3D Analog	-
Screw(Sea)	AC-MUIT5

### Straumann Bonelevel SRA

\*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Nm)
				Torx 15Nm

Description-Library code	0707- MUSR35	0708- MUSR46
Connection Platform Ø	Ø 3.5	Ø 4.6
Product		
SC-Scanbody	SC-MUSR35	SC-MUSR46
LB-Multibase / Bridge	LB-MUSR35	LB-MUSR46
3D Analog	3D-MUSR35	3D-MUSR46
Screw(Sea)	AC-MUSR5	AC-MUSR5

# Compatible to Implant Systems

## MultiUnit (Multibase)

### Nobel MultiUnit 4.8

\*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Nm)
			 	 unigrip 15Nm

Description-Library code	0812- MUTB
Product	Connection Platform Ø
SC-Scanbody	Ø 4.8
LB-Multibase / Bridge(4.0mm)	SC-MUTB / SC2-MUTB
LB-Multibase / Bridge(5.5mm)	LB-MUTB2
3D Analog	LB-MUTBL2
Screw(5ea)	3D-MUTB
	AC-MUTB5

## Ankylos Multi-Unit

\*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Link	3D-Analog	Screw	Driver (Nm)
		-		 1.0Hex 10Nm

Description-Library code	0907- AKLM	0908- AKLR
Product	Connection Platform Ø	
SC-Scanbody	Mini	Regular
LB-Multibase / Bridge	SCB-16	SCB-16
3D Analog	LB-AKLM	LB-AKLR
Screw(5ea)	-	-
	AC-AKL5	AC-AKL5

## Zimmer Tapered Screw-vent Multi-Unit(TSV)

\*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Nm)
				 1.27Hex 20Nm



Description-Library code	1206- MUTSV
Product	Connection Platform Ø
SC-Scanbody	Universel
LB-Multibase / Bridge	SC-MUTSV
3D Analog	LB-MUTSV
Screw(5ea)	3D-MUTSV
	AC-MUTSV5

# GTS Transfer Scanbody

## GTS Transfer Scanbody

### Digital implant technology

Convenient & prompt digital implant impression  
Ideal subgingiva designing referenced from digital library



Differentiate the Gingival height by appearance



- Biocompatible Titanium (Ti-6Al-4V ELI) from CARPENTER Co guarantees durability which can be autoclavable and reusable.
- Provides various size with **different gingiva heights**.  
– 5mm(GH:2mm), 7mm(GH:4mm), 9mm(GH:6mm)
- Osstem TS and Dentium Super line connection available.
- **Highly precise Geo library** is provided for CAD s/w (3shape & exocad)
- **Special tool** allows to find right position and prevent detachment of screw
- **Customized Abutment Design** with CNC or Premilled abutment is available
- **Zirconia prothesis design** with Geo Link.

## Guide Tool



Guide Tool protects the scan body from any damage and helps to find dental implant connection easily.




GTS Scanbody is multifunctional that can be also used as a healing abutment and impression coping.

The groove on the scan body allows users to define different gingiva types and provide a bite index when taking an impression.

## Healing Abutment & Bite Index


# GTS Transfer scanbody

Direct Impression → Digital process



GTS Transfer Scanbody List

GTS Transfer scanbody



DER – Series (Dentium Super-line / Megagen AnyOne / Dentis S-clean / Neo IS)

Diameter GH / Height	4,0mm	4,5mm	5,5mm	6,5mm
2,0mm / 5,0mm	GTS-DER4020	GTS-DER4520	GTS-DER5520	GTS-DER6520
4,0mm / 7,0mm	GTS-DER4040	GTS-DER4540	GTS-DER5540	GTS-DER6540
6,0mm / 9,0mm	GTS-DER4060	GTS-DER4560	GTS-DER5560	GTS-DER6560
Guide Tool	GTST40	GTST45	GTST55	GTST65

\* Laser mark [DER55: Diameter 5,5mm]  
Define height and diameter by a different groove on the scan body.

SUM – Series (Osstem Implant TS Mini)

Diameter GH / Height	4,0mm	4,5mm
2,0mm / 5,0mm	GTS-SUM4020	GTS-SUM4520
4,0mm / 7,0mm	GTS-SUM4040	GTS-SUM4540
6,0mm / 9,0mm	GTS-SUM4060	GTS-SUM4560
Guide Tool	GTST40	GTST45

\* Laser mark [SUM45: Diameter 4,5mm]  
Define height and diameter by a different groove on the scan body.

SURO – Series (Osstem Implant TS Standard / Dio UF Regular)

Diameter GH / Height	4,5mm	5,0mm	6,0mm	7,0mm
2,0mm / 5,0mm	GTS-SURO4520	GTS-SURO5020	GTS-SURO6020	GTS-SURO7020
4,0mm / 7,0mm	GTS-SURO4540	GTS-SURO5040	GTS-SURO6040	GTS-SURO7040
6,0mm / 9,0mm	GTS-SURO4560	GTS-SURO5060	GTS-SURO6060	GTS-SURO7060
Guide Tool	GTST45	GTST50	GTST60	GTST70

\* Laser mark [SURO50: Diameter 5,0mm]  
Define height and diameter by a different groove on the scan body.

# Sirona Ti-base

Sirona Ti-base

## Sirona Multibase



### Description

	SHADE S2	SHADE S3
	S-MEZ0-S2	S-MEZ0-S3

### Scan Cap

Scan Cap (short)	Scan Cap (long)
SCSN-CAP	SCSNL-CAP

### Osstem Implant TS

		*set (Sirona Multibase + screw)									
		Platform Ø	Connection	TS Mini		Ø3.0/3.5	TS Standard		Ø4.0/4.5/5.0	Driver (Nm)	1.2Hex
	Single	Single [GH:1.2mm]		SRN2-SUM-H2			SRN2-SURO-H2			20Nm	Ø3.0/3.5
	Single	Single [GH:2.6mm]		-			SRN2-SUR026-H2			30Nm	Ø4.0/4.5/5.0

### Osstem Implant US / Biomet 3i Osseotite 4.1

		*set (Sirona Multibase + screw)						
		Platform Ø	Connection	US Regular		Ø4.1	Driver (Nm)	1.2Hex
	Single			SRN2-EXR-H2			30Nm	

# Sirona Ti-base

Sirona Ti-base

## Sirona Multibase

### Description

#### Dentium Implantium / SuperLine

		Platform Ø			Driver (Nm)	
		Connection	universal	Ø3.7/4.0/4.5	1.27Hex	
		Single [GH:1.2mm]	SRN2-DER-H2		25~30Nm	
		Single [GH:2.6mm]	SRN2-DER26-H2			

\*set (Sirona Multibase + screw)

#### Megagen AnyRidge

		Platform Ø			Driver (Nm)	
		Connection	Universal	Ø3.5/4.0/4.5	1.2Hex	
		Single	SRN2-MARR-H2		35Nm	

\*set (Sirona Multibase + screw)

#### Megagen AnyOne

		Platform Ø			Driver (Nm)	
		Connection	Universal	Ø3.5/4.0/4.5/5.0	1.2Hex	
		Single [GH:1.2mm]	SRN2-DERA-H2		25Nm	
		Single [GH:2.6mm]	SRN2-DERA26-H2			

\*set (Sirona Multibase + screw)

#### Neo Biotech IS

		Platform Ø			Driver (Nm)	
		Connection	Ø3.5/4.0/4.5	Driver (Nm)	1.2Hex	
		Single [GH:1.2mm]	SRN2-DERN-H2		30Nm	
		Single [GH:2.6mm]	SRN2-DERN26-H2			

\*set (Sirona Multibase + screw)

#### Straumann Bone Level

		Platform Ø			Driver (Nm)	
		Connection	NC Ø3.3	RC Ø4.1/4.8	Torx	
		Single	SRN2-BLM-R2		35Nm	
			SRN2-BLR-R2			

\*set (Sirona Multibase + screw)

# Sirona Ti-base

## Sirona Multibase

### Description

Sirona Ti-base

#### Straumann Tissue Level

Platform Ø Connection	RN	Ø4.8	WN	Ø6.5	Driver (Nm)  Torx	*set (Sirona Multibase + screw)
	Single	SRN2-ITR-02	SRN2-ITW-02	35Nm		
						

#### Nobel Biocare Active

Platform Ø Connection	NP	Ø3.5	RP	Ø4.3/5.0	Driver (Nm)  unigrip	*set (Sirona Multibase + screw)
	Single	SRN2-ACN-H2	SRN2-ACR-H2	35Nm		
						

#### Nobel Biocare Branemark

Platform Ø Connection	RP Ø4.1		Driver (Nm)  unigrip	*set (Sirona Multibase + screw)
	Single	SRN2-EXRB-H2		
				

#### Nobel Biocare Replace

Platform Ø Connection	NP	Ø3.5	RP	Ø4.3	WP	Ø5.0	Driver (Nm)  unigrip	*set (Sirona Multibase + screw)
	Single	SRN2-RPM-E2	SRN2-RPR-E2	SRN2-RPW-E2	35Nm			
								

#### Dentsply Xive / Frialit2

Platform Ø Connection		Ø3.4	Ø3.8	Ø4.5	Driver (Nm)  1.22Hex	*set (Sirona Multibase + screw)
	Single	SRN2-XVM-H2	SRN2-XVR-H2	SRN2-XVW-H2		
						

# Sirona Ti-base

## Sirona Multibase

### Description

#### Dentsply Ankylos

	Platform Ø Connection	universal Ø 3.4/4.5/5.5	*set (Sirona Multibase + screw) Driver (Nm) 1.0Hex
Single      Bridge	Single [GH:1.0mm]	SRN2-AKLL1-T	15Nm
	Single [GH:2.0mm]	SRN2-AKLL2-T	
	Bridge [GH:1.0mm]	SRN2-AKLL1-N	
	Bridge [GH:2.0mm]	SRN2-AKLL2-N	

Sirona Ti-base

#### Zimmer Tapered Screw-vent

	Platform Ø Connection	Ø3.5	Ø4.5	Ø5.7	*set (Sirona Multibase + screw) Driver (Nm) 1.25Hex
Single	Single	SRN2-TSV3-H2	SRN2-TSV4-H2	SRN2-TSV5-H2	30Nm

#### Astra Tech Osseospeed

	Platform Ø Connection	Aqua Ø3.5/4.0	Lilac Ø4.5/5.0	*set (Sirona Multibase + screw) Driver (Nm) 1.27Hex
Single	Single	SRN2-ASR-H2	SRN2-ASW-H2	20Nm Ø 3.5/4.0
				25Nm Ø 4.5/5.0

# CCM Casting Abutment

CCM Casting Abutment

## CCM Casting Abutment

Asia's first developed product in 2009  
 Six years of accumulated clinical data  
 Excellent biocompatibility  
 Incl. screw  
 Optimum platform design for emergence profile  
 Could make all type of implant restorations



### Description

#### Osstem Implant TS

		Platform Ø Connection		TS Mini  Ø3.0/3.5		TS Standard  Ø4.0/4.5/5.0		Driver (Ncm)  1.2Hex							
Single	Bridge	SUM-H2	GH : 1mm	SURO-H2		20Ncm	Ø3.0/3.5								
			GH : 2.2mm	SURO22-H2											
Bridge		SUM-N2	GH : 3mm	SURO30-H2		30Ncm	Ø4.0/4.5/5.0								
			GH : 1mm	SURO-N2											
			GH : 2.2mm	SURO22-N2											
				SURO30-N2											

#### Osstem Implant US

		Platform Ø Connection		US Mini  Ø3.5	US Regular  Ø4.1	US Wide  Ø5.1	Driver (Ncm)  1.2Hex	
Single	Bridge	Single	EXM-H2	EXR-H2	EXW-H2	30Ncm	1.2Hex	
		Bridge	EXM-N2	EXR-N2	EXW-N2			

#### Osstem Implant SS

		Platform Ø Connection		SS Regular  Ø4.8	SS Wide  Ø6.0	Driver (Ncm)  1.2Hex	
Single	Bridge	Single	ISR-02	ISW-02	30Ncm	1.2Hex	
		Bridge	ISR-N2	ISW-N2			

# CCM Casting Abutment

## Description

### Dentium SuperLine / Implantium

	Platform Ø Connection	Standard		Driver (Ncm)	*set (CCM Abt + screw)
		Ø3.7/4.0/4.5			
	Single	GH : 1mm	DER-H2	25~30Ncm	
		GH : 2mm	DER20-H2		
		GH : 3mm	DER30-H2		
	Bridge	GH : 1mm	DER-N2		
		GH : 2mm	DER20-N2		
		GH : 3mm	DER30-N2		

### Megagen Anyridge

	Platform Ø Connection	Standard		Driver (Ncm)	*set (CCM Abt + screw)
		Ø3.5/4.0/4.5			
	Single	MARR-H2		35Ncm	
	Bridge	MARR-N2			

### Megagen EZ-plus

	Platform Ø Connection	Ø 3.3	Ø 4.0/4.5/5.0	Driver (Ncm)	*set (CCM Abt + screw)
	Single	MEPM-E	MEPR-E	35Ncm	

### Megagen Ex-feel [Internal]

	Platform Ø Connection	Ø 4.8		Driver (Ncm)	*set (CCM Abt + screw)
			Ø 4.8		
	Single	ITRK-02		35Ncm	
	Bridge	ITRK-N2			

### DIO SM

	Platform Ø Connection	Ø 3.8/4.1	Ø 4.5/5.0	Driver (Ncm)	*set (CCM Abt + screw)
	Single	DSMM-T	DSMR-T	30Ncm	
	Bridge	DSMM-N	DSMR-N		

# CCM Casting Abutment

## Description

### DIO UF / UFII

		Platform Ø Connection	Ø 3.8/4.0/4.5/5.0	*set (CCM Abt + screw)
		Single	SURD-H2	30Ncm
		Bridge	SURD-N2	
Single	Bridge			

### Neo Biotech IS

		Platform Ø Connection	Ø 3.5/4.0/4.5	*set (CCM Abt + screw)
		Single	DERN-H2	30Ncm
		Bridge	DERN-N2	
Single	Bridge			

### Straumann Bone Level

		Platform Ø Connection	NC Ø3.3	RC Ø4.1/4.8	*set (CCM Abt + screw)
		Single	BLM-R2	BLR-R2	35Ncm
		Bridge	BLM-N2	BLR-N2	
Single	Bridge				

### Straumann Tissue Level

		Platform Ø Connection	RN Ø4.8	WN Ø 6.5	*set (CCM Abt + screw)
		Single	ITR-02	ITW-02	35Ncm
		Bridge	ITR-N2	ITW-N2	
Single	Bridge				

### Nobel Biocare Active

		Platform Ø Connection	NP Ø3.5	RP Ø 4.3/5.0	*set (CCM Abt + screw)
		Single	ACN-H2	ACR-H2	35Ncm
		Bridge	ACN-N2	ACR-N2	
Single	Bridge				

### Nobel Biocare Branemark

		Platform Ø Connection	NP Ø3.5	RP Ø 4.1	WP Ø5.1	*set (CCM Abt + screw)
		Single	EXMB-H2	EXRB-H2	EXWB-H2	35Ncm
		Bridge	EXMB-N2	EXRB-N2	EXWB-N2	
Single	Bridge					

# CCM Casting Abutment

## Description

### Nobel Biocare Replace

		Platform Ø Connection	NP		Ø3.5	RP		Ø4.3	WP		Ø5.0	Driver (Nm)	Driver unigrip	*set (CCM Abt + screw)
		Single		RPM-E2		RRP-E2		RPW-E2		RPW-N2		35Nm		
		Bridge		RPM-N2		RRP-N2		RPW-N2		RPW-N2				

### Nobel Biocare Multi Unit

		Platform Ø Connection			Ø4.8		Driver (Nm)	Driver unigrip	*set (CCM Abt + screw)
		Bridge		MUTB-N2			15Nm		

### Dentsply Xive / Frialit2

		Platform Ø Connection		Ø3.4		Ø3.8		Ø4.5	Driver (Nm)	Driver 1.22Hex	*set (CCM Abt + screw)
		Single	XVM-H2		XVR-H2		XVW-H2		24Nm		
		Bridge	XVM-N2		XVR-N2		XVW-N2				

### Dentsply Ankylos

		Platform Ø Connection		A/B/C/D		Ø3.4/4.5/5.5		Driver (Nm)	Driver 1.0Hex	*set (CCM Abt + screw)
		Single		-				25Nm		
		Bridge		DNF3-N						

### Biomet Bi Certain

		Platform Ø Connection		Ø3.4		Ø4.1		Driver (Nm)	Driver 1.2Hex	*set (CCM Abt + screw)
		Single	3ICM-H2		3ICR-H2			20Nm		
		Bridge	3ICM-N2		3ICR-N2					

### Biomet Bi Osseotite

		Platform Ø Connection		Ø4.1		Ø5.0		Driver (Nm)	Driver 1.2Hex	*set (CCM Abt + screw)
		Single	EXR-H2		EXT-H2			20Nm		
		Bridge	EXR-N2		EXT-N2					

### Zimmer Tapered Screw-vent

		Platform Ø Connection		Ø3.5		Ø4.5		Ø5.7	Driver (Nm)	Driver 1.25Hex	*set (CCM Abt + screw)
		Single	TSV3-H2		TSV4-H2		TSV5-H2		30Nm		
		Bridge	TSV3-N2		TSV4-N2		TSV5-N2				

# CCM Casting Abutment

## Description

### Zimmer Spline

		Platform Ø Connection	 Ø 3.7 / 4.0	Driver (Nm)
		Single	SP40-H	1.25Hex
				30Nm

### Astra Tech Osseospeed

		Platform Ø Connection	X-Small  Ø 3.0	Aqua  Ø 3.5/4.0	Lilac  Ø 4.5/5.0	Driver (Nm)
		Single	AS3-T2	ASR-H2	ASW-H2	1.27Hex
Single	Bridge		-	ASR-N2	ASW-N2	15Nm Ø 3.0 20Nm Ø 3.5/4.0 25Nm Ø 4.5/5.0

### MIS C1

		Platform Ø Connection	C1  Ø 3.3	C1  Ø 3.75/4.2	C1  Ø 5.0	Driver (Nm)
		Single	C1M-T	C1R-T	C1W-T	1.25Hex
Single	Bridge		C1M-N	C1R-N	C1W-N	30Nm

### SIC ace / max

		Platform Ø Connection	 Ø 3.3	 Ø 4.2	Driver (Nm)
		Single	SICM-H	SICR-H	1.2Hex
Single					20Nm

### BEGO

		Platform Ø Connection	 Ø 3.7	 Ø 4.1	Driver (Nm)
		Single	BEG037-H	BEG041-H	1.25Hex
Single	Bridge		BEG037-N	BEG041-N	30Nm

### Anthogyr AXIOM

		Platform Ø Connection	universal  Ø 3.4/4.0/4.6/5.2	Driver (Nm)
		Single	ANG-E	1.18Hex
Single	Bridge		ANG-N	25Nm

### Tekka In-Kone

		Platform Ø Connection	Standard  Ø 3.5/4.0/4.5/5.5	Driver (Nm)
		Single	TEK-H	1.2Hex
Single				15Nm

# Geo Screw Driver



## Angled Driver

Square Type / Ratchet Type

Use with LS Multibase Abutment

Maximum Angle 25°

Recommended torque allowance : below 30Ncm.(reference manufacture's guide)

Caution: unstable connection may bend or destroy the driver

Recommended usage (1.2hex/1.27hex)-10 times

(Universal Type)- 20 times

Connect with Square type Torque wrench

May use Round Type Torque wrench  
by connecting HDEC - Short (Clinic)Handle



Connect Straight & Angle driver  
with the Torque Wrench



## Straight Driver

Square Type / Ratchet Type

Connect by manufacturer's recommended torque allowance

Damage on the tip or usage of long period of time need to exchange driver

Connect with Square type Torque wrench

May use Round Type Torque wrench  
by connecting HDEC - Short (Clinic)Handle

## Universal - Round hexagon

Square Type / Ratchet Type

Connect by manufacturer's recommended torque allowance

Damage on the tip or usage of long period of time need to exchange driver

Connect with Square type Torque wrench

May use Round Type Torque wrench  
by connecting HDEC - Short (Clinic)Handle

## Angled Driver

### 3i & Osstem - 1.2mm Hex

	ASD3IL	Square -20mm		ASD3IR	Square -13mm		ALD3IL	Ratchet -20mm		ALD3IR	Ratchet -13mm
--	--------	--------------	--	--------	--------------	--	--------	---------------	--	--------	---------------

### Astra & Zimmer & Dentium - 1.27mm Hex

	ASDZIL	Square -20mm		ASDZIR	Square -13mm		ALDZIL	Ratchet -20mm		ALDZIR	Ratchet -13mm
--	--------	--------------	--	--------	--------------	--	--------	---------------	--	--------	---------------

### Universal - Round hexagon

	ASDUVL	Square -20mm		ASDUVR	Square -13mm		ALDUVL	Ratchet -20mm		ALDUVR	Ratchet -13mm
--	--------	--------------	--	--------	--------------	--	--------	---------------	--	--------	---------------

# Geo Screw Driver

## Straight Driver

### 3i - 1.2mm Hex

			
SSD3IL Square -20mm	SSD3IR Square -13mm	SLD3IL Ratchet -20mm	SLD3IR Ratchet -13mm
	SSD3IEL Square -35mm		

### Zimmer - 1.27mm Hex

			
SSDZIL Square -20mm	SSDZIR Square -13mm	SLDZIL Ratchet -20mm	SLDZIR Ratchet -13mm

### Astra - 1.27mm Hex

			
SSDATL Square -20mm	SSDATR Square -13mm	SLDATL Ratchet -20mm	SLDATR Ratchet -13mm
	SSDATEL Square -35mm		

### Nobel - unigrip

			
SSDNBL Square -20mm	SSDNBR Square -13mm	SLDNBL Ratchet -20mm	SLDNBR Ratchet -13mm
	SSDNBEL Square -35mm		

### Straumann - Torx

			
SSDSTL Square -20mm	SSDSTR Square -13mm	SLDSTL Ratchet -20mm	SLDSTR Ratchet -13mm
	SSDSTEL Square -35mm		

# Geo Screw Driver

## Driver Tool

HDEL	Long (Lab) Handle		
HDEC	Short (Clinic) Handle		
GTW	Torque Wrench		 <b>Round Type</b> Torque Wrench & Handle Short & Driver   <b>Square Type</b> Torque Wrench & Driver

## Lab Driver Set

Description	Product		
	Size	code	Implant System
 1.2 Hex	Blue	Ossstem Implant / Dio Implant Biomet 3i / Xive	
 1.27 Hex	Red	Dentium /Astra Tech Zimmer	
 Unigrip	Yellow	Nobel Biocare	
 Torx	Green	Straumann	
 1.0Hex	Purple	Ankylos	



Lab Driver Set

# Instructions of CAD/CAM Abutment

< Procedure for intra-oral scan >



1. Scanbody insert



2. Scanning



3. Matching [Scanbody and Library]



4. Design - 3shape / Exocad



5. Sintering




6. Setting

< Procedure for model scan >




1. Scanbody insert






2. Scanning



3. Matching [Scanbody and Library]




4. Design - 3shape / Exocad




5. Sintering and model setting

# Instructions of CCM Casting Abutment

## < Screw Type manufacturing >




## < Customized Type manufacturing >




# Instructions of CCM Casting Abutment



## < Hybrid Type manufacturing>



## < Bar Overdenture manufacturing>



# Geo Library



Free Geo Library available for 3Shape and exocad  
High quality of stability and accuracy  
One scanbody can be used to design all type of implant restorations  
Only implant library including  
Ti-base, implant blank, implant bridge



## Benefits :

The Geo Library enables you to design individual prosthetic solutions for almost all current implant systems.  
The Geo Library is available 3Shape and Exocad libraries, providing 71 different models of titanium bases.  
The Geo Library is the dental optimized for CAD was designed as a function.  
You are able to fabricate individual milled structures in your own lab, for glued or one piece structures with maximum precision and individuality.



## Application :

Individual ceramic abutments on titanium bases

Individual one-piece titanium abutments

Occlusal screwable crowns and bridges

Bars, telescopes

Screw retained abutment

Customized Zirconia abutment

Customized titanium abutment



# Selection of Geo Library

## Selection of Library

**Library for Implant System**

- Geo\_0101\_SUM
- Geo\_0102\_SURO
- Geo\_0103\_EXM
- Geo\_0104\_EXR

**Geo product code for different type**

**Cement gap option with correspond Scanbody**

Cement [gap]-Zirconia
027: SC
035: SC
042: SC
052: SC
027: GMS-S
035: GMS-S
045: GMS-S
052: GMS-S
027: GMS-L
035: GMS-L
042: GMS-L
052: GMS-L
08: PMMA-SC
10: PMMA-SC
08: PMMA-GMS-S
10: PMMA-GMS-S
08: PMMA-GMS-L
10: PMMA-GMS-L

**Cement [gap]-PMMA,Metal**

Cement [gap]-PMMA,Metal
027: SC
034: SC
042: SC
052: SC
027: GMS_S
034: GMS_S
042: GMS_S
052: GMS_S
027: GMS_L
042: GMS_L
052: GMS_L
08: PMMA_SC
10: PMMA_SC
08: PMMA_GMS_S
10: PMMA_GMS_S
08: PMMA_GMS_L
10: PMMA_GMS_L

**Wizard Detect Implant Position**

ONE-CLICK REGISTRATION  
SELECT IMPLANT FOR TOOTH 21  
Geo\_0101 SUM  
CMFit  
10mm SC

Click marked point on scan  
Best fit matching

Back Next

Wizard Detect Implant Position  
ONE-CLICK REGISTRATION  
SELECT IMPLANT FOR TOOTH 21  
Geo\_0101 SUM  
LL\_C20\_H43  
027 SC


Click marked point  
Best fit matching

Back Next

Wizard Detect Implant Position  
ONE-CLICK REGISTRATION  
SELECT IMPLANT FOR TOOTH 21  
Geo\_0101 SUM  
LL\_C20\_H43  
027 SC

Click marked point  
Best fit matching

Back Next



# Geo BS Library- Bridge&Single

**GEOMEDI**

**FOLLOW-ME!**  
TECHNOLOGY GROUP  
**hyperDENT**  
COMPACT CLASSIC OPTIONS



Brand	Implant System	System code	Bridge(Bar)	Single	Scanbody	Screw	
					code	Geo screw code	Compatible with original
Osstem	TS M	SUM	√	√	SC-SUM	SUMC	√
	TS Standard	SURO	√	√	SC-SURO	SUROC	√
	US M	EXM	√	√	SC-EXM	EXMC	√
	US R	EXR	√	√	SC-EXR	EXRC	√
	US W	EXW	√	√	SC-EXW	EXWC	√
Dentium	Superline	DER	√	√	SC-DER	DERC	√
	Screw-Abutment	MUDE45	√		SC-MUDE	MUDEC	√
		MUDE55	√		SC-MUDE	MUDEC	√
Straumann	Bonelevel 3.3	BLM	√	√	SC-BLM	BLMC	√
	Bonelevel 4.1	BLR	√	√	SC-BLR	BLMC	√
	Bonelevel 4.8	BLR	√	√	SC-BLR	BLRC	√
	Tissue level 4.8	ITR	√	√	SC-ITRW	bridge - ITTC	Temporary screw
						single - ITC	√
	Tissue level 6.5	ITW	√	√	SC-ITRW	bridge - ITTC	Temporary screw
						single - ITC	√
	SynOcta coping	MUIT	√		SC-MUIT	MUITC	√
	SRA Abutment	MUSR35	√		SC-MUSR35	MUSRC	√
		MUSR46	√		SC-MUSR46	MUSRC	√
Nobel	Multi Unit	MUTB	√		SC-MUTB	MUTBC	√
	Active 3.5	ACN	√	√	SC-ACN	ACNC	√
	Active 4.3	ACR	√	√	SC-ACR	ACRC	√
	Active 5.0	ACR	√	√	SC-ACR	ACRC	√
	Active 5.5	ACW	√	√	SC-ACW	ACRC	√
	Branemark 3.4	EXM	√	√	SC-EXM(B)	EXMBC	√
	Branemark 4.1	EXR	√	√	SC-EXR(B)	EXRBC	√
	Branemark 5.1	EXW	√	√	SC-EXW(B)	EXWBC	√
	Replace 3.5 rpm	RPM	√	√	SC-RPM	RPMC	√
	Replace 4.3 rpm	RPR	√	√	SC-RPR	RPRC	√
	Replace 5.0 rpm	RPW	√	√	SC-RPW	RPRC	√
	Replace 6.0 rpm	RP60	√	√	SC-RP60	RPRC	√
Dentsply	Xive 3.0	XV3	√	√	SC-XV3	XV3C	√
	Xive 3.4	XVM	√	√	SC-XVM	XVC	√
	Xive 3.8	XVR	√	√	SC-XVR	XVC	√
	Xive 4.5	XVW	√	√	SC-XVW	XVC	√
	Xive 5.5	XV55	√	√	SC-XVW	XVC	√
Ankylos	AKII	AKLL	√		OS-AKLL	DNF3C	special screw

# Geo BS Library- Bridge&Single

Brand	Implant System	System code	Bridge(Bar)	Single	Scanbody	Screw	
					code	GEO screw code	Compatible with original
3i	Certain 3.4	BICM	✓	✓	SC-BICM	3ICC	v
	Certain 4.1	BICR	✓	✓	SC-BICR	3ICC	v
	Certain 5.0	BICW	✓	✓	SC-BICR	3ICC	v
	Certain 6.0	BIC60	✓	✓	SC-BICR	3ICC	v
	External 3.4	EX3M	✓	✓	SC-EX3M	EXRC	v
	External 4.1	EXR	✓	✓	SC-EXR	EXRC	v
	External 5.0	EXT	✓	✓	SC-EXR	EXRC	v
	External 6.0	EXTW	✓	✓	SC-EXR	EXRC	v
Zimmer	TSV 3.5	TSV3	✓	✓	SC-TSV3	TSVC	v
	TSV 4.5	TSV4	✓	✓	SC-TSV4	TSVC	v
	TSV 5.7	TSV4	✓	✓	SC-TSV5	TSVC	v
Astra	Astra EV 3.0	AEV30	✓		SC-AEV30	AEV30C	v
	Astra EV 3.6	AEV36	✓		SC-AEV36	AEV36C	v
	Astra EV 4.2	AEV42	✓		SC-AEV42	AEV42C	v
	Astra EV 4.8	AEV48	✓		SC-AEV48	AEV4854C	v
	Astra EV 5.4	AEV54	✓		SC-AEV54	AEV4854C	v
	Astra 3.5	ASR	✓	✓	SC-ASR	ASRC	v
	Astra 4.0	ASR	✓	✓	SC-ASR	ASRC	v
	Astra 4.5	ASW	✓	✓	SC-ASW	ASWC	v
	Astra 5.0	ASW	✓	✓	SC-ASW	ASWC	v
Camlog	Camlog 3.3	CAM33	✓		SC-CAM33	CAMC	v
	Camlog 3.8	CAM38	✓		SC-CAM38	CAMC	v
	Camlog 4.3	CAM43	✓		SC-CAM43	CAMC	v
	Camlog 5.0	CAM50	✓		SC-CAM5060	CAM5060C	v
	Camlog 6.0	CAM60	✓		SC-CAM5060	CAM5060C	v
Magazen	AnyOne	MAO	✓	✓	SC-DER	MAOC	v
	Anyridge	MARR	✓	✓	SC-MARR	MARRC	v
Dentis	S-clean	DET	✓	✓	SC-DER	DERC	v
Dio	UF	SURD	✓	✓	SC-SURO	SURDC	v
Neo	IS	DERN	✓	✓	SC-SURO	DERNC	v
MIS	C1 3.3	C1M	✓		SC-C1M	C1MC	v
	C1 3.75 / 4.2	C1R	✓		SC-C1R	C1C	v
	C1 5.0	C1W	✓		SC-C1W	C1C	v

## Geo Product List (Scanbody / Multibase / CMFit / 3D Analog / CCM Abutment / Sirona Ti-base)

No	Company	Implant system	Platform	Library Code	Type	sc(os) Scanbody	GMS Scanbody	LL-Link	LS-Link	CMFit 10mm	CMFit 14mm	3D Analog	Screw(5ea)	CCM Abutment (Co-Cr base)	Sirona Ti-base
101			3.0/3.5	0101-SUM	Single	SC-SUM	LL2-SUM-M12-H2	LS2-SUM-MH2	CMFIT10-SUM-H2	CMFIT14-SUM-H2	3-D-SUM	AC-SUM5	SUM-H2	SUM-N2	SRN2-SUM-H2
102	TS	4.0/4.5/5.0-	0102-SURO	Single	Bridge	GMS-SUML	LL2-SUM20-H2	LS2-SUM-N2	CMFIT10-SUM-N2	CMFIT14-SUM-N2	-	-	-	-	-
103	Oststem	3.5	0103-EXM	Single	Bridge	GMS-SUROS	LL2-SUR020-H2	LS2-SURO-H2	CMFIT10-SURO-H2	CMFIT14-SURO-H2	-	-	SURO22-H2	SURO30-H2	SURO26-H2
104	US	4.1	0104-EXR	Single	Bridge	GMS-EXMS	LL2-SUR013-H2	LS2-SURO-N2	CMFIT10-SURO-N2	CMFIT14-SURO-N2	3-D-SURO	AC-SUROS	SURO22-N2	SURO30-N2	SURO26-N2
105		5.1	0105-EXW	Single	Bridge	GMS-EXW	LL2-SUR020-H2	LS2-EXW-H2	CMFIT10-EXW-H2	CMFIT14-EXW-H2	-	-	-	-	-
106		4.8	0106-4SR	Single	Bridge	GMS-4SR	LL2-EXW-N2	LS2-EXW-N2	CMFIT10-EXW-N2	CMFIT14-EXW-N2	3-D-EXW	AC-EXW5	EXW-H2	EXW-N2	-
107	SS	6.0	0107-4SW	Single	Bridge	GMS-4SW	LL2-4SR-N2	-	CMFIT10-4SR-O2	CMFIT14-4SR-O2	3-D-4SR	AC-EXR5	ISR-H2	ISR-N2	-
201	Dentium	3.7/4.0/4.5 (Universal)	0201-DER	Single	Bridge	GMS-DERS	LL2-4SW-N2	-	CMFIT10-4SW-O2	CMFIT14-4SW-O2	3-D-4SW	AC-EXR5	ISW-H2	ISW-N2	-
202	NR-Line	3.1/3.6/4.3 (Standard)	0202-DNR	Single	Os-DNR	GMS-MARRS	LL2-4SW-N2	LS2-DER-H2	CMFIT10-4SW-H2	CMFIT14-4SW-H2	3-D-DER	AC-DERS	DER-H2	DER-N2	SRN2-4SW-H2
301	AnyRidge	3.5/4.0/4.5 (Universal)	0301-MARR	Single	Bridge	GMS-MARR	LL2-MARR-H2	LS2-MARR-H2	CMFIT10-MARR-H2	CMFIT14-MARR-H2	3-D-MARR	AC-MARR5	MARR-H2	MARR-N2	SRN2-4MARR-H2
302	Megagen	3.0/4.0/4.5/5.0 (Universal)	0302-DERA	Single	SO-DER	GMS-DERS	LL2-MARR20-H2	LL2-4MARR20-H2	-	-	-	-	-	-	SRN2-DERA-H2
303	EZ-Plus	3.3	0303-WEPM	Single	Bridge	GMS-DERL	LL2-4DERA13-H2	LL2-4DERA20-H2	CMFIT10-4DERA-H2	CMFIT14-4DERA-H2	3-D-DER	AC-MAR5	DER-H2	DER-N2	SRN2-DERA26-H2
304	EX-feel	4.0/4.5/5.0	0304-MEPR	Single	SC-MEPM	GMS-MARR	LL2-4DERA13-N2	LL2-4DERA20-N2	CMFIT10-4DERA-N2	CMFIT14-4DERA-N2	-	-	AC-MEPMS	MEPM-E	-
305		4.8 (Regular)	0305-ITRK	Single	Bridge	SC-ITRW	LL2-ITRK<20	LL2-ITRK<22	CMFIT10-ITRK<20	CMFIT14-ITRK<22	3-D-ITR	AC-EXR5	ITRK-C02	ITRK-N2	-
306		3.0/3.25	0306-MMN	Single	SC-MMN	GMS-MMN	LL2-4DERA30-H2	LL2-4DERA30-N2	CMFIT10-MMN-H	CMFIT14-MMN-N2	-	-	AC-MMNS	-	-
401	Dio	SM	3.8/4.1	0401-DSMM	Single	SC-DSMM	-	-	CMFIT10-DSMM-T	CMFIT14-DSMM-T	-	-	DSMM-T	DSMM-N	-
402		4.5/5.0	0402-DSMR	Single	SC-DSMR	GMS-DSMR	-	-	CMFIT10-DSMR-T	CMFIT14-DSMR-T	-	-	DSMR-T	DSMM-N	-

No	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Link	LS-Link	CMFit 10mm	CMFit 14mm	3D Analog	Screw(5ea)	CCM Abutment (Co-Cr base)	Sirona Ti-base
403	Dio	3.0/3.3 (Narrow)	0403-DN30	Single	SC-MMN	-	-	-	-	CMFIT10-DN30-H	-	-	-	AC-DN305	-
404	UF / UFII	3.8/4.0/4.5/5.0 (Standard)	0404-SURD	Single	SC-SURO	GMS-SUROS	LL2-SURD-H2	CMFIT10-SURD-H2	CMFIT14-SURD-H2	CMFIT10-SURD-H2	CMFIT14-SURD-H2	3D-SURD	AC-SURD5	SURD+H2	-
501	NeoBiotech	IS	0501-DERN	Single	SC-SURO	GMS-DERN	LL2-SURD20-N2	-	-	CMFIT10-SURD-N2	CMFIT14-SURD-N2	3D-SURD	AC-SURD5	SURD+N2	-
601	Dentis	S-Clean/OneQ	0601-DERT	Single	SC-DER	GMS-DERS	LL2-DERNH2	-	-	CMFIT10-DERNH2	CMFIT14-DERNH2	3D-DER	AC-DERNS	DERNH2	SRN2-DERNH2
709		2.9 (SC)	0700-BLS	Single	-	GMS-BLSL	LL2-BLS-R	-	-	CMFIT10-BLS-R	CMFIT14-BLS-R	3D-BLS	AC-BLS5	-	-
701	Straumann	BoneLevel	3.3 (NC)	0701-BLM	Single	SC-BLM	GMS-BLMS	LL2-BLM-N2	LS2-BLM-R2	CMFIT10-BLM-R2	CMFIT14-BLM-R2	3D-BLM	AC-BLM5	BLM-R2	SRN2-BLM-R2
702		4.1/4.8 (RC)	0702-BLR	Single	SC-BLR	GMS-BLRL	LL2-BLR-N2	LL2-BLR20-R2	LS2-BLR-R2	CMFIT10-BLR-R2	CMFIT14-BLR-R2	3D-BLR	AC-BLR5	BLR-R2	SRN2-BLR-R2
703		3.5 (NN)	0703-IT35	Single	SC-IT35	-	-	-	LS2-IT35-O	CMFIT10-IT35-O	-	-	AC-IT355	-	-
704	Straumann	Tissue Level	4.8 (RN)	0704-ITR	Single	GMS-ITRS	LL2-ITR20-O2	LS2-ITR-O2	CMFIT10-ITR-O2	CMFIT14-ITR-O2	3D-ITR	-	-	ITR-O2	SRN2-ITR-O2
705		6.5 (WN)	0705-ITW	Single	SC-ITRW	GMS-ITRW	LL2-ITR20-N2	LS2-ITR-N2	CMFIT10-ITR-N2	CMFIT14-ITR-N2	3D-ITR	AC-IT5	ITR-N2	-	-
710		BLX	3.6 - 6.5	0710-BLX	Single	GMS-BLXL	LL2-BLX-T	-	-	CMFIT10-BLX-T	CMFIT14-BLX-T	3D-BLX	AC-BLX5	-	-
801	Nobel Biocare	3.0	0801-AC30	Single	OS-AC30	-	LL2-AC30-T2	LS2-AC30-T2	CMFIT10-AC30-T2	CMFIT14-AC30-T2	3D-AC30	AC-AC305	-	-	
802		3.5 (NP)	0802-ACN	Single	SC-ACN	GMS-ACNS	LL2-ACNH2	LS2-ACNH2	CMFIT10-ACNH2	CMFIT14-ACNH2	3D-ACN	AC-ACN5	ACN-H2	SRN2-ACN-H2	
803		4.3/5.0 (RP)	0803-ACR	Bridge	SC-ACR	GMS-ACRS	LL2-ACRN2	LS2-ACRN2	CMFIT10-ACR-H2	CMFIT14-ACR-H2	3D-ACR	AC-ACR5	ACR-H2	SRN2-ACR-H2	
804		5.5 (WP)	0804-ACW	Single	SC-ACW	GMS-ACW	LL2-ACNW2	LS2-ACNW2	CMFIT10-ACNW2	CMFIT14-ACNW2	3D-ACW	-	-	-	
805	Branemark	3.5 (NP)	0805-EXMB	Single	SC-EXMB	-	LL2-EXMB-H2	LS2-EXMB-H2	CMFIT10-EXMB-H2	CMFIT14-EXMB-H2	3D-EXM	AC-EXMB5	EXMB-H2	-	
806		4.1 (RP)	0806-EXRB	Bridge	SC-EXRB	-	LL2-EXRB-H2	LS2-EXRB-H2	CMFIT10-EXRB-H2	CMFIT14-EXRB-H2	3D-EXR	AC-EXRB5	EXRB-H2	SRN2-EXRB-H2	
807		5.1 (WP)	0807-EXWB	Bridge	SC-EXWB	-	LL2-EXWB-H2	LS2-EXWB-H2	CMFIT10-EXWB-H2	CMFIT14-EXWB-H2	3D-EXW	AC-EXWB5	EXWB-H2	-	

NO	Company	Implant system	Platform	Library Code	Type	SC(Os) Scanbody	GMS Scanbody	LL-Link	LS-Link	CMFIT 10mm	CMFIT 14mm	3D Analog	Screw(5ea)	CCM Abutment (Co-Cr base)	Sirona Ti-base	
808			3.5 (NP)	0808-RPM	Single	GMS-RPMS	LL2-RPM-E2	CMFIT10-RPM-E2	CMFIT14-RPM-E2	3D-RPM	AC-RPM5	RPM-E2		SRN2-RPM-E2		
809	Nobel Biocare	Replace Select	4.3 (RP)	0809-RPR	Bridge	SC-RPR	GMS-RPRL	LL2-RPR-M-N2	LS2-RPR-E2	CMFIT10-RPR-E2	CMFIT14-RPR-E2	3D-RPR	RPR-E2	RPR-N2	SRN2-RPR-E2	
810			5 (WP)	0810-RPW	Bridge	SC-RPW	GMS-RPWS	LL2-RPW-M-N2	LS2-RPW-E2	CMFIT10-RPW-E2	CMFIT14-RPW-E2	3D-RPW	RPW-E2	RPW-N2	SRN2-RPW-E2	
811			6.0	0811-RP60	Bridge	SC-RP60	GMS-RP60S	LL2-RP60-H-E2	LS2-RP60-E2	CMFIT10-RP60-E2	CMFIT14-RP60-E2	3D-RP60	-	-		
901			3.0	0901-XV3	Bridge	SC-XV3	-	LL2-XV3-H2	LS2-XV3-H2	CMFIT10-XV3-H2	CMFIT14-XV3-H2	3D-XV3	AC-XV35	-		
902	Xive		3.4	0902-XVM	Bridge	SC-XVM	-	LL2-XVM-N2	LS2-XVM-H2	CMFIT10-XVM-H2	CMFIT14-XVM-H2	3D-XVM	XVM-H2	XVM-N2	SRN2-XVM-H2	
903			3.8	0903-XVR	Bridge	SC-XVR	-	LL2-XVR-H2	LS2-XVR-H2	CMFIT10-XVR-H2	CMFIT14-XVR-H2	3D-XVR	XVR-H2	XVR-N2	SRN2-XVR-H2	
904	Dentsply		4.5	0904-XW	Bridge	SC-XW	-	LL2-XW-H2	LS2-XW-H2	CMFIT10-XW-H2	CMFIT14-XW-H2	3D-XW	AC-XW5	XW-H2	SRN2-XW-H2	
905			5.5	0905-XV55	Bridge	SC-XW	-	LL2-XV55-H2	LS2-XV55-H2	CMFIT10-XV55-H2	CMFIT14-XV55-H2	3D-XV5	-	XW-N2		
906	Antylos	3.4/4.5/5.5 (Universal)	0906-AKLL	Bridge	OS-AKLL	-	LL2-AKLL-T	CMFIT10-AKLL-T	-	3D-AKLL	-	-	-	SRN2-AKLL1-T		
1001	Certain		3.4	1001-3ICM	Bridge	SC-3CM	GMS-3ICMS	LL2-3ICM-H2	LS2-3ICM-H2	CMFIT10-3ICM-H2	CMFIT14-3ICM-H2	3D-3ICM	-	3ICM-N2		
1002			4.1	1002-3ICR	Bridge	SC-3CR	GMS-3ICRS	LL2-3ICR-H2	LS2-3ICR-H2	CMFIT10-3ICR-H2	CMFIT14-3ICR-H2	3D-3ICR	-	3ICR-H2		
1003			5.0	1003-3ICW	Bridge	SC-3ICR	GMS-3ICRS	LL2-3ICW-H2	LS2-3ICW-H2	CMFIT10-3ICW-H2	CMFIT14-3ICW-H2	3D-3ICW	-	3ICRN2		
1004	Biomet 3i		6.0	1004-3IC60	Bridge	SC-3ICR	GMS-3ICRS	LL2-3IC60-H2	LS2-3IC60-H2	CMFIT10-3IC60-H2	CMFIT14-3IC60-H2	3D-3IC60	-	-		
1005			3.4	1005-EX3M	Bridge	SC-EX3M	-	LL2-EX3M-H2	LS2-EX3M-H2	CMFIT10-EX3M-H2	CMFIT14-EX3M-H2	-	-	-		
1006	External / Osseotite		4.1	1006-EXR	Bridge	SC-EXR	GMS-EXRS	LL2-EXR-H2	LS2-EXR-H2	CMFIT10-EXR-H2	CMFIT14-EXR-H2	3D-EXR	-	EXR-H2	SRN2-EXR-H2	
1007			5.0	1007-EXT	Bridge	SC-EXR	GMS-EXRS	LL2-EXT-H2	LS2-EXT-H2	CMFIT10-EXT-H2	CMFIT14-EXT-H2	-	-	EXT-H2	SRN2-EXT-H2	
1008			6.0	1008-EXTW	Bridge	SC-EXTW	GMS-EXRS	LL2-EXTW-H2	-	CMFIT10-EXTW-H2	CMFIT14-EXTW-H2	-	-	EXT-N2		
1101	Implant Direct	Legacy	3.0	1101-TSM	Single	OS-TSM	-	LL2-TSM-H2	-	CMFIT10-TSM-H2	-	3D-TSM	AC-TSM5	-	TSV3-H2	
1201	Zimmer	Tapered Screw-Vent	3.5	1201-TSV3	Bridge	SC-TSV3	GMS-TSV3S	LL2-TSV3-H2	LS2-TSV3-H2	CMFIT10-TSV3-H2	CMFIT14-TSV3-H2	3D-TSV3	-	TSV3-H2	SRN2-TSV3-H2	
1202			4.5	1202-TSV4	Bridge	SC-TSV4	GMS-TSV4S	LL2-TSV4-H2	LS2-TSV4-H2	CMFIT10-TSV4-H2	CMFIT14-TSV4-H2	3D-TSV4	-	TSV4-H2	SRN2-TSV4-H2	
1203-B 11102; Q150			5.7	1203-TSV5	Bridge	SC-TSV5	GMS-TSV5S	LL2-TSV5-H2	LS2-TSV5-H2	CMFIT10-TSV5-H2	CMFIT14-TSV5-H2	3D-TSV5	-	TSV5-H2	SRN2-TSV5-H2	

No	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	L1-Link	Ls-Link	CMFit 10mm	3D Analog	Screw(Sea)	CCM Abutment (Co-Cr base)	Sirona Ti-base
1204	Zimmer	Spline	3.7/4.0	1204-SP40	Single	SC-SP40	-	L1-2-SP40-H	-	-	-	CMFit14-SP40-H	SP40-H	-
1205			5.0	1205-SP50	Single	-	-	L1-2-SP50-H	-	-	-	CMFit14-SP50-H	-	-
1301			3.0 (Yellow)	1301-AS3	Single	SC-AS3	-	L1-2-AS3-T2	-	-	-	CMFit10-AS3-T2	-	-
1302	Osseospeed		3.5/4.0 (Aqua)	1302-ASR	Single	SC-ASR	GMS-ASRS	L1-2-ASR2-H2	L1-2-ASR20-H2	CMFit10-ASR-H2	3D-ASR	AC-ASR5	ASR-H2	SRN2-ASR-H2
1303	Astra Tech		4.5/5.0 (Lilac)	1303-ASW	Bridge	GMS-ASRL	GMS-ASWS	L1-2-ASW-H2	L1-2-ASW20-H2	CMFit14-ASW-H2	3D-ASW	AC-ASW5	ASW-H2	SRN2-ASW-H2
1304			3 (Green)	1304-AEV30	Single	SC-ASW	GMS-ASWL	L1-2-ASW-H2	L1-2-ASW20-N2	CMFit10-AEV30-T2	-	3D-AEV30	AC-AEV305	-
1305			3.6 (Purple)	1305-AEV36	Bridge	SC-AEV30	-	L1-2-AEV30-T2	-	-	-	CMFit10-AEV36-T2	AC-AEV365	-
1306			4.2 (Yellow)	1306-AEV42	Single	SC-AEV42	-	L1-2-AEV30-N2	-	-	-	CMFit10-AEV42-T2	AC-AEV425	-
1307			4.8 (Blue)	1307-AEV48	Bridge	SC-AEV42	-	L1-2-AEV36-T2	-	-	-	CMFit10-AEV48-T2	AC-AEV485	-
1308			5.4 (Brown)	1308-AEV54	Single	SC-AEV54	-	L1-2-AEV42-N2	-	-	-	CMFit10-AEV54-T2	AC-AEV545	-
1401			3.3	1401-CAM33	Bridge	SC-CAM33	-	L1-2-CAM33-T2	-	-	-	CMFit10-CAM33-T2	3D-CAM33	-
1402			3.8	1402-CAM38	Single	SC-CAM38	-	L1-2-CAM38-T2	-	-	-	CMFit10-CAM38-T2	3D-CAM38	-
1403	Camlog		4.3	1403-CAM43	Bridge	SC-CAM43	-	L1-2-CAM43-T2	-	-	-	CMFit10-CAM43-T2	3D-CAM43	-
1404	Camlog		5.0	1404-CAM50	Single	SC-CAM50	-	L1-2-CAM50-T2	-	-	-	CMFit10-CAM50-T2	3D-CAM50	-
1405			6.0	1405-CAM60	Bridge	SC-CAM560	-	L1-2-CAM60-T2	-	-	-	CMFit10-CAM60-T2	3D-CAM60	-
1406			3.3	1406-CN33	Single	GMS-CN33L	-	-	-	-	-	CMFit10-CN33-E	3D-CN33	-
1407	Coneleg		3.8/4.3	1407-CN843	Single	-	GMS-CN843L	-	-	-	-	CMFit14-CN843-E	3D-CN843	-
1408			5.0	1408-CN50	Single	-	GMS-CN50L	-	-	-	-	CMFit10-CN50-E	3D-CN50	-
1501			3.3 (Narrow)	1501-MS33	Bridge	SC-MS33	-	L1-2-MS33-H	L1-2-MS33-H	-	-	CMFit10-MS33-H	AC-MS33	-
1502	Seven	3.7/4.2 (Standard)	1502-MSR	Single	SC-TSV3	GMS-TSV3	L1-2-MSR-H	CMFit10-MSR-H	CMFit14-MSR-H	3D-TSV3	-	-	-	-
1503	MIS	5.0/5.0 (Wide)	1503-MSW	Single	SC-TSV4	GMS-TSV4S	L1-2-MSR-N	CMFit10-MSW-H	CMFit14-MSW-H	3D-TSV4	-	AC-MSS	-	-
1504		3.3	1504-C1M	Single	SC-C1M	GMS-C1MS	L1-2-C1M-T	CMFit10-C1M-T	CMFit14-C1M-T	3D-C1M	AC-C1M5	C1M-T	C1M-N	-
1505		3.75/4.2	1505-C1R	Single	SC-C1R	GMS-C1RS	L1-2-C1R-T	CMFit10-C1R-T	CMFit14-C1R-T	3D-C1R	AC-C15	C1R-T	C1R-N	-
1506	C1	5.0	1509-C1W	Single	SC-C1W	GMS-C1WS	L1-2-C1W-T	CMFit10-C1W-T	CMFit14-C1W-T	3D-C1W	-	C1W-T	C1W-N	-

No	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Link	LS-Link	CMFit 10mm	CMFit 14mm	3D Analog	Screw(5ea)	CCM Abutment (Co-Cr base)	Strona Ti-base
1601	SIC	ace / SIC max	3.3	1601-SICM	Single	-	-	LL2-SICM-H	-	CMFIT10-SICM-H	CMFIT14-SICM-H	3D-SICM	SICM-H	SICM-H	-
1602			4.2	1602-SICR	Bridge	SC-SIC	-	LL2-SICR-N	-	CMFIT10-SICR-H	CMFIT14-SICR-H	3D-SICR	SICR-H	SICR-H	-
1701	Thommen	SPI	3.5	1701-SPI35	Bridge	SC-SPI35	-	LL2-SPI35-H	LS2-SPI35-H	CMFIT10-SPI35-H	-	3D-SPI35	AC-SPI35	-	-
1702			4.0	1702-SPI40	Bridge	SC-SPI40	-	LL2-SPI40-H	LS2-SPI40-H	CMFIT10-SPI40-H	CMFIT14-SPI40-H	3D-SPI40	-	-	-
1703	Thommen	5.0	4.5	1703-SPI45	Bridge	SC-SPI45	-	LL2-SPI45-H	LS2-SPI45-H	CMFIT10-SPI45-H	CMFIT14-SPI45-H	3D-SPI45	-	-	-
1704			5.0	1704-SPI50	Bridge	SC-SPI50	-	LL2-SPI50-H	LS2-SPI50-H	CMFIT10-SPI50-H	CMFIT14-SPI50-H	3D-SPI50	AC-SPI5	-	-
1705	Bego	Bego	6.0	1705-SPI60	Bridge	SC-SPI60	-	LL2-SPI60-H	LS2-SPI60-H	CMFIT10-SPI60-H	CMFIT14-SPI60-H	-	-	-	-
1801			3.7	1801-BEGO37	Bridge	SC-BEGO	-	LL2-BEGO37-H	LS2-BEGO37-H	CMFIT10-BEGO37-H	CMFIT14-BEGO37-H	3D-BEGO37	BEGO37-H	BEGO37-N	-
1802	Bego	Bego	4.1	1802-BEGO41	Bridge	SC-BEGO	-	LL2-BEGO41-H	LS2-BEGO41-H	CMFIT10-BEGO41-H	CMFIT14-BEGO41-H	-	-	BEGO41-H	-
1803			4.5	1803-BEGO45	Single	SC-BEGO	-	LL2-BEGO45-H	-	CMFIT10-BEGO45-H	CMFIT14-BEGO45-H	-	-	BEGO41-N	-
1804	Mytis	Arrow	5.5	1804-BEGO55	Single	SC-BEGO	-	LL2-BEGO55-H	-	CMFIT10-BEGO55-H	CMFIT14-BEGO55-H	-	-	-	-
1901			3.74.0	1901-AR3740	Single	OS-AR3740	-	-	-	CMFIT10-AR3740-H	-	-	-	-	-
1902	Biohorizon	POI	4.7/6.0	1902-AR4760	Single	OS-AR4760	-	-	-	CMFIT10-AR4760-H	-	-	-	-	-
2001			3.7	2001-PQ37	Single	SC-PQ37	-	LL2-PQ37-H	-	CMFIT10-PQ37-H	-	-	-	AC-PQ375	-
2002	KYOCERA	Laser Lok	4.2	2002-PQ42	Single	SC-PQ42	-	LL2-PQ42-H	-	CMFIT10-PQ42-H	-	-	-	AC-PQ425	-
2101			3.0	2101-BH30	Single	OS-BH30	-	LL2-BH30-H	-	CMFIT10-BH30-H	-	-	-	AC-BH305	-
2201	Anatomic	AXIOM Bone Level	3.4/4.0/4.6/5.2 (Universal)	2201-ANG	Single	SC-ANG	-	LL2-ANG-E	-	CMFIT10-ANG-E	-	-	-	3D-ANG	ANG-E
2301			4.3/4.8/5.3 (Standard)	2301-WOP	Bridge	SC-WOP	-	-	-	CMFIT10-WOP-O	-	-	-	AC-WOP5	ANG-N
2401	Teeka	In-Kone	3.5/4.0/4.5/5.5 (Standard)	2401-TEK	Bridge	SC-TEK	-	-	LS2-TEK-H	CMFIT10-TEK-H	-	-	3D-TEK	AC-TEK5	TEK-H
2501			3.5	2501-KEP35	Bridge	SC-KEP35	-	LL2-KEP35-T	LS2-KEP35-T (SLKEP35)	CMFIT10-KEP35-T	CMFIT14-KEP35-T	3D-KEP35	-	-	-
2502	Keystone	PrimaComex / Max-Tri-Lobe	4.1	2502-KEP41	Bridge	SC-KEP41	-	LL2-KEP41-T	LS2-KEP41-T (SLKEP41)	CMFIT10-KEP41-T	CMFIT14-KEP41-T	3D-KEP41	-	-	-
2503			5.0/6.0	2503-KEP50	Bridge	SC-KEP50	-	LL2-KEP50-T	LS2-KEP50-T (SLKEP50)	CMFIT10-KEP50-T	CMFIT14-KEP50-T	3D-KEP50	AC-KEP5	-	-
2504	Neodent	GM	7.0	2504-KEX50	Bridge	SC-NCM	-	LL2-KEX50-T	LS2-KEX50-T (SLKEX50)	CMFIT10-KEX50-T	CMFIT14-KEX50-T	3D-KEX50	-	-	-
2505			8.0	2505-KEX65	Bridge	SC-NCM	-	LL2-KEX65-T	LS2-KEX65-T (SLKEX65)	CMFIT10-KEX65-T	CMFIT14-KEX65-T	3D-KEX65	-	-	-
2601	Neodent	GM	3.5/4.3/5.0 (Standard)	2601-NCM	Single	SC-NCG	-	LL2-NCG-H	-	CMFIT10-NCM-H	-	-	3D-NCM	-	-
2602			3.5/3.7/4.4/3.5/5.0 (Standard)	2602-NCG	Single	SC-NCG	-	LL2-NCG-H	-	CMFIT10-NCG-H	-	-	3D-NCG	AC-NCG5	-

No	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Link	LS-Link	CMFit 10mm	CMFit 14mm	3D Analog	Screw(See)	CCM Abutment (Co-Cr base)	Sirona TiBase
2701	BICON		2.0	2701-B-BC20	Bridge	SC-BC20	-	LL-2-B-BC20-N	-	-	-	-	-	-	-
2702	BICON	BICON Implant	2.5	2702-B-BC25	Bridge	SC-BC25	-	LL-2-B-BC25-N	-	-	-	-	-	-	-
2703			3.0	2703-B-BC30	Bridge	SC-BC30	-	LL-2-B-BC30-N	-	-	-	-	-	-	-
2801	Medenika	IPS	3.5/4.0/4.5/5.0	2801-MTR	Single	SC-MTR	-	-	-	-	-	-	-	AC-MTR5	-
2901	GC	Adeva Implant (Regular / Wide)	4.0/5.0	2901-GAD	Single	SC-GAD	-	-	-	-	-	-	-	AC-GAD5	-
3001	BIOTECH	KONTACT	3.0	3001-BK30	Single	SC-BK30	-	LL-2-BK30-T	LS2-BK30-T	-	-	3D-BK30	-	-	-
3002			3.6/4.2/4.8/5.4	3002-BK40	Single	SC-BK40	-	LL-2-BK40-T	LS2-BK40-T	-	-	3D-BK40	-	AC-BK405	-
3101	Medentis medical	ICX	3.75/4.14/4.8	3101-ICX	Single	SC-ICX	-	LL-2-ICX-H	LS2-ICX-H	CMFit14-ICX-H	3D-ICX	AC-ICX5	-	-	-
3201	B&B Dental	Dura-vit	4.0/3.5/5.0	3201-BNB	Single	-	GMS-BNBL	LL-2-BNB-H	-	-	-	3D-BNB	AC-BNBS	-	-
3301	Bredent Medical	Sky	3.5 / 4.0	3301-SKR	Single	-	GMS-SKRL	LL-2-SKR-T	-	CMFit14-SKR-T	3D-SKR	AC-SKR5	-	-	-
3401	Neoss	ProActive	3.5/4.0/4.5/5.0/5.5	3401-NEP	Single	-	GMS-NEPL	-	-	-	CMFit14-NEP-T	3D-NEP	AC-NEP5	-	-
3501	Champions	(Revolution	3.5/4.0/4.5/5.5	3501-CHR	Single	-	GMS-CHRL	-	-	CMFit14-CHR-H	3D-CHR	AC-CHR5	-	-	-

## Geo MultiUnit List (Scanbody / Cylinder / 3D Analog)

No	Compatible to implant Systems	Platform	Library Code	Type	SC(OS) Scanbody	LB-short	LB-long	Temporary Cylinder	3D Analog	Screw(See)	CCM Abutment (Co-Cr base)
203	Dentium Screw Abutment	4.5	0703-MUDE5	Bridge	SC-MUDE	LB-MUDE45	LB-MUDE55				AC-MUDE5
204		5.5	0704-MUDE5	Bridge	SC-MUFT	LB-MUFT					AC-MUFT5
706	Straumann Tissue Level Sincta coping	4.8	0706-MUFT	Bridge	SC-MUFT	LB-MUFT35	LB-MUFT35				AC-MUFT5
707	Straumann bonelevel RA	3.4	0707-MUSR35	Bridge	SC-MUSR35	LB-MUSR35	TA-MUSR35				AC-MUSR5
708	Nobel Multi Unit	4.6	0708-MUSR46	Bridge	SC-MUSR46	LB-MUSR46	TA-MUSR46				AC-MUSR6
812	Dentsply Nobel	4.8	0812-MUTB	Bridge	SC-MUTB	LB-MUTB2	TA-MUTB	3D-MUTB	AC-MUTB5	MUTB-H	
907	Dentsply Ankylos Multi-Unit	mini	0907-ANKLM	Bridge	SC-B-16	LB-ANKLM	LB-ANKLR				AC-ANKL5
908	Ankylos Multi-Unit regular		0908-ANKLR	Bridge	SC-MUTSV	LB-MUTSV	TA-MUTSV	3D-MUTSV	AC-MUTSV5		
1206	Zimmer Tapered Screw-Vent Tapered Coping	4.5	1206-MUTSV	Bridge	SC-MUTSV						

## Geo Transfer Scanbody List

No	Compatible to Implant Systems	Platform	Type	Section				
				Diameter GH / Height	4.0mm	4.5mm		
101	Osstem Implant TS Mini	3.0/3.5	Single	2.0 / 5.0	GTS-SUM4020	GTS-SUM4520		
				4.0 / 7.0	GTS-SUM4040	GTS-SUM4540		
				6.0 / 9.0	GTS-SUM4060	GTS-SUM4560		
102	Osstem Implant TS Standard	4.0/4.5/5.0-	Single	Diameter GH / Height	4.5mm	5.0mm	6.0mm	7.0mm
				2.0 / 5.0	GTS-SURO4520	GTS-SURO5020	GTS-SURO6020	GTS-SURO7020
				4.0 / 7.0	GTS-SURO4540	GTS-SURO5040	GTS-SURO6040	GTS-SURO7040
201	Dentium Super-Line Implantum	3.7/4.0/4.5-	Single	Diameter GH / Height	4.0mm	4.5mm	5.5mm	6.5mm
				2.0 / 5.0	GTS-DER4020	GTS-DER4520	GTS-DER5520	GTS-DER6520
				4.0 / 7.0	GTS-DER4040	GTS-DER4540	GTS-DER5540	GTS-DER6540
				6.0 / 9.0	GTS-DER4060	GTS-DER4560	GTS-DER5560	GTS-DER6560

Guide Tool	Diameter	4.0mm	4.5mm	5.0mm	5.5mm	6.0mm	6.5mm	7.0mm
	code	GTST40	GTST45	GTST50	GTST55	GTST60	GTST65	GTST70

**GEOMEDI**

Product Catalog - 2022

GeoMedi Co., Ltd.

**GEOMEDI**

PREMIUM CAD/CAM PRODUCTS

(주) 지오메디

경기도 의왕시 이미로 40 (포일동, 인덕원IT밸리) D동 702~705, 810-2, 811호

Tel : 02-588-2871 / Fax : 02-6280-2849

sales@geomedi.com www.geomedi.com

GeoMedi Co., Ltd.

#811, 810-2, 702-705, Unit D, Indeokwon IT Valley, 40, Imi-ro, Uiwang-si, Gyeonggi-do, Korea

Tel : 82-2-588-2871 Fax : 82-2-6280-2849

sales@geomedi.com www.geomedi.com



**GTS** Transfer Scanbody

**Multi-screw** Abutment