


Product Catalog

2024

Implant Solutions PREMIUM CAD/CAM PRODUCT



Implant Library System
(3Shape, exocad, Dentalwings)

LL/LS [Multibase]

CMFit [Premill, Blank]

3D-Analog

SC / GMS / GTS [Scanbody]

Sirona Multibase

CCM Casting Abutment

Geo Screw Driver

GeOMEDI

www.geomedi.com

We are exploring excellent dental implant solutions

Since 2009, we have been developing patient-customize abutment, and all of our dentists and dental technicians are also trying to make the fittest implant prosthesis for each individual person.

For Ideal Implant Restorations

We have a commitment to provide the best implant solution at a reasonable price, easy and fast way for everyone.

Co-Cr Abutment Wax up & Casting

CAD/CAM Abutment CAD design & CAM milling

For Digital Implant Solution

It is possible to make implant prosthesis by using a scanner, dental CAD, milling machine and 3D printer as developing our own digital solution for dental clinic and laboratory.



Contents

■ Company History / Certification	04	Warantec Oneplant	50
■ Dental CAD/CAM System Workflow	05	Warantec IU	50
■ Design concept	06	Tekka In-Kone	51
GMS Scanbody /SC Scanbody	06	Keystone PrimaConnex / Max-TL	52
TiN screw / Geo LS Multibase	08	Neodent CM	53
Geo LL Multibase	09	Neodent GM	53
Geo CMFit / Reverse Jig	10	Bicon Implant	54
Geo Multi-Holder for CMFit	11	Medentika IPS	55
3D-Analog	12	GC Aadvia Implant	55
Geo CAD/CAM prosthetic system	13	BIOTECH KONTACT	56
with digital library		Medentis medical ICX	57
■ Compatible Implant Systems	14	B&B Dental	58
Osstem Implant TS	14	Bredent Medical Sky	58
Osstem Implant US	15	NEOSS ProActive	59
Osstem Implant SS	16	Champions (R)evolution	59
Dentium Super-Line / Implantium	17		
Dentium NR-Line	18	■ MultiUnit(Multi-base)	60
Megagen AnyRidge	19	■ GTS Transfer Scanbody	62
Megagen AnyOne	20	■ Sirona Multibase	64
Megagen EZ-plus	21	■ CCM Casting Abutment	68
Megagen Ex-Feel [Internal]	22	■ Geo Screw Driver	73
Megagen MiNi	22	■ Instructions of CAD/CAM Abutment	76
DIO SM	23	■ Instructions of CCM Casting Abutment	77
DIO UF / UFII	23	■ Geo Library	79
NeoBiotech IS	24	■ Selection of Library	80
Dentis S-Clean / OneQ	25	■ Geo Bridge-Single Library	81
Straumann Bone Level	26	■ Products List	83
Straumann Tissue Level	27		
Straumann BLX	27		
Nobel Biocare Active	28		
Nobel Biocare Branemark	29		
Nobel Biocare Replace	30		
Dentsply Xive / Frialit2	31		
Dentsply Ankylos	32		
Biomet 3i Certain	33		
Biomet 3i Osseotite	34		
Implant Direct 3.0	35		
Zimmer Tapered Screw-vent	36		
Zimmer Spline	37		
Astra Tech Osseospeed TX	38		
Astra Tech Osseospeed EV	39		
Camlog	40		
Conelog	40		
MIS SEVEN	41		
MIS C1	42		
SIC ace / max	43		
Thommen SPI	44		
BEGO	45		
MYTIS Arrow Implant	46		
Kyocera POI	47		
Kyocera FINESIA	47		
Biohorizons Laser-Lok 3.0	48		
Anthogyr AXIOM	49		

Company History & Certification

- 2009 ● 01 Released Geo CCM Abutment
- 2010 ● 03 Established Geodent.
- 2011 ● 12 Released Geo Link, CMFit, Library, Geo Driver
- 2012 ● 10 The 1st Geo Symposium at Seoul National University
- 2013 ● 08 Obtained KGMP
08 Obtained KFDA – Implant, endosseous, laboratory supplies
10 Obtained KFDA – Geo Link abutment, Geo CCM abutment
12 Obtained KFDA – Geo CMFit abutment
12 Established GeoMedi Co.,Ltd.
- 2014 ● 07 Obtained KFDA – TCP
09 The 2nd Geo Symposium at Seoul Trade Exhibition & Convention
- 2015 ● 02 Obtained ISO13485, International Organization for Standardization' Medical devices Certification
02 Obtained CE, the European quality certification
04 Released Geo CAD/CAM system, Geo Mill
06 Certified as a venture business
07 Expansion & Relocation of corporation
09 Established research institute affiliated with a corporation
- 2016 ● 03 Obtained KFDA – Geo Scanbody, Geo 3D Analog, Geo screw
06 Obtained KFDA – Geo Multibase abutment
09 Obtained CFDA (China) – Geo Scanbody, Geo 3D analog
- 2017 ● 07 Obtained CE – Geo CCM abutment, Geo CMFit abutment *(Withdrawal, 2023.05)
07 Obtained EN ISO13485 : 2012
11 Obtained certification of INNOBIZ – Ministry of SMEs and Startups
- 2018 ● 05 Obtained CE Multibase – Geo Multibase abutment
11 Obtained EN ISO13485 : 2016 (31720001)
12 Achieved 1 Million dollars in export
- 2019 ● 01 Obtained FDA – Keystone connection abutment (Link, Premill)
09 Obtained certification of MainBIZ – Ministry of SMEs and Startups
- 2020 ● 02 Obtained KFDA – Geo GTS transfer scanbody, Geo Multiscrew abutment
10 Obtained PMDA (Japan) – Geo Multibase abutment AR(SUM/SURO/ASR/ASW)
- 2021 ● 03 Obtained PMDA (Japan) – Geo Multibase abutment IR(ITR/ITW/ISR/ISW)
04 Being selected as a Small-Giant enterprises by Ministry of Employment and Labor
04 Obtained Kyrgyzstan – Geo CCM abutment, Geo CMFit abutment, Geo Multibase abutment, Geo Multiscrew abutment, Geo GTS transfer scanbody, Geo healing cap, Geo Scanbody, Geo Impression coping, Geo 3D analog
05 Obtained KFDA – OMEGA System
05 Obtained CE – Geo GTS transfer scanbody, Geo Multiscrew abutment *(Withdrawal, 2023.05)
12 Certified as a venture business – Innovative Growth
- 2022 ● 04 Obtained RZN (Russia) – Geo GMS scanbody
04 Obtained PMDA (Japan) – Geo Multibase abutment BD(BLS/BLM/BLR)
07 Obtained Kazakhstan – Geo CCM abutment, Geo CMFit abutment, Geo Multibase abutment, Geo Multiscrew abutment, Geo healing cap, Geo Scanbody, Geo GMS scanbody, Geo Impression coping, Geo 3D analog, Geo LB Cylinder, Geo Temporary Cylinder
11 Obtained RZN (Russia) – Geo Multiscrew abutment, Geo CCM abutment, Geo CMFit abutment, Geo Multibase abutment
11 Obtained FDA – Biomet 3i OSSEOTITE² Certain², NobelActive, NobelReplace Conical Connection, NobelParallel Conical, Nobel Replace Internal Tri-Channel Connection abutment (Link, Premill)
- 2023 ● 10 Obtained Saudi Arabia – Geo Multibase abutment
- 2024 ● 05 Obtained Thailand – Geo Scanbody, Geo GMS scanbody, Geo Impression coping, Geo 3D analog, Geo screwdriver, Geo Torque wrench, Geo Screw

Dental CAD/CAM System Workflow

1. Scanbody inserting




Model and Intra-oral

2. Scanning



Geo Library System

3. Designing - 3shape / exocad



Geo Library System

Geo
Abutment & Libraries
with CAD/CAM

7. Setting



6. Sintering



5. RP and 3D printing



3D Analog

4. Milling



Multi-jig for CMFit



Titanium blank

Design concept

GMS Scanbody (Intra-oral)

- Biocompatible Titanium (Ti-6Al-4V ELI) from CARPENTER Co. guarantees high-durability.
- Autoclavable.
- Special coating allows spray-free scanning.
- Providing two different sizes - 7.5mm(Short) / 11.5mm(Long)
- Consistent performance & Visibility in X-ray.
- Various implant systems are available.
- Providing highly precise Geo Library for various CAD software. (3Shape, Exocad, Dental Wings)
- Customized abutment design for CNC & Pre-milled blank is available.
- Zirconia prosthesis designing with multi-base.



Design concept

GMS Scanbody (Intra-oral)



GMS Scanbody / SC Scanbody


SC Scanbody(Model)




Design concept

TiN Screw

- Tin-coated (Surface treatment)
- High abrasion resistance
 - Improved fixation of screw



LS Multibase



Geo Multibase & Screw made of Titanium Grade ELI
Two Screws - TiN coating screw & Titanium screw



Design concept

LL Multibase

LL Multibase height is editable on Geo Library.



Grooves - Good adhesive strength.
No need for sandblasting.

Natural curve to make emergence profile

Tin-coating - High esthetic.
LL Multibase for Posterior.



Connection

Gingival height-1.3mm --> height-5.5mm
Gingival height-2.0mm --> height-6.0mm
Two types of gingival height are available.



Zirconia - Screw retained restorations
Customized Zirconia abutment
Geo Multibase & Screw made of Titanium Grade ELI
Two Screws - TiN coating screw & Titanium screw



Design concept

Geo CMFit [Premill Blank]



Two types of diameter size are available ($\Phi 10$, $\Phi 14$)

Laser marking - Geo implant system code

Two screws - Tin coating screw & Titanium screw

Qualified material - Ti-6Al-4V ELI made by CARPENTER co.

Reverse Jig



The capability of occlusal milling


No support pin is required.

Original screws are used to fix premill blank.

* GEO milling holder & machine flexibility are required.

Design concept

Geo Multi-Holder for CMFit



● Datron



● imes-icore

Multi-holder and Reverse holder

Geo Multi-Holder for CMFit



Geo multi holder shows compatibility with various milling machines.
Available with all geo premill blanks.

● VHF Z4




● VHF N4



Design concept

Geo 3D Analog



3D Analog has grooves for additional retention.
And extra space allows material-filling to make higher stability.



Pioneer of providing 3D digital Analog
Can be used both stone model & 3D printed model
Reasonable price with excellent durability
Provided all type of implant system
Stainless Steel



Cross section
3D Analog / LL -Link / Zirconia



Geo CMFit for Titan blank



Geo Link for Zirconia

Design concept

Geo CAD/CAM prosthetic system with digital library

Dental CAD s/w



Dental CAM s/w



Geo Scanbody



Scanbody that provides the precise library matching according to the user's environment.

CMFit Abutment

Geo libraries provide stable and accurate sitting of customized abutments



Multibase Abutment




Geo Ti-base abutment is well known for its emphasis on aesthetic and functional design.

*Various gingival heights are available with dedicated library.

3D Analog

Multi-functional analog can be used for both plaster(stone) model and 3D printing model.



Retention groove on 3D analog and extra space to fill up wax provides stable fix

Bridge & Single Library



By cooperation with dental CAM s/w companies, precise direct interface milling from Geo-geometry is available.

Compatible Implant Systems

Osstem Implant TS

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				
				20Ncm ∅3.0 / 3.5 30Ncm ∅4.0 / 4.5 / 5.0

Description-Library code		0101 - SUM	0102 - SURO
Connection Platform ∅		TS Mini  ∅3.0 / 3.5	TS Standard  ∅4.0/4.5/5.0
Product			
SC-Scanbody		SC-SUM	SC-SURO
GMS-Scanbody(Short)		GMS-SUMS	GMS-SUROS
GMS-Scanbody(Long)		GMS-SUML	GMS-SUROL
LL-Multibase (Single)	GH : 1.3	LL2-SUM13-H2	LL2-SURO13-H2
	GH : 2.0	LL2-SUM20-H2	LL2-SURO20-H2
	GH : 3.0	-	LL2-SURO30-H2
	GH : 4.0	-	LL2-SURO40-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-SUM13-N2	LL2-SURO13-N2
	GH : 2.0	LL2-SUM20-N2	LL2-SURO20-N2
	GH : 3.0	-	LL2-SURO30-N2
	GH : 4.0	-	LL2-SURO40-N2
LS-Multibase / Single		LS2-SUM-H2	LS2-SURO-H2
LS-Multibase / Bridge		LS2-SUM-N2	LS2-SURO-N2
CMFit Abutment(S) / ∅10		CMFit10-SUM-H2	CMFit10-SURO-H2
CMFit Abutment(B) / ∅10		CMFit10-SUM-N2	CMFit10-SURO-N2
CMFit Abutment(S) / ∅14		CMFit14-SUM-H2	CMFit14-SURO-H2
CMFit Abutment(B) / ∅14		-	CMFit14-SURO-N2
3D-Analog		3D-SUM	3D-SURO
Screw(5ea)		AC-SUM5	AC-SURO5

Compatible Implant Systems

Osstem Implant US

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 12Hex
				30Ncm

Description-Library code	0103 - EXM	0104 - EXR	0105 - EXW
Connection Platform ∅	US Mini  ∅3.5	US Regular  ∅4.1	US Wide  ∅5.1
Product			
SC-Scanbody	SC-EXM	SC-EXR	SC-EXW
GMS-Scanbody(Short)	GMS-EXMS	GMS-EXRS	GMS-EXWS
GMS-Scanbody(Long)	GMS-EXML	GMS-EXRL	GMS-EXWL
LL-Multibase / Single	LL2-EXM-H2	LL2-EXR-H2	LL2-EXW-H2
LL-Multibase / Bridge	LL2-EXM-N2	LL2-EXR-N2	LL2-EXW-N2
LS-Multibase / Single	LS2-EXM-H2	LS2-EXR-H2	LS2-EXW-H2
LS-Multibase / Bridge	-	-	-
CMFit Abutment(S) / ∅10	CMFit10-EXM-H2	CMFit10-EXR-H2	CMFit10-EXW-H2
CMFit Abutment(B) / ∅10	-	CMFit10-EXR-N2	-
CMFit Abutment(S) / ∅14	CMFit14-EXM-H2	CMFit14-EXR-H2	CMFit14-EXW-H2
3D-Analog	3D-EXM	3D-EXR	3D-EXW
Screw(5ea)	AC-EXM5	AC-EXR5	AC-EXW5

Compatible Implant Systems

Osstem Implant SS

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				30Ncm

Description—Library code	0106 - ISR	0107 - ISW
Connection Platform ∅	SS Regular  ∅4.8	SS Wide  ∅6.0
Product		
SC-Scanbody	SC-ISRW	SC-ISRW
GMS-Scanbody(Short)	GMS-ISRS	GMS-ISWS
GMS-Scanbody(Long)	GMS-ISRL	GMS-ISWL
LL-Multibase / Single	LL2-ISR-02	LL2-ISW-02
LL-Multibase / Bridge	LL2-ISR-N2	LL2-ISW-N2
LS-Multibase / Single	-	-
CMFit Abutment / ∅10	CMFit10-ISR-02	CMFit10-ISW-02
CMFit Abutment / ∅14	CMFit14-ISR-02	CMFit14-ISW-02
3D-Analog	3D-ISR	3D-ISW
Screw(5ea)	AC-EXR5	AC-EXR5

Compatible Implant Systems

Dentium Super-Line / Implantium

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.27Hex
				25-30Ncm

Description-Library code		0201 - DER	
Product	Connection Platform ∅	Universal	 ∅3.7/4.0/4.5
	SC-Scanbody	SC-DER	
GMS-Scanbody(Short)		GMS-DERS	
GMS-Scanbody(Long)		GMS-DERL	
LL-Multibase (Single)	GH : 1.3	LL2-DER13-H2	
	GH : 2.0	LL2-DER20-H2	
	GH : 3.0	LL2-DER30-H2	
LL-Multibase (Bridge)	GH : 1.3	LL2-DER13-N2	
	GH : 2.0	LL2-DER20-N2	
	GH : 3.0	LL2-DER30-N2	
LS-Multibase / Single		LS2-DER-H2	
LS-Multibase / Bridge		LS2-DER-N2	
CMFit Abutment(S) / ∅10		CMFit10-DER-H2	
CMFit Abutment(B) / ∅10		CMFit10-DER-N2	
CMFit Abutment(S) / ∅14		CMFit14-DER-H2	
CMFit Abutment(B) / ∅14		CMFit14-DER-N2	
3D-Analog		3D-DER	
Screw(5ea)		AC-DER5	

Compatible Implant Systems

Dentium NR-Line

*Set (Scanbody/Multibase/CMFit Abutment + screw)


OS-Scanbody(Intra-Oral)		Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
		-	-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-	-		 1.0square
				20Ncm

Description—Library code	0202 - DNR	
Product	Standard	 ∅3.1/3.6/4.3
OS-Scanbody	OS-DNR	
LL-Multibase / Single	-	
LL-Multibase / Bridge	-	
LS-Multibase / Single	-	
CMFit Abutment / ∅10	CMFit10-DNR-R	
CMFit Abutment / ∅14	-	
3D-Analog	-	
Screw(5ea)	AC-DNR5	

Compatible Implant Systems

Megagen AnyRidge

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 12Hex
				35Ncm

Description-Library code		0301 - MARR
Product	Connection Platform ∅	Universal  ∅3.5/4.0/4.5
	SC-Scanbody	SC-MARR
	GMS-Scanbody(Short)	GMS-MARRS
	GMS-Scanbody(Long)	GMS-MARRL
LL-Multibase (Single)	GH : 1.3	LL2-MARR-H2
	GH : 2.0	LL2-MARR20-H2
	GH : 3.0	LL2-MARR30-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-MARR-N2
	GH : 2.0	LL2-MARR20-N2
	GH : 3.0	-
	LS-Multibase / Single	LS2-MARR-H2
	CMFit Abutment / ∅10	CMFit10-MARR-H2
	CMFit Abutment / ∅14	CMFit14-MARR-H2
	3D-Analog	3D-MARR
	Screw(5ea)	AC-MARR5

Compatible Implant Systems

Megagen AnyOne

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral) *Using 1.27Hex Driver*	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				25Ncm

Description—Library code		0302 - DERA
Product	Connection Platform ∅	Universal  ∅ 3.5/4.0/4.5/5.0
SC-Scanbody		SC-DER
GMS-Scanbody(Short)		GMS-DERS
GMS-Scanbody(Long)		GMS-DERL
LL-Multibase (Single)	GH : 1.3	LL2-DERA13-H2
	GH : 2.0	LL2-DERA20-H2
	GH : 3.0	LL2-DERA30-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-DERA13-N2
	GH : 2.0	LL2-DERA20-N2
	GH : 3.0	LL2-DERA30-N2
LS-Multibase / Single		LS2-DERA-H2
LS-Multibase / Bridge		-
CMFit Abutment(S) / ∅10		CMFit10-DERA-H2
CMFit Abutment(B) / ∅10		CMFit10-DERA-N2
CMFit Abutment(S) / ∅14		CMFit14-DERA-H2
CMFit Abutment(B) / ∅14		CMFit14-DERA-N2
3D-Analog		3D-DER
Screw(5ea)		AC-MA05

Compatible Implant Systems

Megagen EZ-plus

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-	-	-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-	-		 12Hex
				35Ncm


Megagen EZ-plus

Description-Library code	0303 - MEPM	0304 - MEPR
Connection Platform ∅	 ∅3.3	 ∅4.0/4.5/5.0
Product		
SC-Scanbody	SC-MEPM	SC-MEPR
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Multibase / Single	-	-
LL-Multibase / Bridge	-	-
LS-Multibase / Single	-	-
CMFit Abutment / ∅10	CMFit10-MEPM-E	CMFit10-MEPR-E
CMFit Abutment / ∅14	-	-
3D-Analog	-	-
Screw(5ea)	AC-MEPM5	AC-MEPR5

Compatible Implant Systems

Megagen ExFeel Internal

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				30Ncm

Description—Library code	0305 - ITRK	
Connection Platform ∅	Regular	 ∅4.8
Product		
SC-Scanbody	SC-ITRK	
GMS-Scanbody(Short)	-	
GMS-Scanbody(Long)	-	
LL-Multibase / Single	LL2-ITRK-02	
LL-Multibase / Bridge	LL2-ITRK-N2	
LS-Multibase / Single	-	
CMFit Abutment(S) / ∅10	CMFit10-ITRK-02	
CMFit Abutment(B) / ∅10	CMFit10-ITRK-N2	
CMFit Abutment(S) / ∅14	CMFit14-ITRK-02	
3D-Analog	3D-ITR	
Screw(5ea)	AC-EXR5	

Megagen MiNi

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	Single LL-Multibase	∅10 CMFit Abutment	Screw	Driver (Ncm)
	-			 1.2Hex
				15Ncm Mini 3.0/3.25

Description—Library code	0306 - MMN	
Connection Platform ∅	MINI	 ∅3.0/3.25
Product		
SC-Scanbody	SC-MMN	
LL-Multibase / Single	-	
CMFit Abutment / ∅10	CMFit10-MMN-H	
Screw(5ea)	AC-MMN5	

Compatible Implant Systems

DIO SM

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	∅10 CMFit Abutment	Single LS-Multibase	Screw	Driver (Ncm)
		-		 1.7 Torx
				30Ncm

Description-Library code	0401 - DSMM	0402 - DSMR
Connection Platform ∅	Mini  ∅3.8 / 4.1	Regular  ∅4.5 / 5.0
Product		
Scanbody	SC-DSMM	SC-DSMR
LS-Multibase / Single	-	-
CMFit Abutment ∅10	CMFit10-DSMM-T	CMFit10-DSMR-T
Screw(5ea)	AC-DSMM5	AC-DSMR5

DIO UF / UFII

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				20Ncm ∅ 3.0/3.3
				30Ncm ∅ 3.8/4.0/4.5/5.0

Description-Library code	0403 - DN30	0404 - SURD
Connection Platform ∅	UFII Narrow  ∅ 3.0/3.3	Standard  ∅ 3.8/4.0/4.5/5.0
Product		
SC-Scanbody	SC-MMN	SC-SURO
GMS-Scanbody(Short)	-	GMS-SUROS
GMS-Scanbody(Long)	-	GMS-SUROL
LL-Multibase (Single)	GH : 1.3	-
	GH : 2.0	-
	GH : 3.0	-
LL-Multibase (Bridge)	GH : 1.3	-
	GH : 2.0	-
	GH : 3.0	-
CMFit Abutment(S) / ∅10	CMFit10-DN30-H	CMFit10-SURD-H2
CMFit Abutment(B) / ∅10	-	CMFit10-SURD-N2
CMFit Abutment(S) / ∅14	-	CMFit14-SURD-H2
CMFit Abutment(B) / ∅14	-	CMFit14-SURD-N2
3D-Analog	-	3D-SURO
Screw(5ea)	AC-DN305	AC-SURD5

Compatible Implant Systems

NeoBiotech IS

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				30Ncm ∅ 3.5/4.0/4.5/5.0

Description–Library code		0501 - DERN
Connection Platform ∅	Standard	 ∅ 3.5/4.0/4.5/5.0
Product		
SC-Scanbody		SC-SURO
GMS-Scanbody(Short)		GMS-DERNS
GMS-Scanbody(Long)		GMS-DERNL
LL-Multibase (Single)	GH : 1.3	LL2-DERN-H2
	GH : 2.0	LL2-DERN20-H2
LL-Multibase (Single)	GH : 1.3	LL2-DERN-N2
	GH : 2.0	LL2-DERN20-N2
LS-Multibase / Single		-
CMFit Abutment(S) / ∅10		CMFit10-DERN-H2
CMFit Abutment(B) / ∅10		CMFit10-DERN-N2
CMFit Abutment(S) / ∅14		CMFit14-DERN-H2
CMFit Abutment(B) / ∅14		CMFit14-DERN-N2
3D-Analog		3D-DER
Screw(Sea)		AC-DERN5

Compatible Implant Systems

Dentis S-Clean / One Q

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25Hex
				30Ncm

Description–Library code	0601- DERT
Product	Regular  ∅ 4.2/4.7/5.2/6.0
SC-Scanbody	SC-DER
GMS-Scanbody(Short)	GMS-DERS
GMS-Scanbody(Long)	GMS-DERL
LL-Multibase (Single)	GH : 1.3 LL2-DERT13-H2
	GH : 2.0 -
LL-Multibase (Bridge)	GH : 1.3 LL2-DERT13-N2
	GH : 2.0 -
LS-Multibase / Single	-
CMFit Abutment(S) / ∅10	CMFit10-DERT-H2
CMFit Abutment(S) / ∅14	CMFit14-DERT-H2
3D-Analog	3D-DER
Screw(5ea)	AC-DERT5

Compatible Implant Systems

Straumann Bone Level

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 Torx
				35Ncm

Description—Library code	0700 - BLS		0701 - BLM		0702 - BLR	
Product	SC	 ∅2.9	NC	 ∅3.3	RC	 ∅4.1/4.8
SC-Scanbody	-		SC-BLM		SC-BLR	
GMS-Scanbody(Short)	-		GMS-BLMS		GMS-BLRS	
GMS-Scanbody(Long)	GMS-BLSL		GMS-BLML		GMS-BLRL	
LL-Multibase (Single)	GH : 1.3	LL2-BLS-R2	LL2-BLM-R2		LL2-BLR-R2	
	GH : 2.0	-	LL2-BLM20-R2		LL2-BLR20-R2	
	GH : 3.0	-	-		LL2-BLR30-R2	
LL-Multibase (Bridge)	GH : 1.3	-	LL2-BLM-N2		LL2-BLR-N2	
	GH : 2.0	-	LL2-BLM20-N2		LL2-BLR20-N2	
	GH : 3.0	-	-		-	
LS-Multibase / Single	-		LS2-BLM-R2		LS2-BLR-R2	
CMFit Abutment / ∅10	CMFit10-BLS-R2		CMFit10-BLM-R2		CMFit10-BLR-R2	
CMFit Abutment / ∅14	-		CMFit14-BLM-R2		CMFit14-BLR-R2	
3D-Analog	3D-BLS		3D-BLM		3D-BLR	
Screw(5ea)	AC-BLS5		AC-BLM5		AC-BLR5	

Compatible Implant Systems

Straumann Tissue Level

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 Torx
				35Ncm

Description–Library code	0703 - IT35	0704 - ITR	0705 - ITW
Connection Platform ∅	NN  ∅3.5	RN  ∅4.8	WN  ∅6.5
Product			
SC-Scanbody	SC-IT35	SC-ITRW	SC-ITRW
GMS-Scanbody(Short)	-	GMS-ITRS	GMS-ITWS
GMS-Scanbody(Long)	-	GMS-ITRL	GMS-ITWL
LL-Multibase (Single)	GH : 1.3	-	LL2-ITR-02
	GH : 2.0	-	LL2-ITR20-02
LL-Multibase (Bridge)	GH : 1.3	-	LL2-ITW-02
	GH : 2.0	-	LL2-ITW20-02
LS-Multibase / Single	LS2-IT35-0	LS2-ITR-02	LS2-ITW-02
CMFit Abutment / ∅10	CMFit10-IT35-0	CMFit10-ITR-02	CMFit10-ITW-02
CMFit Abutment / ∅14	-	CMFit14-ITW-02	CMFit14-ITW-02
3D-Analog	-	3D-ITR	3D-ITW
Screw(5ea)	AC-IT355	AC-IT5	AC-IT5

Straumann BLX

*Set (Scanbody/Multibase/CMFit Abutment + screw)

GMS-Scanbody (Intra-Oral)	Single LL-Multibase	∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw
					
(10mm)					Driver (Ncm)
					 Torx
					35Ncm
Description–Library code	0710 - BLX		0711 - BLXW		
Connection Platform ∅	 ∅ 3.75/4.5		 ∅ 5.5/6.5		
Product					
GMS-Scanbody(Long)	GMS-BLXL		GMS-BLXL		
LL-Multibase / Single	LL2-BLX-T2		-		
CMFit Abutment / ∅10	CMFit10-BLX-T2		CMFit10-BLXW-T2		
CMFit Abutment / ∅14	CMFit14-BLX-T2		CMFit14-BLXW-T2		
3D-Analog	3D-BLX		-		
Screw(5ea)	AC-BLX5		AC-BLX5		

Compatible Implant Systems

Nobel Biocare Active

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 unigrip
				15Ncm ∅3.0 35Ncm ∅3.5/4.3/5.0/6.0

Description—Library code	0801- AC30	0802- ACN	0803- ACR	0804- ACW
Connection Platform ∅	 ∅3.0	NP  ∅3.5	RP  ∅4.3/5.0	WP  ∅5.5
Product				
Scanbody	OS-AC30	SC-ACN	SC-ACR	SC-ACW
GMS-Scanbody(Short)	-	GMS-ACNS	GMS-ACRS	GMS-ACWS
GMS-Scanbody(Long)	-	GMS-ACNL	GMS-ACRL	GMS-ACWL
LL-Multibase (Single)	GH: 1.3 LL2-AC30-T2	GH: 1.3 LL2-ACN-H2	GH: 1.3 LL2-ACR-H2	GH: 1.3 LL2-ACW-H2
	GH: 2.0 -	GH: 2.0 LL2-ACN20-H2	GH: 2.0 LL2-ACR20-H2	GH: 2.0 -
LL-Multibase (Bridge)	GH: 1.3 -	GH: 1.3 LL2-ACN-N2	GH: 1.3 LL2-ACR-N2	GH: 1.3 LL2-ACW-N2
	GH: 2.0 -	GH: 2.0 LL2-ACN20-N2	GH: 2.0 LL2-ACR20-N2	GH: 2.0 -
LS-Multibase / Single	LS2-AC30-T2	LS2-ACN-H2	LS2-ACR-H2	-
CMFit Abutment / ∅10	CMFit10-AC30-T2	CMFit10-ACN-H2	CMFit10-ACR-H2	CMFit10-ACW-H2
CMFit Abutment / ∅14	-	CMFit14-ACN-H2	CMFit14-ACR-H2	CMFit14-ACW-H2
3D-Analog	3D-AC30	3D-ACN	3D-ACR	3D-ACW
Screw(Sea)	AC-AC305	AC-ACN5	AC-ACR5	AC-ACR5

Compatible Implant Systems

Nobel Biocare Branemark

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 unigrip
				35Ncm

Description-Library code	0805 - EXMB	0806 - EXRB	0807 - EXWB
Connection Platform ∅	NP ∅3.5	RP ∅4.1	WP ∅5.1
Product			
SC-Scanbody	SC-EXMB	SC-EXRB	SC-EXWB
GMS-Scanbody(Short)	-	-	-
GMS-Scanbody(Long)	-	-	-
LL-Multibase / Single	LL2-EXMB-H2	LL2-EXRB-H2	LL2-EXWB-H2
LL-Multibase / Bridge	LL2-EXMB-N2	LL2-EXRB-N2	LL2-EXWB-N2
LS-Multibase / Single	LS2-EXMB-H2	LS2-EXRB-H2	LS2-EXWB-H2
LS-Multibase / Bridge	-	-	-
CMFit Abutment(S) / ∅10	CMFit10-EXMB-H2	CMFit10-EXRB-H2	CMFit10-EXWB-H2
CMFit Abutment(S) / ∅14	CMFit14-EXMB-H2	CMFit14-EXRB-H2	CMFit14-EXWB-H2
3D-Analog	3D-EXM	3D-EXR	3D-EXW
Screw(5ea)	AC-EXMB5	AC-EXRB5	AC-EXWB5

Compatible Implant Systems

Nobel Biocare Replace

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 unigrip
				35Ncm

Description—Library code	0808- RPM	0809- RPR	0810- RPW	0811- RP60
Connection Platform ∅	NP  ∅3.5	RP  ∅4.3	WP  ∅5.0	 ∅6.0
Product				
Scanbody	SC-RPM	SC-RPR	SC-RPW	SC-RP60
GMS-Scanbody(Short)	GMS-RPMS	GMS-RPRS	GMS-RPWS	GMS-RP60S
GMS-Scanbody(Long)	GMS-RPML	GMS-RPRL	GMS-RPWL	GMS-RP60L
LL-Multibase / Single	LL2-RPM-E2	LL2-RPR-E2	LL2-RPW-E2	LL2-RP60-E2
LL-Multibase / Bridge	LL2-RPM-N2	LL2-RPR-N2	LL2-RPW-N2	LL2-RP60-N2
LS-Multibase / Single	LS2-RPM-E2	LS2-RPR-E2	LS2-RPW-E2	LS2-RP60-E2
CMFit Abutment / ∅10	CMFit10-RPM-E2	CMFit10-RPR-E2	CMFit10-RPW-E2	CMFit10-RP60-E2
CMFit Abutment / ∅14	CMFit14-RPM-E2	CMFit14-RPR-E2	CMFit14-RPW-E2	CMFit14-RP60-E2
3D-Analog	3D-RPM	3D-RPR	3D-RPW	3D-RP60
Screw(5ea)	AC-RPM5	AC-RPR5	AC-RPR5	AC-RPR5

Compatible Implant Systems

Dentsply Xive/Frialit2

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				
				24Ncm

Description—Library code	0901 - XV3	0902 - XVM	0903 - XVR	0904 - XVW	0905 - XV55
Connection Platform ∅	 ∅3.0	 ∅3.4	 ∅3.8	 ∅4.5	 ∅5.5
Product					
SC-Scanbody	SC-XV3	SC-XVM	SC-XVR	SC-XVW	SC-XVW
GMS-Scanbody(Short)	-	GMS-XVMS	GMS-XVRS	GMS-XVWS	GMS-XVWS
GMS-Scanbody(Long)	GMS-XV3L	GMS-XVML	GMS-XVRL	GMS-XVWL	GMS-XVWL
LL-Multibase / Single	LL2-XV3-H2	LL2-XVM-H2	LL2-XVR-H2	LL2-XVW-H2	LL2-XV55-H2
LL-Multibase / Bridge	LL2-XV3-N2	LL2-XVM-N2	LL2-XVR-N2	LL2-XVW-N2	LL2-XV55-N2
LS-Multibase / Single	LS2-XV3-H2	LS2-XVM-H2	LS2-XVR-H2	LS2-XVW-H2	LS2-XV55-H2
CMFit Abutment / ∅10	CMFit10-XV3-H2	CMFit10-XVM-H2	CMFit10-XVR-H2	CMFit10-XVW-H2	CMFit10-XV55-H2
CMFit Abutment / ∅14	-	CMFit14-XVM-H2	CMFit14-XVR-H2	CMFit14-XVW-H2	CMFit14-XV55-H2
3D-Analog	3D-XV3	3D-XVM	3D-XVR	3D-XVW	3D-XV55
Screw(5ea)	AC-XV35	AC-XV5	AC-XV5	AC-XV5	AC-XV5

Compatible Implant Systems

Dentsply Ankylos

*Set (Scanbody/Multibase/CMFit Abutment + screw)

OS-Scanbody(Intra-Oral)		Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				15Ncm

Description—Library code		0906 - AKLL	
Connection Platform ∅		universal	∅3.4/4.5/5.5
Product			
OS-Scanbody		OS-AKLL	
LL-Multibase (Single)	GH : 1.3	LL2-AKLL-T	
	GH : 2.0	LL2-AKLL20-T	
LL-Multibase (Bridge)	GH : 1.3	LL2-AKLL-N	
	GH : 2.0	LL2-AKLL20-N	
LS-Multibase / Single		-	
CMFit Abutment(S) / ∅10		CMFit10-AKLL-T	
CMFit Abutment(B) / ∅10		CMFit10-AKLL-N	
CMFit Abutment(S) / ∅14		CMFit14-AKLL-T	
3D-Analog		3D-AKLL	
Screw		-	

Compatible Implant Systems

Biomet 3i Certain

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				20Ncm

Description-Library code	1001 - 3ICM	1002 - 3ICR	1003 - 3ICW	1004 - 3IC60
Connection Platform ∅	 ∅3.4	 ∅4.1	 ∅5.0	 ∅6.0
Product				
SC-Scanbody	SC-3ICM	SC-3ICR	SC-3ICR	SC-3ICR
GMS-Scanbody(Short)	GMS-3ICMS	GMS-3ICRS	GMS-3ICRS	GMS-3ICRS
GMS-Scanbody(Long)	GMS-3ICML	GMS-3ICRL	GMS-3ICRL	GMS-3ICRL
LL-Multibase / Single	LL2-3ICM-H2	LL2-3ICR-H2	LL2-3ICW-H2	LL2-3IC60-H2
LL-Multibase / Bridge	LL2-3ICM-N2	LL2-3ICR-N2	LL2-3ICW-N2	LL2-3IC60-N2
LS-Multibase / Single	LS2-3ICM-H2	LS2-3ICR-H2	LS2-3ICW-H2	LS2-3IC60-H2
CMFit Abutment / ∅10	CMFit10-3ICM-H2	CMFit10-3ICR-H2	CMFit10-3ICW-H2	CMFit10-3IC60-H2
CMFit Abutment / ∅14	-	CMFit14-3ICR-H2	CMFit14-3ICW-H2	CMFit14-3IC60-H2
3D-Analog	3D-3ICM	3D-3ICR	3D-3ICW	3D-3IC60
Screw(5ea)	AC-3IC5	AC-3IC5	AC-3IC5	AC-3IC5

Compatible Implant Systems

Biomet 3i Osseotite

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				20Ncm

Description—Library code	1005- EX3M	1006- EXR	1007- EXT	1008- EXTW
Connection Platform ∅	 ∅3.4	 ∅4.1	 ∅5.0	 ∅6.0
Product				
Scanbody	SC-EX3M	SC-EXR	SC-EXR	SC-EXR
GMS-Scanbody(Short)	-	GMS-EXRS	GMS-EXRS	GMS-EXRS
GMS-Scanbody(Long)	-	GMS-EXRL	GMS-EXRL	GMS-EXRL
LL-Multibase / Single	LL2-EX3M-H2	LL2-EXR-H2	LL2-EXT-H2	LL2-EXTW-H2
LL-Multibase / Bridge	LL2-EX3M-N2	LL2-EXR-N2	LL2-EXT-N2	LL2-EXTW-N2
LS-Multibase / Single	-	LS2-EXR-H2	LS2-EXT-H2	-
LS-Multibase / Bridge	-	-	-	-
CMFit Abutment(S) / ∅10	CMFit10-EX3M-H2	CMFit10-EXR-H2	CMFit10-EXT-H2	CMFit10-EXTW-H2
CMFit Abutment(B) / ∅10	-	CMFit10-EXR-N2	-	-
CMFit Abutment(S) / ∅14	-	CMFit14-EXR-H2	CMFit14-EXT-H2	CMFit14-EXTW-H2
3D-Analog	-	3D-EXR	-	-
Screw(5ea)	AC-EXR5	AC-EXR5	AC-EXR5	AC-EXR5

Compatible Implant Systems

Implant Direct 3.0

*Set (Scanbody/Multibase/CMFit Abutment + screw)

OS-Scanbody(Intra-Oral)		Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
			-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.25Hex
				30Ncm

Implant Direct 3.0

Description—Library code	1101 - TSVM
Connection Platform ∅	Implant Direct  ∅3.0
Product	
OS-Scanbody	OS-TSVM
LL-Multibase / Single	LL2-TSVM-H2
LL-Multibase / Bridge	-
LS-Multibase / Single	-
CMFit Abutment / ∅10	CMFit10-TSVM-H2
CMFit Abutment / ∅14	-
3D-Analog	3D-TSVM
Screw(5ea)	AC-TSVM5

Compatible Implant Systems

Zimmer Tapered Screw-vent

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25Hex
				30Ncm

Description-Library code	1201- TSV3	1202- TSV4	1203- TSV5	
Connection Platform ∅	 ∅3.5	 ∅4.5	 ∅5.7	
Product				
Scanbody	SC-TSV3	SC-TSV4	SC-TSV5	
GMS-Scanbody(Short)	GMS-TSV3S	GMS-TSV4S	GMS-TSV5S	
GMS-Scanbody(Long)	GMS-TSV3L	GMS-TSV4L	GMS-TSV5L	
LL-Multibase (Single)	GH : 1.0	LL2-TSV3-H2	LL2-TSV4-H2	LL2-TSV5-H2
	GH : 2.0	LL2-TSV320-H2	LL2-TSV420-H2	-
LL-Multibase (Bridge)	GH : 1.0	LL2-TSV3-N2	LL2-TSV4-N2	LL2-TSV5-N2
	GH : 2.0	LL2-TSV320-N2	LL2-TSV420-N2	-
LS-Multibase / Single	LS2-TSV3-H2	LS2-TSV4-H2	LS2-TSV5-H2	
CMFit Abutment / ∅10	CMFit10-TSV3-H2	CMFit10-TSV4-H2	CMFit10-TSV5-H2	
CMFit Abutment / ∅14	CMFit14-TSV3-H2	CMFit14-TSV4-H2	CMFit14-TSV5-H2	
3D-Analog	3D-TSV3	3D-TSV4	3D-TSV5	
Screw(5ea)	AC-TSV5	AC-TSV5	AC-TSV5	

Compatible Implant Systems

Zimmer Spline

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-		-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
-		-		 1.25Hex
				30Ncm

Description—Library code	1204 - SP40	1205- SP50
Connection Platform ∅	 ∅ 3.7 / 4.0	 ∅ 5.0
Product		
SC-Scanbody	SC-SP40	SC-SP40
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Multibase / Single	LL2-SP40-H	LL2-SP50-H
LL-Multibase / Bridge	-	-
LS-Multibase / Single	-	-
CMFit Abutment / ∅10	-	-
CMFit Abutment / ∅14	CMFit14-SP40-H	CMFit14-SP50-H
3D-Analog	-	-
Screw(5ea)	AC-TSV5	AC-TSV5

Compatible Implant Systems

Astra Tech Osseospeed TX

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.27Hex
				15Ncm ∅3.0 20Ncm ∅3.5/4.0 25Ncm ∅4.5/5.0

Description-Library code	1301- AS3	1302- ASR	1303- ASW
Connection Platform ∅	Yellow  ∅3.0	Aqua  ∅3.5/4.0	Lilac  ∅4.5/5.0
Product			
Scanbody	SC-AS3	SC-ASR	SC-ASW
GMS-Scanbody(Short)	-	GMS-ASRS	GMS-ASWS
GMS-Scanbody(Long)	GMS-AS3L	GMS-ASRL	GMS-ASWL
LL-Multibase (Single)	GH : 1.3 LL2-AS3-T2	LL2-ASR-H2	LL2-ASW-H2
	GH : 2.0 -	LL2-ASR20-H2	LL2-ASW20-H2
LL-Multibase (Bridge)	GH : 1.3 -	LL2-ASR-N2	LL2-ASW-N2
	GH : 2.0 -	LL2-ASR20-N2	LL2-ASW20-N2
LS-Multibase / Single	-	LS2-ASR-H2	LS2-ASW-H2
CMFit Abutment / ∅10	CMFit10-AS3-T2	CMFit10-ASR-H2	CMFit10-ASW-H2
CMFit Abutment / ∅14	-	CMFit14-ASR-H2	CMFit14-ASW-H2
3D-Analog	3D-AS3	3D-ASR	3D-ASW
Screw(5ea)	AC-AS35	AC-ASR5	AC-ASW5

Compatible Implant Systems

Astra Tech Osseospeed EV

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.27Hex
				25Ncm

Description—Library code	1304 - AEV30	1305 - AEV36	1306 - AEV42	1307 - AEV48	1308 - AEV54	
Product	Green ∅3.0	Purple ∅3.6	Yellow ∅4.2	Blue ∅4.8	Brown ∅5.4	
Connection Platform ∅						
SC-Scanbody	SC-AEV30	SC-AEV36	SC-AEV42	SC-AEV48	SC-AEV54	
GMS-Scanbody(Short)	GMS-AEV30S	GMS-AEV36S	GMS-AEV42S	GMS-AEV48S	GMS-AEV54S	
GMS-Scanbody(Long)	GMS-AEV30L	GMS-AEV36L	GMS-AEV42L	GMS-AEV48L	GMS-AEV54L	
LL-Multibase (Single)	GH : 1.3	LL2-AEV30-T2	LL2-AEV36-T2	LL2-AEV42-T2	LL2-AEV48-T2	LL2-AEV54-T2
	GH : 2.0	-	LL2-AEV3620-T2	LL2-AEV4220-T2	LL2-AEV4820-T2	LL2-AEV5420-T2
LL-Multibase (Bridge)	GH : 1.3	LL2-AEV30-N2	LL2-AEV36-N2	LL2-AEV42-N2	LL2-AEV48-N2	LL2-AEV54-N2
	GH : 2.0	-	LL2-AEV3620-N2	LL2-AEV4220-N2	LL2-AEV4820-N2	LL2-AEV5420-N2
LS-Multibase / Single	-	-	-	-	-	
CMFit Abutment / ∅10	CMFit10-AEV30-T2	CMFit10-AEV36-T2	CMFit10-AEV42-T2	CMFit10-AEV48-T2	CMFit10-AEV54-T2	
CMFit Abutment / ∅14	-	CMFit14-AEV36-T2	CMFit14-AEV42-T2	CMFit14-AEV48-T2	CMFit14-AEV54-T2	
3D-Analog	3D-AEV30	3D-AEV36	3D-AEV42	3D-AEV48	3D-AEV54	
Screw(Sea)	AC-AEV305	AC-AEV365	AC-AEV425	AC-AEV48545	AC-AEV48545	

Compatible Implant Systems

Camlog

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.28Hex
				20Ncm

Description-Library code	1401- CAM33	1402- CAM38	1403- CAM43	1404- CAM50	1405- CAM60
Connection Platform ∅	 ∅3.3	 ∅3.8	 ∅4.3	 ∅5.0	 ∅6.0
Product					
Scanbody	SC-CAM33	SC-CAM38	SC-CAM43	SC-CAM5060	SC-CAM5060
LL-Multibase / Single	LL2-CAM33-T2	LL2-CAM38-T2	LL2-CAM43-T2	LL2-CAM50-T2	LL2-CAM60-T2
LL-Multibase / Bridge	LL2-CAM33-N2	LL2-CAM38-N2	LL2-CAM43-N2	LL2-CAM50-N2	LL2-CAM60-N2
LS-Multibase / Single	-	-	-	-	-
CMFit Abutment / ∅10	CMFit10-CAM33-T2	CMFit10-CAM38-T2	CMFit10-CAM43-T2	CMFit10-CAM50-T2	CMFit10-CAM60-T2
CMFit Abutment / ∅14	-	CMFit14-CAM38-T2	CMFit14-CAM43-T2	CMFit14-CAM50-T2	CMFit14-CAM60-T2
3D-Analog	3D-CAM33	3D-CAM38	3D-CAM43	3D-CAM5060	3D-CAM5060
Screw(Sea)	AC-CAM5	AC-CAM5	AC-CAM5	AC-CAM50605	AC-CAM50605

Conelog

*Set (Scanbody/Multibase/CMFit Abutment + screw)

GMS-Scanbody (Intra-Oral)	∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
 (11.5mm)					 1.28Hex
					20Ncm

Description-Library code	1406- CN33	1407- CN3843	1408- CN50
Connection Platform ∅	 ∅3.3	 ∅3.8/4.3	 ∅5.0/6.0
Product			
GMS-Scanbody(Long)	GMS-CN33L	GMS-CN3843L	GMS-CN50L
CMFit Abutment / ∅10	CMFit10-CN33-E	-	-
CMFit Abutment / ∅14	-	CMFit14-CN3843-E	CMFit14-CN50-E
3D-Analog	3D-CN33	3D-CN3843	3D-CN50
Screw(Sea)	AC-CN5	AC-CN5	AC-CN505

Compatible Implant Systems

MIS SEVEN

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25Hex
				30Ncm

Description-Library code	1501 - MS33	1502 - MSR	1503 - MSW
Connection Platform ∅	Narrow  ∅3.3	Standard  ∅ 3.75 / 4.2	Wide  ∅ 5.0 / 6.0
Product			
SC-Scanbody	SC-MS33	SC-TSV3	SC-TSV4
GMS-Scanbody(Short)	-	GMS-TSV3S	GMS-TSV4S
GMS-Scanbody(Long)	-	GMS-TSV3L	GMS-TSV4L
LL-Multibase / Single	LL2-MS33-H	LL2-MSR-H	LL2-MSW-H
LL-Multibase / Bridge	LL2-MS33-N	LL2-MSR-N	LL2-MSW-N
LS-Multibase / Single	-	-	-
CMFit Abutment / ∅10	CMFit10-MS33-H	CMFit10-MSR-H	CMFit10-MSW-H
CMFit Abutment / ∅14	-	CMFit14-MSR-H	CMFit14-MSW-H
3D-Analog	3D-MS33	3D-TSV3	3D-TSV4
Screw(5ea)	AC-MS335	AC-MS5	AC-MS5

Compatible Implant Systems

MIS C1

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25Hex
				30Ncm

Description—Library code	1504- C1M	1505- C1R	1506- C1W
Connection Platform ∅	C1  ∅3.3	C1  ∅ 3.75/4.2	C1  ∅ 5.0
Product			
Scanbody	SC-C1M	SC-C1R	SC-C1W
GMS-Scanbody(Short)	GMS-C1MS	GMS-C1RS	GMS-C1WS
GMS-Scanbody(Long)	GMS-C1ML	GMS-C1RL	GMS-C1WL
LL-Multibase / Single	LL2-C1M-T	LL2-C1R-T	LL2-C1W-T
LL-Multibase / Bridge	LL2-C1M-N	LL2-C1R-N	LL2-C1W-N
LS-Multibase / Single	-	-	-
CMFit Abutment / ∅10	CMFit10-C1M-T	CMFit10-C1R-T	CMFit10-C1W-T
CMFit Abutment / ∅14	CMFit14-C1M-T	CMFit14-C1R-T	CMFit14-C1W-T
3D-Analog	3D-C1M	3D-C1R	3D-C1W
Screw(Sea)	AC-C1M5	AC-C15	AC-C15

Compatible Implant Systems

SIC ace / max

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				20Ncm

SIC ace / max

Description—Library code	1601- SICM	1602- SICR
Connection Platform ∅	 ∅ 3.3	 ∅ 4.2
Product		
Scanbody	SC-SIC	SC-SIC
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Multibase / Single	LL2-SICM-H	LL2-SICR-H
LL-Multibase / Bridge	LL2-SICM-N	LL2-SICR-N
LS-Multibase / Single	-	-
CMFit Abutment / ∅10	CMFit10-SICM-H	CMFit10-SICR-H
CMFit Abutment / ∅14	CMFit14-SICM-H	CMFit14-SICR-H
3D-Analog	3D-SICM	3D-SICR
Screw(5ea)	AC-SIC5	AC-SIC5

Compatible Implant Systems

Thommen SPI

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase				
	-			-				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)				
				<div style="text-align: center;">  Four lobe </div> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">15Ncm</td> <td style="width: 50%;">∅3.5</td> </tr> <tr> <td>25Ncm</td> <td>∅ 4.0/4.5/5.0/6.0</td> </tr> </table>	15Ncm	∅3.5	25Ncm	∅ 4.0/4.5/5.0/6.0
15Ncm	∅3.5							
25Ncm	∅ 4.0/4.5/5.0/6.0							

Description—Library code	1701- SPI35	1702- SPI40	1703- SPI45	1704- SPI50	1705- SPI60
Connection Platform ∅	 ∅3.5	 ∅4.0	 ∅4.5	 ∅5.0	 ∅6.0
Product					
Scanbody	SC-SPI35	SC-SPI40	SC-SPI45	SC-SPI45	SC-SPI45
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-SPI35-H	LL2-SPI40-H	LL2-SPI45-H	LL2-SPI50-H	LL2-SPI60-H
LL-Multibase / Bridge	LL2-SPI35-N	LL2-SPI40-N	LL2-SPI45-N	LL2-SPI50-N	LL2-SPI60-N
LS-Multibase / Single	-	-	-	-	-
LS-Multibase / Bridge	-	-	-	-	-
CMFit Abutment / ∅10	CMFit10-SPI35-H	CMFit10-SPI40-H	CMFit10-SPI45-H	CMFit10-SPI50-H	CMFit10-SPI60-H
CMFit Abutment / ∅14	-	CMFit14-SPI40-H	CMFit14-SPI45-H	CMFit14-SPI50-H	CMFit14-SPI60-H
3D-Analog	3D-SPI35	3D-SPI40	3D-SPI45	3D-SPI50	-
Screw(5ea)	AC-SPI355	AC-SPI5	AC-SPI5	AC-SPI5	AC-SPI5

Compatible Implant Systems

BEGO

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25Hex
				30Ncm

Description—Library code	1801 - BEG037	1802 - BEG041	1803 - BEG045	1804 - BEG055
Connection Platform ∅	 ∅3.7	 ∅4.1	 ∅4.5	 ∅5.5
Product				
SC-Scanbody	SC-BEG0	SC-BEG0	SC-BEG0	SC-BEG0
GMS-Scanbody(Short)	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-
LL-Multibase / Single	LL2-BEG037-H	LL2-BEG041-H	LL2-BEG045-H	LL2-BEG055-H
LL-Multibase / Bridge	LL2-BEG037-N	LL2-BEG041-N	-	-
LS-Multibase / Single	-	-	-	-
CMFit Abutment / ∅10	CMFit10-BEG037-H	CMFit10-BEG041-H	CMFit10-BEG045-H	CMFit10-BEG055-H
CMFit Abutment / ∅14	-	-	CMFit14-BEG045-H	CMFit14-BEG055-H
3D-Analog	3D-BEG037	3D-BEG041	3D-BEG045	3D-BEG055
Screw(Sea)	AC-BEG05	AC-BEG05	AC-BEG05	AC-BEG05

Compatible Implant Systems

MYTIS Arrow Implant

*set (Scanbody/Link/CMFit Abutment + screw)

OS-Scanbody(Intra Oral)	GMS-Scanbody(Intra-Oral)	∅10 CMFit Abutment	Screw	Driver (Ncm)
	 (11.5mm)			 1.3Hex
				30Ncm

Description-Library code	1901 - AR3740	1902 - AR4760
Connection Platform ∅	 ∅ 3.7/4.0	 ∅ 4.7/6.0
Product		
OS-Scanbody	-	OS-AR4760
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	GMS-AR3740L	-
CMFit Abutment / ∅10	CMFit10-AR3740-H	CMFit10-AR4760-H
Screw(5ea)	AC-AR5	AC-AR5

Compatible Implant Systems

Kyocera POI

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	Single LL-Multibase	Bridge LL-Multibase	Ø10 CMFit Abutment	Screw	Driver (Ncm)
		-			 1.1Hex 20Ncm

Description—Library code	2001 - POI37	2002 - POI42
Connection Platform Ø	 Ø3.7	 Ø4.2
Product		
SC-Scanbody	SC-POI37	SC-POI42
LL-Multibase / Single	LL2-POI37-H	LL2-POI42-H
LL-Multibase / Bridge	-	-
CMFit Abutment / Ø10	CMFit10-POI37-H	CMFit10-POI42-H
Screw(5ea)	AC-POI375	AC-POI425

Kyocera FINESIA

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 Torx 20Ncm

Description—Library code	3601- FBLN	3602- FBLR
Connection Platform Ø	 Ø 3.2 / 3.4	 Ø 3.7 / 4.2
Product		
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	GMS-FBLNL	GMS-FBLRL
LL-Multibase / Single	-	LL2-FBLR-H
LL-Multibase / Bridge	-	-
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	CMFit10-FBLN-H	CMFit10-FBLR-H
CMFit Abutment / Ø14	-	-
3D-Analog	3D-FBLN	3D-FBLR
Screw(5ea)	AC-FBLN5	AC-FBLR5

Compatible Implant Systems

Biohorizons Laser Lok 3.0

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.25Hex
				30Ncm

Description-Library code	2101 - BH30
Connection Platform ∅	 ∅ 3.0
Product	
OS-Scanbody	OS-BH30
LL-Multibase / Single	LL2-BH30-H2
LL-Multibase / Bridge	LL2-BH30-N2
LS-Multibase / Single	-
CMFit Abutment / ∅10	CMFit10-BH30-H2
CMFit Abutment / ∅14	-
3D-Analog	3D-BH30
Screw(5ea)	AC-BH305

Compatible Implant Systems

Anthogyr AXIOM

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.18Hex
				25Ncm

Description-Library code	2201 - ANG
Product	universal  ∅ 3.4/4.0/4.6/5.2
SC-Scanbody	SC-ANG
GMS-Scanbody(Short)	GMS-ANGS
GMS-Scanbody(Long)	GMS-ANGL
LL-Multibase / Single	LL2-ANG-E
LL-Multibase / Bridge	LL2-ANG-N
LS-Multibase / Single	-
CMFit Abutment / ∅10	CMFit10-ANG-E
CMFit Abutment / ∅14	-
3D-Analog	3D-ANG
Screw(5ea)	AC-ANG5

Compatible Implant Systems

Warantec Oneplant

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	∅10 CMFit Abutment	Single LS-Multibase	Screw	Driver (Ncm)
		-		1.2Hex 30Ncm

Description-Library code	2301 - WOP	
Connection Platform ∅	Standard	 ∅ 4.3/4.8/5.3
Product		
Scanbody	SC-WOP	
LS-Multibase / Single	-	
CMFit Abutment ∅10	CMFit10-WOP-0	
Screw(5ea)	AC-WOP5	

Warantec IU

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)		-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			1.2Hex 15Ncm ∅ 3.5/4.0/4.5/5.5

Description-Library code	2302 - SUMW	
Connection Platform ∅	Mini	 ∅ 3.3/3.6
Product		
SC-Scanbody	-	
GMS-Scanbody(Short)	-	
GMS-Scanbody(Long)	GMS-SUMWL	
LL-Multibase / Single	LL2-SUMW-H2	
CMFit Abutment / ∅10	CMFit10-SUMW-H2	
CMFit Abutment / ∅14	-	
3D-Analog	3D-SUMW	
Screw(5ea)	AC-SUMW5	

Compatible Implant Systems

Tekka In-Kone

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-	-	-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.2Hex
				15Ncm ∅ 3.5/4.0/4.5/5.5

Description-Library code	2401 - TEK	
Connection Platform ∅	Standard	 ∅3.5/4.0/4.5/5.5
Product		
SC-Scanbody		SC-TEK
GMS-Scanbody(Short)		-
GMS-Scanbody(Long)		-
LS-Multibase / Single		-
CMFit Abutment / ∅10		CMFit10-TEK-H
CMFit Abutment / ∅14		-
3D-Analog		3D-TEK
Screw(5ea)		AC-TEK5

Compatible Implant Systems

Keystone PrimaConnex / Max-TL

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25 Square 30Ncm

Description—Library code	2501 - KEP35	2502 - KEP41	2503 - KEP50	2504 - KEX50	2505 - KEX60
Connection Platform ∅	 ∅ 3.5	 ∅ 4.1	 ∅ 5.0 / 6.0	 ∅ 7.0	 ∅ 8.0
Product					
SC-Scanbody	SC-KEP35	SC-KEP41	SC-KEP50	SC-KEP50	SC-KEP50
GMS-Scanbody(Short)	GMS-KEP35S	GMS-KEP41S	GMS-KEP50S	GMS-KEP50S	GMS-KEP50S
GMS-Scanbody(Long)	GMS-KEP35L	GMS-KEP41L	GMS-KEP50L	GMS-KEP50L	GMS-KEP50L
LL-Multibase / Single	LL2-KEP35-T2	LL2-KEP41-T2	LL2-KEP50-T2	LL2-KEX50-T2	LL2-KEX65-T2
LL-Multibase / Bridge	LL2-KEP35-N2	LL2-KEP41-N2	LL2-KEP50-N2	LL2-KEX50-N2	LL2-KEX65-N2
LS-Multibase / Single	LS2-KEP35-T2	LS2-KEP41-T2	LS2-KEP50-T2	LS2-KEX50-T2	LS2-KEX65-T2
CMFit Abutment / ∅10	CMFit10-KEP35-T2	CMFit10-KEP41-T2	CMFit10-KEP50-T2	-	-
CMFit Abutment / ∅14	CMFit14-KEP35-T2	CMFit14-KEP41-T2	CMFit14-KEP50-T2	CMFit14-KEX50-T2	CMFit14-KEX65-T2
3D-Analog	3D-KEP35	3D-KEP41	3D-KEP50	3D-KEX50	3D-KEX65
Screw(Sea)	AC-KEPT5	AC-KEPT5	AC-KEPT5	AC-KEPT5	AC-KEPT5

Compatible Implant Systems

Neodent CM


*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-	-	-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 0.9Hex 15Ncm

Description-Library code	2601 - NCM	
Product	Standard	 ∅ 3.5/4.3/5.0
SC-Scanbody	SC-NCM	
LS-Multibase / Single	-	
CMFit Abutment / ∅10	CMFit10-NCM-H	
CMFit Abutment / ∅14	-	
3D-Analog	3D-NCM	
Screw	-	

Neodent GM

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm)			-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 unigrip (not compatible with Nobel) 20Ncm

Description-Library code	2602 - NCG	
Product	Standard	 ∅ 3.5/3.75/4.0/4.3/5.0
SC-Scanbody	SC-NCG	
GMS-Scanbody(Long)	GMS-NCGL	
CMFit Abutment / ∅10	CMFit10-NCG-H2	
CMFit Abutment / ∅14	CMFit14-NCG-H2	
LL-Multibase (Single)	GH : 1.3	LL2-NCG-H2
	GH : 2.0	LL2-NCG20-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-NCG-N2
	GH : 2.0	LL2-NCG20-N2
3D-Analog	3D-NCG	
Screw(5ea)	AC-NCG5	

Neodent CM

Neodent GM

Compatible Implant Systems

Bicon Implant

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	Bridge LL-Multibase	Bridge LS-Multibase	∅10 CMFit Abutment	∅14 CMFit Abutment
		-		-

Description-Library code	2701- BC20	2702- BC25	2703- BC30
Connection Platform ∅	○ ∅2.0	○ ∅2.5	○ ∅3.0
Product			
SC-Scanbody	SC-BC20	SC-BC25	SC-BC30
LL-Multibase / Bridge	LL2-BC20-N	LL2-BC25-N	LL2-BC30-N
LS-Multibase / Bridge	-	-	-
CMFit Abutment / ∅10	CMFit10-BC20-N	CMFit10-BC25-N	CMFit10-BC30-N
CMFit Abutment / ∅14	-	-	-

Compatible Implant Systems

Medentika IPS

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	∅10 CMFit Abutment	∅14 CMFit Abutment	Screw	Driver (Ncm)
		-		 1.3 Hex <hr/> 25Ncm

Description–Library code	2801- MTR
Connection Platform ∅	 ∅3.5/4.3/5.0
Product	
SC-Scanbody	SC-MTR
CMFit Abutment / ∅10	CMFit10-MTR-R
Screw(5ea)	AC-MTR5

GC Aadvia Implant

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	∅10 CMFit Abutment	∅14 CMFit Abutment	Screw	Driver (Ncm)
		-		 Torx <hr/> 20Ncm

Description–Library code	2901- GAD
Connection Platform ∅	Regular / Wide  ∅4.0 / 5.0
Product	
SC-Scanbody	SC-GAD
CMFit Abutment / ∅10	CMFit10-GAD-H
Screw(5ea)	AC-GAD5

Compatible Implant Systems

BIOTECH KONTACT

BIOTECH KONTACT

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.2 Hex 20Ncm

Description–Library code	3001- BK30	3002- BK40
Connection Platform ∅	 ∅ 3.0	 ∅ 3.6 / 4.2 / 4.8 / 5.4
Product		
SC-Scanbody	SC-BK30	SC-BK40
LL-Multibase / Single	LL2-BK30-T	LL2-BK40-T
LL-Multibase / Bridge	LL2-BK30-N	LL2-BK40-N
LS-Multibase / Single	-	-
CMFit Abutment / ∅10	CMFit10-BK30-T	CMFit10-BK40-T
3D-Analog	3D-BK30	3D-BK40
Screw(5ea)	AC-BK305	AC-BK405

Compatible Implant Systems

Medentis medical ICX-System

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.42Hex
				30Ncm

Description- Library code	3101 - ICX
Connection Platform ∅	universal  ∅3.75/4.1/4.8
Product	
SC-Scanbody	SC-ICX
LL-Multibase / Single	LL2-ICX-H
LL-Multibase / Bridge	LL2-ICX-N
LS-Multibase / Single	-
CMFit Abutment / ∅10	CMFit10-ICX-H
CMFit Abutment / ∅14	CMFit14-ICX-H
3D-Analog	3D-ICX
Screw(5ea)	AC-ICX5

Compatible Implant Systems

B&B Dental

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)		-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.27 Hex
				25Ncm

Description—Library code	3201- BNB
Connection Platform ∅	 ∅4.0/4.5/5.0
Product	
GMS-Scanbody(long)	GMS-BNBL
LL-Multibase / Single	LL2-BNB-H
LL-Multibase / Bridge	-
CMFit Abutment / ∅10	CMFit10-BNB-H
3D-Analog	3D-BNB
Screw(5ea)	AC-BNB5

Bredent Medical Sky

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)		-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 Torx
				25Ncm

Description—Library code	3301- SKR
Connection Platform ∅	 ∅3.5(NP) /4.0-6.0
Product	
GMS-Scanbody(long)	GMS-SKRL
LL-Multibase / Single	LL2-SKR-T
LL-Multibase / Bridge	-
CMFit Abutment / ∅10	CMFit10-SKR-T
CMFit Abutment / ∅14	CMFit14-SKR-T
3D-Analog	3D-SKR
Screw(5ea)	AC-SKR5

Compatible Implant Systems

NEOSS ProActive

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)		-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
-				 Unigrip (not compatible with Nobel)
32Ncm				

Description–Library code	3401- NEP
Connection Platform ∅	 ∅3.5/4.0/4.5/5.0/5.5
Product	
GMS-Scanbody(long)	GMS-NEPL
LL-Multibase / Single	LL2-NEP-T
LL-Multibase / Bridge	-
CMFit Abutment / ∅10	-
CMFit Abutment / ∅14	CMFit14-NEP-T
3D-Analog	3D-NEP
Screw(5ea)	AC-NEP5

Champions (R)evolution

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)	-	-	-
∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
-				 1.2 Hex
30Ncm				

Description–Library code	3501- CHR
Connection Platform ∅	 ∅ 3.5/4.5/5.5
Product	
GMS-Scanbody(long)	GMS-CHRL
LL-Multibase / Single	-
LL-Multibase / Bridge	-
CMFit Abutment / ∅10	-
CMFit Abutment / ∅14	CMFit14-CHR-H
3D-Analog	3D-CHR
Screw(5ea)	AC-CHR5

Compatible Implant Systems

MultiUnit(Multibase)

Dentium Screw Abutment

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Ncm)
		-		 1.27Hex
				10Ncm

Description-Library code	0203- MUDE45	0204- MUDE55
Connection Platform Ø	Ø 4.5	Ø 5.5
Product		
SC-Scanbody	SC-MUDE	SC-MUDE
LB-Multibase / Bridge	LB2-MUDE45	LB2-MUDE55
3D Analog	-	-
Screw(Sea)	AC-MUDE5	AC-MUDE5

Straumann Tissue Level SynOcta coping

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Ncm)
		-		 Torx
				35Ncm

Description-Library code	0706- MUIT
Connection Platform Ø	Ø 4.8
Product	
SC-Scanbody	SC-MUIT
LB-Multibase / Bridge	LB2-MUIT
3D Analog	-
Screw(Sea)	AC-MUIT5

Straumann Bonelevel SRA

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Ncm)
				 Torx
				15Ncm

Description-Library code	0707- MUSR35	0708- MUSR46
Connection Platform Ø	Ø 3.5	Ø 4.6
Product		
SC-Scanbody	SC-MUSR35	SC-MUSR46
LB-Multibase / Bridge	LB2-MUSR35	LB2-MUSR46
3D Analog	3D-MUSR35	3D-MUSR46
Screw(Sea)	AC-MUSR5	AC-MUSR5

Compatible Implant Systems

MultiUnit (Multibase)

Nobel MultiUnit 4.8

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Ncm)
				 unigrip
				15Ncm

Description—Library code	0812- MUTB
Product	Connection Platform Ø Ø 4.8
SC-Scanbody	SC-MUTB
LB-Multibase / Bridge(4.0mm)	LB2-MUTB2
LB-Multibase / Bridge(5.5mm)	LB2-MUTBL2
3D Analog	3D-MUTB
Screw(5ea)	AC-MUTB5

Ankylos Multi-Unit

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Link	3D-Analog	Screw	Driver (Ncm)
		-		 1.0Hex
				10Ncm

Description—Library code	0907- AKLM	0908- AKLR
Product	Connection Platform Ø Mini	Regular
SC-Scanbody	SCB-16	SCB-16
LB-Multibase / Bridge	LB2-AKLM	LB2-AKLR
3D Analog	-	-
Screw(5ea)	AC-AKL5	AC-AKL5

Zimmer Tapered Screw-vent Multi-Unit(TSV)

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Ncm)
				 1.27Hex
				20Ncm



Description—Library code	1206- MUTSV
Product	Connection Platform Ø Universel
SC-Scanbody	SC-MUTSV
LB-Multibase / Bridge	LB2-MUTSV
3D Analog	3D-MUTSV
Screw(5ea)	AC-MUTSV5

GTS Transfer Scanbody

GTS Transfer Scanbody

Digital implant technology

Convenient & prompt digital implant impression
Ideal subgingiva designing referenced from digital library



Differentiate the Gingival height by appearance



- Biocompatible Titanium (Ti-6Al-4V ELI) from CARPENTER Co. guarantees high-durability.
- Autoclavable.
- Providing various sizes.
- Osstem TS And Dentium Super Line Connections are available.
- Providing High precise Geo library for CAD s/w(3Shape, Exocad, Dental wings)
- Special tool allows finding right position and preventing detachment of screw.
- Zirconia prosthesis designing with multi-base

Guide Tool



Guide Tool protects the scan body from any damage and helps to find dental implant connection easily.

Healing Abutment & Bite Index



GTS Scanbody is multifunctional that can be also used as a healing abutment and impression coping.

The groove on the scan body allows users to define different gingiva types and provide a bite index when taking an impression.

GTS Transfer scanbody

Direct Impression → Digital process



GTS scanbody in patient's mouth



Impression taking and stone model casting



Designing process with scan data from model scanner



Designing process with scan data from model scanner



Final prosthesis setting to the patient's mouth

GTS Transfer scanbody

GTS Transfer Scanbody List



Guide Tool



DER – Series (Dentium Super-line / Megagen AnyOne / Dentis S-clean / Neo IS)

Diameter / GH / Height	4.0mm	4.5mm	5.5mm	6.5mm
2.0mm / 5.0mm	GTS-DER4020	GTS-DER4520	GTS-DER5520	GTS-DER6520
4.0mm / 7.0mm	GTS-DER4040	GTS-DER4540	GTS-DER5540	GTS-DER6540
6.0mm / 9.0mm	GTS-DER4060	GTS-DER4560	GTS-DER5560	GTS-DER6560
Guide Tool	GTST-D40-127H	GTST-D45-127H	GTST-D55-127H	GTST-D65-127H

* Laser mark [DER55: Diameter 5.5mm]
Easy verifying of gingival heights by differentiated appearance.

SUM – Series (Osstem Implant TS Mini)

Diameter / GH / Height	4.0mm	4.5mm
2.0mm / 5.0mm	GTS-SUM4020	GTS-SUM4520
4.0mm / 7.0mm	GTS-SUM4040	GTS-SUM4540
6.0mm / 9.0mm	GTS-SUM4060	GTS-SUM4560
Guide Tool	GTST-S40-120H	GTST-S45-120H

* Laser mark [SUM45: Diameter 4.5mm]
Easy verifying of gingival heights by differentiated appearance.

SURO – Series (Osstem Implant TS Standard / Dio UF Regular)

Diameter / GH / Height	4.5mm	5.0mm	6.0mm	7.0mm
2.0mm / 5.0mm	GTS-SURO4520	GTS-SURO5020	GTS-SURO6020	GTS-SURO7020
4.0mm / 7.0mm	GTS-SURO4540	GTS-SURO5040	GTS-SURO6040	GTS-SURO7040
6.0mm / 9.0mm	GTS-SURO4560	GTS-SURO5060	GTS-SURO6060	GTS-SURO7060
Guide Tool	GTST-S45-120H	GTST-S50-120H	GTST-S60-120H	GTST-S70-120H

* Laser mark [SURO50: Diameter 5.0mm]
Easy verifying of gingival heights by differentiated appearance.

Sirona Ti-base

Sirona Multibase



Scan Cap for Sirona Multibase

Scan Cap

Scan Cap (short)	Scan Cap (long)
SCSN-CAP	SCSNL-CAP

Description

Osstem Implant TS

*Set (Sirona Multibase + screw)

 Single	Platform Ø	TS Mini  Ø3.0/3.5	TS Standard  Ø4.0/4.5/5.0	Driver (Ncm)  1.2Hex	
	Connection				
Single [GH:1.2mm]		SRN2-SUM-H2	SRN2-SURO-H2	20Ncm	Ø3.0/3.5
Single [GH:2.6mm]		-	SRN2-SURO26-H2	30Ncm	Ø4.0/4.5/5.0

Osstem Implant US / Biomet 3i Osseotite 4.1

*Set (Sirona Multibase + screw)

 Single	Platform Ø	US Regular  Ø4.1	Driver (Ncm)  1.2Hex
	Connection		
Single		SRN2-EXR-H2	30Ncm

Sirona Ti-base

Sirona Multibase

Description

Dentium Implantium / SuperLine

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	universal  Ø3.7/4.0/4.5		Driver (Ncm)  1.27Hex
	Connection			
	Single [GH:1.2mm]	SRN2-DER-H2		25-30Ncm
Single [GH:2.6mm]	SRN2-DER26-H2			

Megagen AnyRidge

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	Universal  Ø3.5/4.0/4.5		Driver (Ncm)  1.2Hex
	Connection			
Single	SRN2-MARR-H2		35Ncm	

Megagen AnyOne

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	Universal  Ø3.5/4.0/4.5/5.0		Driver (Ncm)  1.2Hex
	Connection			
	Single [GH:1.2mm]	SRN2-DERA-H2		25Ncm
Single [GH:2.6mm]	SRN2-DERA26-H2			

Neo Biotech IS

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	 Ø3.5/4.0/4.5		Driver (Ncm)  1.2Hex
	Connection			
	Single [GH:1.2mm]	SRN2-DERN-H2		30Ncm
Single [GH:2.6mm]	SRN2-DERN26-H2			

Straumann Bone Level

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	NC  Ø3.3	RC  Ø4.1/4.8	Driver (Ncm)  Torx
	Connection			
Single	SRN2-BLM-R2	SRN2-BLR-R2	35Ncm	

Sirona Ti-base

Sirona Multibase

Description

Straumann Tissue Level

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	RN  Ø4.8	WN  Ø6.5	Driver (Ncm)  Torx
	Connection			
	Single	SRN2-ITR-02	SRN2-ITW-02	35Ncm

Nobel Biocare Active

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	NP  Ø3.5	RP  Ø4.3/5.0	Driver (Ncm)  unigrip
	Connection			
	Single	SRN2-ACN-H2	SRN2-ACR-H2	35Ncm

Nobel Biocare Branemark

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	RP  Ø4.1	Driver (Ncm)  unigrip
	Connection		
	Single	SRN2-EXRB-H2	35Ncm

Nobel Biocare Replace

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	NP  Ø3.5	RP  Ø4.3	WP  Ø5.0	Driver (Ncm)  unigrip
	Connection				
	Single	SRN2-RPM-E2	SRN2-RPR-E2	SRN2-RPW-E2	35Ncm

Dentsply Xive / Frialit2

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	 Ø3.4	 Ø3.8	 Ø4.5	Driver (Ncm)  1.22Hex
	Connection				
	Single	SRN2-XVM-H2	SRN2-XVR-H2	SRN2-XVW-H2	24Ncm

Sirona Ti-base

Sirona Multibase

Description

Dentsply Ankylos

*Set (Sirona Multibase + screw)

 <p>Single Bridge</p>	Platform Ø	universal  Ø3.4/4.5/5.5			Driver (Ncm)  1.0Hex
	Connection				15Ncm
	Single [GH:1.0mm]	SRN2-AKLL1-T			
	Single [GH:2.0mm]	SRN2-AKLL2-T			
	Bridge [GH:1.0mm]	SRN2-AKLL1-N			
Bridge [GH:2.0mm]	SRN2-AKLL2-N				

Zimmer Tapered Screw-vent

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	 Ø3.5	 Ø4.5	 Ø5.7	Driver (Ncm)  1.25Hex
	Connection				30Ncm
Single	SRN2-TSV3-H2	SRN2-TSV4-H2	SRN2-TSV5-H2		

Astra Tech Osseospeed

*Set (Sirona Multibase + screw)

 <p>Single</p>	Platform Ø	Aqua  Ø3.5/4.0	Lilac  Ø4.5/5.0	Driver (Ncm)  1.27Hex		
	Connection				20Ncm	Ø 3.5/4.0
	Single	SRN2-ASR-H2	SRN2-ASW-H2		25Ncm	Ø 4.5/5.0

CCM Casting Abutment

Description

CCM Casting Abutment

Reliable product with accumulated clinical data since 2009
 Excellent bio-compatibility
 Optimized design for emergence profile.
 All type of implant restorations are enable.



Description


Osstem Implant TS

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø		Connection		Driver (Ncm)	
	TS Mini	Ø3.0/3.5	TS Standard	Ø4.0/4.5/5.0	1.2Hex	
Single	SUM-H2	GH : 1mm	SURO-H2		20Ncm	Ø3.0/3.5
		GH : 2.2mm	SURO22-H2			
		GH : 3mm	SURO30-H2			
Bridge	SUM-N2	GH : 1mm	SURO-N2		30Ncm	Ø4.0/4.5/5.0
		GH : 2.2mm	SURO22-N2			
		GH : 3mm	SURO30-N2			

Osstem Implant US

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø			Connection		Driver (Ncm)	
	US Mini	Ø3.5	US Regular	Ø4.1	US Wide	Ø5.1	1.2Hex
Single	EXM-H2	EXR-H2	EXW-H2	30Ncm			
Bridge	EXM-N2	EXR-N2	EXW-N2				

Osstem Implant SS

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø		Connection		Driver (Ncm)	
	SS Regular	Ø4.8	SS Wide	Ø6.0	1.2Hex	
Single	ISR-02	ISW-02	30Ncm			
Bridge	ISR-N2	ISW-N2				

CCM Casting Abutment

Description


Dentium SuperLine / Implantium *Set (CCM Abt + screw)

	Platform Ø	Standard  Ø3.7/4.0/4.5		Driver (Ncm)  1.27Hex	
	Connection				
	Single	GH : 1mm	DER-H2		25-30Ncm
		GH : 2mm	DER20-H2		
		GH : 3mm	DER30-H2		
	Bridge	GH : 1mm	DER-N2		
GH : 2mm		DER20-N2			
GH : 3mm		DER30-N2			

Megagen Anyridge *Set (CCM Abt + screw)

	Platform Ø	Standard  Ø3.5/4.0/4.5		Driver (Ncm)  1.2Hex	
	Connection				
	Single	MARR-H2		35Ncm	
Bridge	MARR-N2				

Megagen Ex-feel [Internal] *Set (CCM Abt + screw)

	Platform Ø	 Ø 4.8		Driver (Ncm)  1.2Hex	
	Connection				
	Single	ITRK-02		35Ncm	
Bridge	ITRK-N2				



- Qualified Raw Material
Cobalt-chrome (CoCr) alloys
- Single and bridge case
Available single and bridge types of abutments
- Two screws
TiN coating screw & Titanium screw
- Compatibility
Various compatible Implant systems

CCM Casting Abutment

CCM Casting Abutment

CCM Casting Abutment

Description

DIO UF / UFII

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	 Ø3.8/4.0/4.5/5.0		Driver (Ncm)	1.2Hex
	Connection				
	Single	SURD-H2		30Ncm	
Bridge	SURD-N2				

Neo Biotech IS

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	 Ø3.5/4.0/4.5		Driver (Ncm)	1.2Hex
	Connection				
	Single	DERN-H2		30Ncm	
Bridge	DERN-N2				

Straumann Bone Level

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	NC	 Ø3.3	RC	 Ø4.1/4.8	Driver (Ncm)	Torx
	Connection						
	Single	BLM-R2	BLR-R2		35Ncm		
Bridge	BLM-N2	BLR-N2					

Straumann Tissue Level

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	RN	 Ø4.8	WN	 Ø6.5	Driver (Ncm)	Torx
	Connection						
	Single	ITR-O2	ITW-O2		35Ncm		
Bridge	ITR-N2	ITW-N2					

Nobel Biocare Active

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	NP	 Ø3.5	RP	 Ø4.3/5.0	Driver (Ncm)	unigrip
	Connection						
	Single	ACN-H2	ACR-H2		35Ncm		
Bridge	ACN-N2	ACR-N2					

Nobel Biocare Branemark

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	NP	 Ø3.5	RP	 Ø4.1	WP	 Ø5.1	Driver (Ncm)	unigrip
	Connection								
	Single	EXMB-H2	EXRB-H2	EXWB-H2		35Ncm			
Bridge	EXMB-N2	EXRB-N2	EXWB-N2						

CCM Casting Abutment


Description


Nobel Biocare Replace						*Set (CCM Abt + screw)
 Single Bridge	Platform Ø	NP  Ø3.5	RP  Ø4.3	WP  Ø5.0	Driver (Ncm)	 unigrip
	Connection					
	Single	RPM-E2	RPR-E2	RPW-E2		35Ncm
Bridge	RPM-N2	RPR-N2	RPW-N2			


Nobel Biocare Multi Unit						*Set (CCM Abt + screw)
 Bridge	Platform Ø	 Ø4.8			Driver (Ncm)	 unigrip
	Connection	Bridge			MUTB-N2	15Ncm

Dentsply Xive / Frialit2						*Set (CCM Abt + screw)
 Single Bridge	Platform Ø	 Ø3.4	 Ø3.8	 Ø4.5	Driver (Ncm)	 1.22Hex
	Connection					
	Single	XVM-H2	XVR-H2	XVW-H2		24Ncm
Bridge	XVM-N2	XVR-N2	XVW-N2			

Dentsply Ankylos						*Set (CCM Abt + screw)
 Bridge	Platform Ø	A/B/C/D  Ø3.4/4.5/5.5			Driver (Ncm)	 1.0Hex
	Connection	Single			-	25Ncm
	Bridge	DNF3-N				

Biomet 3i Certain						*Set (CCM Abt + screw)
 Single Bridge	Platform Ø	 Ø3.4	 Ø4.1	Driver (Ncm)	 1.2Hex	
	Connection					
	Single	3ICM-H2	3ICR-H2		20Ncm	
Bridge	3ICM-N2	3ICR-N2				

Biomet 3i Osseotite						*Set (CCM Abt + screw)
 Single Bridge	Platform Ø	 Ø4.1	 Ø5.0	Driver (Ncm)	 1.2Hex	
	Connection					
	Single	EXR-H2	EXT-H2		20Ncm	
Bridge	EXR-N2	EXT-N2				

Zimmer Tapered Screw-vent						*Set (CCM Abt + screw)
 Single Bridge	Platform Ø	 Ø3.5	 Ø4.5	 Ø5.7	Driver (Ncm)	 1.25Hex
	Connection					
	Single	TSV3-H2	TSV4-H2	TSV5-H2		30Ncm
Bridge	TSV3-N2	TSV4-N2	TSV5-N2			

CCM Casting Abutment

Description

Zimmer Spline

*Set (CCM Abt + screw)

 Single	Platform Ø	 Ø 3.7 / 4.0		Driver (Ncm)	 1.25Hex
	Connection	Single	SP40-H		30Ncm

Astra Tech Osseospeed

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	X-Small  Ø3.0	Aqua  Ø3.5/4.0	Lilac  Ø4.5/5.0	Driver (Ncm)	 1.27Hex	
	Connection	Single	AS3-T2	ASR-H2	ASW-H2	15Ncm	Ø3.0
	Bridge	-	ASR-N2	ASW-N2	20Ncm	Ø3.5/4.0	
					25Ncm	Ø4.5/5.0	

MIS C1

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	C1  Ø3.3	C1  Ø 3.75/4.2	C1  Ø 5.0	Driver (Ncm)	 1.25Hex
	Connection	Single	C1M-T	C1R-T	C1W-T	30Ncm
	Bridge	C1M-N	C1R-N	C1W-N		

SIC ace / max

*Set (CCM Abt + screw)

 Single	Platform Ø	 Ø3.3	 Ø4.2	Driver (Ncm)	 1.2Hex
	Connection	Single	SICM-H	SICR-H	20Ncm

Anthogyr AXIOM

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	universal  Ø 3.4/4.0/4.6/5.2		Driver (Ncm)	 1.18Hex
	Connection	Single	ANG-E		25Ncm
	Bridge	ANG-N			

Tekka In-Kone

*Set (CCM Abt + screw)

 Single	Platform Ø	Standard  Ø 3.5/4.0/4.5/5.5		Driver (Ncm)	 1.2Hex
	Connection	Single	TEK-H		15Ncm

Geo Screw Driver



Angled Driver

Square Type / Ratchet Type
 Use with LS Multibase Abutment
 Maximum Angle 25°
 Recommended torque allowance : below 30Ncm.(reference manufacture's guide)
 Caution: unstable connection may bend or destroy the driver
 Recommended usage (1.2hex/1.27hex)-10 times
 (Universal Type)- 20 times
 Connect with Square type Torque wrench
 May use Round Type Torque wrench
 by connecting HDEC - Short (Clinic)Handle



Connect Straight & Angle driver with the Torque Wrench



Straight Driver

Square Type / Ratchet Type
 Connect by manufacture's recommended torque allowance
 Damage on the tip or usage of long period of time need to exchange driver
 Connect with Square type Torque wrench
 May use Round Type Torque wrench
 by connecting HDEC - Short (Clinic)Handle

Universal - Round hexagon

Square Type / Ratchet Type
 Connect by manufacture's recommended torque allowance
 Damage on the tip or usage of long period of time need to exchange driver
 Connect with Square type Torque wrench
 May use Round Type Torque wrench
 by connecting HDEC - Short (Clinic)Handle

Angled Driver

3i & Osstem - 1.2mm Hex

			
ASD3IL Square -20mm	ASD3IR Square -13mm	ALD3IL Ratchet -20mm	ALD3IR Ratchet -13mm

Astra & Zimmer & Dentium - 1.27mm Hex

			
ASDZIL Square -20mm	ASDZIR Square -13mm	ALDZIL Ratchet -20mm	ALDZIR Ratchet -13mm

Universal - Round hexagon

			
ASDUVL Square -20mm	ASDUVR Square -13mm	ALDUVL Ratchet -20mm	ALDUVR Ratchet -13mm

Geo Screw Driver

Straight Driver

3i - 1.2mm Hex

 SSD3IL Square -20mm	 SSD3IR Square -13mm	 SLD3IL Ratchet -20mm	 SLD3IR Ratchet -13mm
 SSD3IEL Square -35mm			

Zimmer - 1.27mm Hex

 SSDZIL Square -20mm	 SSDZIR Square -13mm	 SLDZIL Ratchet -20mm	 SLDZIR Ratchet -13mm
--	--	--	---

Astra - 1.27mm Hex

 SSDATL Square -20mm	 SSDATR Square -13mm	 SLDATL Ratchet -20mm	 SLDATR Ratchet -13mm
 SSDATEL Square -35mm			

Nobel - unigrip

 SSDNBL Square -20mm	 SSDNBR Square -13mm	 SLDNBL Ratchet -20mm	 SLDNBR Ratchet -13mm
 SSDNBEL Square -35mm			

Straumann - Torx

 SSDSTL Square -20mm	 SSDSTR Square -13mm	 SLDSTL Ratchet -20mm	 SLDSTR Ratchet -13mm
 SSDSTEL Square -35mm			

Geo Screw Driver

Driver Tool

 <p>HDEL Long (Lab) Handle</p>	 <p>Long (Lab) Handle & Driver</p>
 <p>HDEC Short (Clinic) Handle</p>	 <p>Short (Clinic) Handle & Driver</p>
 <p>GTW Torque Wrench</p>	 <p>Round Type Torque Wrench & Handle Short & Driver</p>  <p>Square Type Torque Wrench & Driver</p>

Lab Driver Set



Description		
Product		
Size	Code	Implant System
 1.2 Hex	Blue	Osstem Implant / Dio Implant Biomet 3i / Xive
 1.27 Hex	Red	Dentium / Astra Tech Zimmer
 Unigrip	Yellow	Nobel Biocare
 Torx	Green	Straumann
 1.0 Hex	Purple	Ankylos
 Universal	Black	Dynamic

Instructions of CAD/CAM Abutment

< Procedure for intra-oral scan >



1. Scanbody inserting



2. Scanning



3. Matching [Scanbody and Library]



4. Designing - 3shape / Exocad



5. Sintering



6. Setting

< Procedure for model scan >



1. Scanbody inserting



2. Scanning



3. Matching [Scanbody and Library]



4. Designing - 3shape / Exocad



5. Sintering and model setting

Instructions of CCM Casting Abutment

< Screw Type manufacturing >



CCM Casting abutment



Status of Wax-up the connection of the metal part has more wax to prevent the occurrence of air bubble in the connection part.



Status after casting of general Porcelain metal.



Investment material and the oxide film is removed by using a Glass bead (50um/ 2-4 bar)



The oxide film can be removed cleanly without damaging the internal side of Hex.



The Junction is checked in the shape of work in the condition where the oxide film is removed.



Upon sanding with alumina oxide, it can be the cause of plaque.



Status after the oxide film is removed.



The formation of oxide film can be observed on the basal area of firing porcelain.



The firing porcelain part is wrapped by using wax.



The oxide film is removed by using the Glass bead.



The oxide film of internal side of Access hole is removed.



Final compatibility of finished upper structure is checked.



< Customized Type manufacturing >



Full contour wax up



Gum sculpture



after Casting



Instructions of CCM Casting Abutment

< Hybrid Type manufacturing >



Wax up on the Model



Separation of Abutment from Model



Sprue installation



After Casting



Screw Type is manufactured by using Denture Resin



< Bar Overdenture manufacturing >



Wax up on the Model




Milling in the desired form on the Wax



Bar Milling after Casting



Check the compliance



- Free Geo Library available for 3Shape and exocad
- High quality of stability and accuracy
- One scanbody can be used to design all type of implant restorations
- Providing library for Ti-base
- Premill blank and implant bridge



Benefits :

- Optimized to show best function in dental CAD
- Enable to design individual prosthetic solutions for various implant systems.
- Available in 3Shape and Exocad with various Ti-bases.
- Flexible to fabricate milled structures with its individuality.



Application :

- Individual ceramic abutments on titanium bases
- Individual one-piece titanium abutments
- Occlusal screwable crowns and bridges
- Bars, telescopes
- Screw retained abutment
- Customized Zirconia abutment
- Customized titanium abutment





Selection of Geo Library

Library for Implant System

- Geo_0101_SUM
- Geo_0102_SURO
- Geo_0103_EXM
- Geo_0104_EXR

Geo product code for different type





Cement [gap] option with correspond Scanbody

Cement [gap]-Zirconia

- 027 SC - Model Scanbody
- 035 SC - Model Scanbody
- 042 SC - Model Scanbody
- 052 SC - Model Scanbody
- 027 GMS-S - Intra-oral Scanbody
- 035 GMS-S - Intra-oral Scanbody
- 045 GMS-S - Intra-oral Scanbody
- 052 GMS-S - Intra-oral Scanbody
- 027 GMS-L - Intra-oral Scanbody
- 035 GMS-L - Intra-oral Scanbody
- 042 GMS-L - Intra-oral Scanbody
- 052 GMS-L - Intra-oral Scanbody
- 08 PMMA-SC - Model Scanbody
- 10 PMMA-SC - Model Scanbody
- 08 PMMA-GMS-S - Intra-oral Scanbody
- 10 PMMA-GMS-S - Intra-oral Scanbody
- 08 PMMA-GMS-L - Intra-oral Scanbody
- 10 PMMA-GMS-L - Intra-oral Scanbody

Cement [gap]-PMMA,Metal



Geo BS Library- Bridge&Single



Geo Bridge/Single Library

Brand	Implant System	System code	Bridge(Bar)	Single	Scanbody	Screw	
					code	Geo screw code	Compatible with original
Osstem	TSM	SUM	✓	✓	SC-SUM	SUMC	v
	TS Standard	SURO	✓	✓	SC-SURO	SUROC	v
	US M	EXM	✓	✓	SC-EXM	EXMC	v
	USR	EXR	✓	✓	SC-EXR	EXRC	v
	US W	EXW	✓	✓	SC-EXW	EXWC	v
Dentium	Superline	DER	✓	✓	SC-DER	DERC	v
	Screw-Abutment	MUDE45	✓		SC-MUDE	MUDEC	v
		MUDE55	✓		SC-MUDE	MUDEC	v
Straumann	Bonelevel 3.3	BLM	✓	✓	SC-BLM	BLMC	v
	Bonelevel 4.1	BLR	✓	✓	SC-BLR	BLMC	v
	Bonelevel 4.8	BLR	✓	✓	SC-BLR	BLRC	v
	Tissue level 4.8	ITR	✓	✓	SC-ITRW	bridge - ITTC	Temporary screw
						single - ITC	v
	Tissue level 6.5	ITW	✓	✓	SC-ITRW	bridge - ITTC	Temporary screw
						single - ITC	v
	SynOcta coping	MUIT	✓		SC-MUIT	MUITC	v
SRA Abutment	MUSR35	✓		SC-MUSR35	MUSRC	v	
	MUSR46	✓		SC-MUSR46	MUSRC	v	
Nobel	Multi Unit	MUTB	✓		SC-MUTB	MUTBC	v
	Active 3.5	ACN	✓	✓	SC-ACN	ACNC	v
	Active 4.3	ACR	✓	✓	SC-ACR	ACRC	v
	Active 5.0	ACR	✓	✓	SC-ACR	ACRC	v
	Active 5.5	ACW	✓	✓	SC-ACW	ACRC	v
	Branemark 3.4	EXM	✓	✓	SC-EXM(B)	EXMBC	v
	Branemark 4.1	EXR	✓	✓	SC-EXR(B)	EXRBC	v
	Branemark 5.1	EXW	✓	✓	SC-EXW(B)	EXWBC	v
	Replace 3.5 rpm	RPM	✓	✓	SC-RPM	RPMC	v
	Replace 4.3 rpm	RPR	✓	✓	SC-RPR	RPRC	v
	Replace 5.0 rpm	RPW	✓	✓	SC-RPW	RPRC	v
Replace 6.0 rpm	RP60	✓	✓	SC-RP60	RPRC	v	
Dentsply	Xive 3.0	XV3	✓	✓	SC-XV3	XV3C	v
	Xive 3.4	XVM	✓	✓	SC-XVM	XVC	v
	Xive 3.8	XVR	✓	✓	SC-XVR	XVC	v
	Xive 4.5	XVW	✓	✓	SC-XVW	XVC	v
	Xive 5.5	XV55	✓	✓	SC-XVW	XVC	v
Ankylos	AK II	AKLL	✓		OS-AKLL	DNF3C	special screw

Geo BS Library- Bridge&Single

Brand	Implant System	System code	Bridge(Bar)	Single	Scanbody	Screw	
					code	Geo screw code	Compatible with original
3i	Certain 3.4	3ICM	✓	✓	SC-3ICM	3ICC	v
	Certain 4.1	3ICR	✓	✓	SC-3ICR	3ICC	v
	Certain 5.0	3ICW	✓	✓	SC-3ICR	3ICC	v
	Certain 6.0	3IC60	✓	✓	SC-3ICR	3ICC	v
	External 3.4	EX3M	✓	✓	SC-EX3M	EXRC	v
	External 4.1	EXR	✓	✓	SC-EXR	EXRC	v
	External 5.0	EXT	✓	✓	SC-EXR	EXRC	v
	External 6.0	EXTW	✓	✓	SC-EXR	EXRC	v
Zimmer	TSV 3.5	TSV3	✓	✓	SC-TSV3	TSVC	v
	TSV 4.5	TSV4	✓	✓	SC-TSV4	TSVC	v
	TSV 5.7	TSV4	✓	✓	SC-TSV5	TSVC	v
Astra	Astra EV 3.0	AEV30	✓		SC-AEV30	AEV30C	v
	Astra EV 3.6	AEV36	✓		SC-AEV36	AEV36C	v
	Astra EV 4.2	AEV42	✓		SC-AEV42	AEV42C	v
	Astra EV 4.8	AEV48	✓		SC-AEV48	AEV4854C	v
	Astra EV 5.4	AEV54	✓		SC-AEV54	AEV4854C	v
	Astra 3.5	ASR	✓	✓	SC-ASR	ASRC	v
	Astra 4.0	ASR	✓	✓	SC-ASR	ASRC	v
	Astra 4.5	ASW	✓	✓	SC-ASW	ASWC	v
	Astra 5.0	ASW	✓	✓	SC-ASW	ASWC	v
Camlog	Camlog 3.3	CAM33	✓		SC-CAM33	CAMC	v
	Camlog 3.8	CAM38	✓		SC-CAM38	CAMC	v
	Camlog 4.3	CAM43	✓		SC-CAM43	CAMC	v
	Camlog 5.0	CAM50	✓		SC-CAM5060	CAM5060C	v
	Camlog 6.0	CAM60	✓		SC-CAM5060	CAM5060C	v
Magazen	AnyOne	MAO	✓	✓	SC-DER	MAOC	v
	Anyridge	MARR	✓	✓	SC-MARR	MARRC	v
Dentis	S-clean	DET	✓	✓	SC-DER	DERC	v
Dio	UF	SURD	✓	✓	SC-SURO	SURDC	v
Neo	IS	DERN	✓	✓	SC-SURO	DERNC	v
MIS	C1 3.3	C1M	✓		SC-C1M	C1MC	v
	C1 3.75 / 4.2	C1R	✓		SC-C1R	C1C	v
	C1 5.0	C1W	✓		SC-C1W	C1C	v

Geo Product List (Scanbody / Multibase / CMFit / 3D Analog / CCM Abutment / Sirona Ti-base)

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMFit 10mm	CMFit 14mm	3D Analog	Screw (Sea Pack)	CCM Abutment (Co-Cr base)	Sirona Ti-base
101			3.0/3.5	0101-SUM	Single	SC-SUM	GMS-SUMS	LL2-SUM13-H2 LL2-SUM20-H2	LS2-SUM-H2	CMFT10-SUM-H2	CMFT14-SUM-H2	3D-SUM	AC-SUM5	SUM-H2	SRN2-SUM-H2
					Bridge		GMS-SUML	LL2-SUM13-N2 LL2-SUM20-N2	LS2-SUM-N2	CMFT10-SUM-N2					
102	Osstem	TS	4.0/4.5/5.0	0102-SURO	Single	SC-SURO	GMS-SUROS	LL2-SURO13-H2 LL2-SURO20-H2 LL2-SURO30-H2 LL2-SURO40-H2	LS2-SURO-H2	CMFT10-SURO-H2	CMFT14-SURO-H2	3D-SURO	AC-SURO5	SURO-H2 SURO22-H2 SURO30-H2	SRN2-SURO-H2 SRN2-SURO26-H2
					Bridge		GMS-SUROL	LL2-SURO13-N2 LL2-SURO20-N2 LL2-SURO30-N2 LL2-SURO40-N2	LS2-SURO-N2	CMFT10-SURO-N2	CMFT14-SURO-N2				
103			3.5	0103-EXM	Single Bridge	SC-EXM	GMS-EXMS GMS-EXML	LL2-EXMH2 LL2-EXMN2	LS2-EXMH2	CMFT10-EXMH2	CMFT14-EXMH2	3D-EXM	AC-EXM5	EXM-H2 EXM-N2	-
104		US	4.1	0104-EXR	Single	SC-EXR	GMS-EXRS GMS-EXRL	LL2-EXRH2 LL2-EXRN2	LS2-EXRH2	CMFT10-EXRH2 CMFT10-EXRN2	CMFT14-EXRH2	3D-EXR	AC-EXR5	EXR-H2 EXR-N2	SRN2-EXR-H2
					Bridge		GMS-EXWS GMS-EXWL	LL2-EXWH2 LL2-EXWN2	LS2-EXWH2	CMFT10-EXWH2	CMFT14-EXWH2				
106		SS	4.8	0106-ISR	Single	SC-ISRW	GMS-ISRS GMS-ISRL	LL2-ISR-O2 LL2-ISR-N2	-	CMFT10-ISR-O2	CMFT14-ISR-O2	3D-ISR	AC-EXR5	ISR-O2 ISR-N2	-
					Bridge		GMS-ISWS GMS-ISWL	LL2-ISW-O2 LL2-ISW-N2	-	CMFT10-ISW-O2	CMFT14-ISW-O2				
107			6.0	0107-ISW	Single	SC-ISRW	GMS-ISRS GMS-ISRL	LL2-ISR-O2 LL2-ISR-N2	-	CMFT10-ISR-O2	CMFT14-ISR-O2	3D-ISR	AC-EXR5	ISR-O2 ISR-N2	-
					Bridge		GMS-ISWS GMS-ISWL	LL2-ISW-O2 LL2-ISW-N2	-	CMFT10-ISW-O2	CMFT14-ISW-O2				
201	Dentium	Super-Line	3.7/4.0/4.5 (Universal)	0201-DER	Single	SC-DER	GMS-DERS	LL2-DER13-H2 LL2-DER20-H2 LL2-DER30-H2	LS2-DER-H2	CMFT10-DER-H2	CMFT14-DER-H2	3D-DER	AC-EXR5	DER-H2 DER20-H2 DER30-H2	SRN2-DER-H2 SRN2-DER26-H2
					Bridge		GMS-DERL	LL2-DER13-N2 LL2-DER20-N2 LL2-DER30-N2	LS2-DER-N2	CMFT10-DER-N2	CMFT14-DER-N2				
202			3.1/3.6/4.3 (Standard)	0202-DNR	Single	OS-DNR	-	-	-	CMFIT10-DNR-R	-	-	AC-DNR5	-	-
301		AnyRidge	3.5/4.0/4.5 (Universal)	0301-MARR	Single	SC-MARR	GMS-MARRS	LL2-MARR-H2 LL2-MARR20-H2 LL2-MARR30-H2	LS2-MARR-H2	CMFT10-MARR-H2	CMFT14-MARR-H2	3D-MARR	AC-MARR5	MARR-H2	SRN2-MARR-H2
					Bridge		GMS-MARRL	LL2-MARR-N2 LL2-MARR20-N2 LL2-MARR30-N2	-	-					
302	Megagen	AnyOne	3.5/4.0/4.5/5.0 (Universal)	0302-DERA	Single	SC-DER	GMS-DERS	LL2-DEA13-H2 LL2-DEA20-H2 LL2-DEA30-H2	LS2-DEA-H2	CMFT10-DEA-H2	CMFT14-DEA-H2	3D-DER	AC-MAO5	-	SRN2-DEA-H2
					Bridge		GMS-DERL	LL2-DEA13-N2 LL2-DEA20-N2 LL2-DEA30-N2	-	CMFT10-DEA-N2	CMFT14-DEA-N2				
303		EZ-Plus	3.3	0303-MEPM	Single	SC-MEPM	-	-	CMFT10-MEPM-E	-	-	AC-MEPM5	-	-	-
304			4.0/4.5/5.0	0304-MEPR	Single	SC-MEPR	-	-	CMFT10-MEPR-E	-	-	AC-MEPR5	-	-	-
305		EX-feel	4.8 (Regular)	0305-ITRK	Single	SC-ITRK	-	LL2-ITRK-O2	-	CMFT10-ITRK-O2	CMFT14-ITRK-O2	3D-ITR	AC-EXR5	ITRK-O2	-
					Bridge		GMS-ITRKN	LL2-ITRKN2	-	CMFT10-ITRKN2					
306		Miri	3.0/3.25	0306-MMN	Single	SC-MMN	-	-	CMFT10-MMNH	-	-	AC-MMN5	-	-	
401	Dio	SM	3.8/4.1	0401-DSMM	Single	SC-DSMM	-	-	-	CMFT10-DSMM-T	-	-	AC-DSMM5	-	-
					Bridge		GMS-DSMRT	LL2-DSMRT	-	CMFT10-DSMRT-T					
402			4.5/5.0	0402-DSMR	Single	SC-DSMR	-	-	CMFT10-DSMR-T	-	-	AC-DSMR5	-	-	

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMFk 10mm	CMFk 14mm	3D Analog	Screw (See Pack)	CCM Abutment (Co-Cr base)	Sirona T-base
403	Dio	UF / UFI	3.0/3.3 (Narrow)	0403-DNG0	Single	SC-MMN	-	-	-	CMFIT10-DN30-H	-	-	AC-DN305	-	-
404					Bridge	SC-SURO	GMS-SUROS	LL2-SURD-H2 LL2-SURD20-H2 LL2-SURD30-H2 LL2-SURD-N2 LL2-SURD20-N2 LL2-SURD30-N2	-	-	-	CMFIT10-SURD-H2 CMFIT10-SURD-N2	CMFIT14-SURD-H2 CMFIT14-SURD-N2	3D-SURO	AC-SURD5
501	NeoBiotech	IS	3.5/4.0/4.5/5.0-	0501-DERN											
601					Dentis	S-Clean/OneQ	3.9/4.2/4.7/5.2	0601-DERT	Single	SC-DER	GMS-DERS	LL2-DERT13-H2 LL2-DERT13-N2	CMFIT10-DERT-H2	CMFIT14-DERT-H2	3D-DER
700	Straumann	BoneLevel	2.9 (SC)	0700-BLS					Single	-	GMS-BLSL	LL2-BLS-R2	CMFIT10-BLS-R2	-	3D-BLS
701					Bridge	SC-BLM	GMS-BLMS	LL2-BLMR2 LL2-BLM20-R2 LL2-BLM-N2 LL2-BLM20-N2	-	-	-	-	CMFIT10-BLM-R2	CMFIT14-BLM-R2	3D-BLM
702	Straumann	BoneLevel	4.1/4.8 (RC)	0702-BLR											
703					Straumann	Tissue Level	3.5 (NN)	0703-IT35	Single	SC-IT35	-	-	LS2-IT35-O	CMFIT10-IT35-O	-
704	Straumann	Tissue Level	4.8 (RN)	0704-ITR					Single	GMS-ITRS	LL2-ITR-O2 LL2-ITR20-O2	LS2-ITR-O2	CMFIT10-ITR-O2	CMFIT14-ITR-O2	3D-ITR
705					Straumann	BLX	6.5 (WN)	0705-ITW	Single	SC-ITRW	GMS-ITWS	LL2-ITW-O2 LL2-ITW20-O2	LS2-ITW-O2	CMFIT10-ITW-O2	CMFIT14-ITW-O2
710	Straumann	BLX	3.75/4.5 (RB)	0710-BLX					Single	-	GMS-BLXL	LL2-BLX-T2	CMFIT10-BLX-T2	CMFIT14-BLX-T2	3D-BLX
711					Straumann	BLX	5.5/6.5 (WB)	0711-BLXW	Single	-	-	-	CMFIT10-BLXW-T2	CMFIT14-BLXW-T2	-
801	Nobel Biocare	Active	3.0	0801-AC30					Single	OS-AC30	-	LS2-AC30-T2	CMFIT10-AC30-T2	-	3D-AC30
802					Nobel Biocare	Active	3.5 (NP)	0802-ACN	Single	SC-ACN	GMS-ACNS	LL2-ACNH2 LL2-ACN20-H2 LL2-ACNH2 LL2-ACN20-N2	LS2-ACN-H2	CMFIT14-ACN-H2	3D-ACN
803	Nobel Biocare	Active	4.3/5.0 (RP)	0803-ACR					Single	SC-ACR	GMS-ACRS	LL2-ACR-H2 LL2-ACR20-H2 LL2-ACR-N2 LL2-ACR20-N2	LS2-ACR-H2	CMFIT14-ACR-H2	3D-ACR
804					Nobel Biocare	Active	5.5 (WP)	0804-ACW	Single	SC-ACW	GMS-ACWS	LL2-ACW-H2 LL2-ACW-N2	LS2-ACW-H2	CMFIT14-ACW-H2	3D-ACW
805	Nobel Biocare	Active	3.5 (NP)	0805-EXMB					Single	SC-EXMB	-	LS2-EXMB-H2	CMFIT10-EXMB-H2	CMFIT14-EXMB-H2	3D-EXM
806					Nobel Biocare	Active	4.1 (RP)	0806-EXRB	Single	SC-EXRB	-	LS2-EXRB-H2	CMFIT10-EXRB-H2	CMFIT14-EXRB-H2	3D-EXR
807	Nobel Biocare	Active	5.1 (WP)	0807-EXWB					Single	SC-EXWB	-	LS2-EXWB-H2	CMFIT10-EXWB-H2	CMFIT14-EXWB-H2	3D-EXW

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMFik 10mm	CMFik 14mm	3D Analog	Screw (see Pack)	CCM Abutment (Co-Cr base)	Sirona Tibase		
808	Nobel Biocare	Replace Select	3.5 (NP)	0808-RPM	Single	SC-RPM	GMS-RPMS	LL2-RPM-E2	LS2-RPM-E2	CMFIT10-RPM-E2	CMFIT14-RPM-E2	3D-RPM	AC-RPM5	RPM-E2	SRN2-RPM-E2		
Bridge						GMS-RPML	LL2-RPM-N2										
Single					SC-RPR	GMS-RPRS	LL2-RPR-E2	LS2-RPR-E2	CMFIT10-RPR-E2	CMFIT14-RPR-E2	3D-RPR					RPR-E2	SRN2-RPR-E2
Bridge						GMS-RPRL	LL2-RPR-N2										
810	Nobel Biocare	Replace Select	5 (WP)	0810-RPW	Single	SC-RPW	GMS-RPWS	LL2-RPW-E2	LS2-RPW-E2	CMFIT10-RPW-E2	CMFIT14-RPW-E2	3D-RPW	AC-RPR5	RPW-E2	SRN2-RPW-E2		
Bridge						GMS-RPWL	LL2-RPW-N2										
Single					SC-RP60	GMS-RP60S	LL2-RP60-E2	LS2-RP60-E2	CMFIT10-RP60-E2	CMFIT14-RP60-E2	3D-RP60						
Bridge						GMS-RP60L	LL2-RP60-N2										
901			3.0	0901-XV3	Single	SC-XV3	GMS-XV3L	LL2-XV3-H2	LS2-XV3-H2	CMFIT10-XV3-H2		3D-XV3	AC-XV35				
Bridge						GMS-XV3M	LL2-XV3-N2										
Single					SC-XVM	GMS-XVMS	LL2-XVM-H2	LS2-XVM-H2	CMFIT10-XVM-H2	3D-XVM						XVM-H2	SRN2-XVM-H2
Bridge						GMS-XVML	LL2-XVM-N2										
902		Xive	3.4	0902-XVM	Single	SC-XVR	GMS-XVRS	LL2-XVR-H2	LS2-XVR-H2	CMFIT10-XVR-H2	CMFIT14-XVR-H2	3D-XVR	AC-XV5	XVR-H2	SRN2-XVR-H2		
Bridge						GMS-XVRL	LL2-XVR-N2										
Single					SC-XVW	GMS-XVWS	LL2-XVW-H2	LS2-XVW-H2	CMFIT10-XVW-H2	3D-XVW						XVW-H2	SRN2-XVW-H2
Bridge						GMS-XVWL	LL2-XVW-N2										
903	Dentsply		4.5	0904-XV5	Single	SC-XV5	GMS-XV5S	LL2-XV5-H2	LS2-XV5-H2	CMFIT10-XV5-H2	CMFIT14-XV5-H2	3D-XV5					
Bridge						GMS-XV5L	LL2-XV5-N2										
Single					OS-AKLL		LL2-AKLL-T		CMFH10-AKLL-T	3D-AKLL							
Bridge							LL2-AKLL20-T										
906		Ankylos	3.4/4.5/5.5 (Universal)	0906-AKLL	Single			LL2-AKLL-N		CMFH10-AKLL-N							
Bridge							LL2-AKLL20-N										
Single					SC-3ICM	GMS-3ICMS	LL2-3ICM-H2	LS2-3ICM-H2	CMFIT10-3ICM-H2	3D-3ICM							
Bridge						GMS-3ICML	LL2-3ICM-N2										
1001	Biomet 3i	Certain	3.4	1001-3ICM	Single			LL2-3ICR-H2	LS2-3ICR-H2	CMFIT10-3ICR-H2	CMFIT14-3ICR-H2	3D-3ICR	AC-3IC5	3ICM-H2			
Bridge							LL2-3ICR-N2										
Single					SC-3ICR	GMS-3ICRS	LL2-3ICW-H2	LS2-3ICW-H2	CMFIT10-3ICW-H2	3D-3ICW							
Bridge						GMS-3ICRL	LL2-3ICW-N2										
1002			4.1	1002-3ICR	Single			LL2-3IC60-H2	LS2-3IC60-H2	CMFIT10-3IC60-H2	CMFIT14-3IC60-H2	3D-3IC60					
Bridge							LL2-3IC60-N2										
Single					SC-EX3M		LL2-EX3M-H2		CMFIT10-EX3M-H2								
Bridge							LL2-EX3M-N2										
1003		External / Osseotite	5.0	1003-3ICW	Single			LL2-EXR-H2	LS2-EXR-H2	CMFIT10-EXR-H2	CMFIT14-EXR-H2	3D-EXR	AC-EXR5	EXR-H2	SRN2-EXR-H2		
Bridge							LL2-EXR-N2										
Single					SC-EXR	GMS-EXRS	LL2-EXT-H2	LS2-EXT-H2	CMFIT10-EXT-H2								
Bridge						GMS-EXRL	LL2-EXT-N2										
1004		Legacy	6.0	1008-EXTW	Single			LL2-EXTW-H2	LS2-EXTW-H2	CMFIT10-EXTW-H2	CMFIT14-EXTW-H2						
Bridge							LL2-EXTW-N2										
Single					OS-TSVM		LL2-TSVM-H2		CMFIT10-TSVM-H2	3D-TSVM							
Bridge							LL2-TSVM-N2										
1101	Zimmer	Tapered Screw-Vent	3.5	1201-TSVM	Single			LL2-TSVM-H2	LS2-TSVM-H2	CMFIT10-TSVM-H2	CMFIT14-TSVM-H2	3D-TSVM	AC-TSVM5				
Bridge							LL2-TSVM-N2										
Single					SC-TSVM	GMS-TSVM5	LL2-TSVM30-H2	LS2-TSVM30-H2	CMFIT10-TSVM30-H2	3D-TSVM30							
Bridge						GMS-TSVM3L	LL2-TSVM30-N2										
1202		Spine	4.5	1202-TSVM	Single			LL2-TSVM4-H2	LS2-TSVM4-H2	CMFIT10-TSVM4-H2	CMFIT14-TSVM4-H2	3D-TSVM4	AC-TSVM5				
Bridge							LL2-TSVM4-N2										
Single					SC-TSVM4	GMS-TSVM4S	LL2-TSVM420-H2	LS2-TSVM420-H2	CMFIT10-TSVM420-H2	3D-TSVM420							
Bridge						GMS-TSVM4L	LL2-TSVM420-N2										
1203		Spine	5.7	1203-TSVM	Single			LL2-TSVM5-H2	LS2-TSVM5-H2	CMFIT10-TSVM5-H2	CMFIT14-TSVM5-H2	3D-TSVM5					
Bridge							LL2-TSVM5-N2										
Single					SC-SP40		LL2-SP40-H		CMFIT14-SP40-H								
Bridge							LL2-SP40-N										
1204			3.7/4.0	1204-SP50	Single			LL2-SP50-H		CMFIT14-SP50-H							
Bridge							LL2-SP50-N										
1205			5.0	1205-SP50	Single			LL2-SP50-H		CMFIT14-SP50-H							

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMFit 10mm	CMFit 14mm	3D Analog	Screw (Sea Pack)	CCM Abutment (Co-Cr base)	Sirona Ti-base
1301			3.0 (Yellow)	1301-AS3	Single	SC-AS3	GMS-AS3L	LL2-AS3-T2	-	CMFIT10-AS3-T2	-	3D-AS3	AC-AS35	AS3-T2	-
1302		Osseospeed	3.5/4.0 (Aqua)	1302-ASR	Single	SC-ASR	GMS-ASRS	LL2-ASR-H2 LL2-ASR20-H2 LL2-ASR-N2 LL2-ASR20-N2	LS2-ASR-H2	CMFIT10-ASR-H2	CMFIT14-ASR-H2	3D-ASR	AC-ASR5	ASR-H2 ASR-N2	SRN2-ASR-H2
1303	Bridge														
1304			4.5/5.0 (Lilac)	1303-ASW	Single	SC-ASW	GMS-ASWS	LL2-ASW-H2 LL2-ASW20-H2 LL2-ASW-N2 LL2-ASW20-N2	LS2-ASW-H2	CMFIT10-ASW-H2	CMFIT14-ASW-H2	3D-ASW	AC-ASW5	ASW-H2 ASW-N2	SRN2-ASW-H2
1305	Astra Tech		3 (Green)	1304-AEV30	Single	SC-AEV30	GMS-AEV30S GMS-AEV30L	LL2-AEV30-T2 LL2-AEV30-N2	-	CMFIT10-AEV30-T2	-	3D-AEV30	AC-AEV305	-	-
1306		Bridge													
1307		EV	3.6 (Purple)	1305-AEV36	Single	SC-AEV36	GMS-AEV36S GMS-AEV36L	LL2-AEV36-T2 LL2-AEV3620-T2 LL2-AEV36-N2 LL2-AEV3620-N2	-	CMFIT10-AEV36-T2	CMFIT14-AEV36-T2	3D-AEV36	AC-AEV365	-	-
1308	Bridge														
1309			4.2 (Yellow)	1306-AEV42	Single	SC-AEV42	GMS-AEV42S GMS-AEV42L	LL2-AEV42-T2 LL2-AEV4220-T2 LL2-AEV42-N2 LL2-AEV4220-N2	-	CMFIT10-AEV42-T2	CMFIT14-AEV42-T2	3D-AEV42	AC-AEV425	-	-
1310			4.8 (Blue)	1307-AEV48	Single	SC-AEV48	GMS-AEV48S GMS-AEV48L	LL2-AEV48-T2 LL2-AEV4820-T2 LL2-AEV48-N2 LL2-AEV4820-N2	-	CMFIT10-AEV48-T2	CMFIT14-AEV48-T2	3D-AEV48	AC-AEV48545	-	-
1311			5.4 (Brown)	1308-AEV54	Single	SC-AEV54	GMS-AEV54S GMS-AEV54L	LL2-AEV54-T2 LL2-AEV5420-T2 LL2-AEV54-N2 LL2-AEV5420-N2	-	CMFIT10-AEV54-T2	CMFIT14-AEV54-T2	3D-AEV54	AC-AEV545	-	-
1312			3.3	1401-CAM33	Single	SC-CAM33	GMS-CAM33S GMS-CAM33L	LL2-CAM33-T2 LL2-CAM33-N2	-	CMFIT10-CAM33-T2	-	3D-CAM33	AC-CAM335	-	-
1313			3.8	1402-CAM38	Single	SC-CAM38	GMS-CAM38S GMS-CAM38L	LL2-CAM38-T2 LL2-CAM38-N2	-	CMFIT10-CAM38-T2	CMFIT14-CAM38-T2	3D-CAM38	AC-CAM385	-	-
1314		Camlog	4.3	1403-CAM43	Single	SC-CAM43	GMS-CAM43S GMS-CAM43L	LL2-CAM43-T2 LL2-CAM43-N2	-	CMFIT10-CAM43-T2	CMFIT14-CAM43-T2	3D-CAM43	AC-CAM435	-	-
1315	Bridge														
1316			5.0	1404-CAM50	Single	SC-CAM50	GMS-CAM50S GMS-CAM50L	LL2-CAM50-T2 LL2-CAM50-N2	-	CMFIT10-CAM50-T2	CMFIT14-CAM50-T2	3D-CAM50	AC-CAM50605	-	-
1317			6.0	1405-CAM60	Single	SC-CAM60	GMS-CAM60S GMS-CAM60L	LL2-CAM60-T2 LL2-CAM60-N2	-	CMFIT10-CAM60-T2	CMFIT14-CAM60-T2	3D-CAM60	AC-CAM60605	-	-
1318			3.3	1406-CN33	Single	SC-CN33	GMS-CN33S GMS-CN33L	LL2-CN33-T2 LL2-CN33-N2	-	CMFIT10-CN33-E	-	3D-CN33	AC-CN335	-	-
1319		Coneolog	3.8/4.3	1407-CN3843	Single	SC-CN3843	GMS-CN3843S GMS-CN3843L	LL2-CN3843-T2 LL2-CN3843-N2	-	CMFIT10-CN3843-E	CMFIT14-CN3843-E	3D-CN3843	AC-CN38435	-	-
1320	Bridge														
1321			5.0	1408-CN50	Single	SC-CN50	GMS-CN50S GMS-CN50L	LL2-CN50-T2 LL2-CN50-N2	-	CMFIT10-CN50-E	CMFIT14-CN50-E	3D-CN50	AC-CN505	-	-
1322			3.3 (Narrow)	1501-MS33	Single	SC-MS33	GMS-MS33S GMS-MS33L	LL2-MS33-H LL2-MS33-N	-	CMFIT10-MS33-H	-	3D-MS33	AC-MS3335	-	-
1323		Seven	3.75/4.2 (Standard)	1502-MSR	Single	SC-TS33	GMS-TS33S GMS-TS33L	LL2-MSR-H LL2-MSR-N	-	CMFIT10-MSR-H	CMFIT14-MSR-H	3D-TS33	AC-MS335	-	-
1324	Bridge														
1325			5.0/6.0 (Wide)	1503-MSW	Single	SC-TS34	GMS-TS34S GMS-TS34L	LL2-MSW-H LL2-MSW-N	-	CMFIT10-MSW-H	CMFIT14-MSW-H	3D-TS34	AC-MS5	-	-
1326			3.3	1504-C1M	Single	SC-C1M	GMS-C1MS GMS-C1ML	LL2-C1M-T LL2-C1M-N	-	CMFIT10-C1M-T	CMFIT14-C1M-T	3D-C1M	AC-C1M5	C1M-T C1M-N	-
1327		C1	3.75/4.2	1505-C1R	Single	SC-C1R	GMS-C1RS GMS-C1RL	LL2-C1R-T LL2-C1R-N	-	CMFIT10-C1R-T	CMFIT14-C1R-T	3D-C1R	AC-C15	-	-
1328	Bridge														
1329			5.0	1509-C1W	Single	SC-C1W	GMS-C1WS GMS-C1WL	LL2-C1W-T LL2-C1W-N	-	CMFIT10-C1W-T	CMFIT14-C1W-T	3D-C1W	AC-C15	C1W-T C1W-N	-

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMFk 10mm	CMFk 14mm	3D Analog	Screw (See Pack)	CCM Abutment (Co-Cr base)	Sirona TI-base							
1601	SIC	ace / SIC max	3.3	1601-SICM	Single Bridge	SC-SIC		LL2-SICM-H		CMFIT10-SICM-H	CMFIT14-SICM-H	3D-SICM	AC-SIC5	SICM-H								
1602				1602-SICR	Single Bridge			LL2-SICR-H		CMFIT10-SICR-H	CMFIT14-SICR-H	3D-SICR		SICR-H								
1701	Thommen	SPI	3.5	1701-SPI35	Single Bridge	SC-SPI35		LL2-SPI35-H		CMFIT10-SPI35-H	-	3D-SPI35	AC-SPI35	-								
1702				1702-SPI40	Single Bridge			LL2-SPI35-N		CMFIT10-SPI40-H	CMFIT14-SPI40-H	3D-SPI40		-								
1703				1703-SPI45	Single Bridge			LL2-SPI40-H		CMFIT10-SPI45-H	CMFIT14-SPI45-H	3D-SPI45		AC-SPI5		-						
1704				1704-SPI50	Single Bridge			LL2-SPI45-N		CMFIT10-SPI50-H	CMFIT14-SPI50-H	3D-SPI50				-						
1705				1705-SPI60	Single Bridge			LL2-SPI60-N		CMFIT10-SPI60-H	CMFIT14-SPI60-H	-		-								
1801				Bego	Bego			3.7		1801-BEGO37	Single Bridge	SC-BEGO				LL2-BEGO37-H		CMFIT10-BEGO37-H	-	3D-BEGO37	AC-BEGO5	-
1802	1802-BEGO41	Single Bridge	LL2-BEGO37-N			CMFIT10-BEGO41-H	3D-BEGO41		-													
1803	1803-BEGO45	Single Bridge	LL2-BEGO41-N			CMFIT10-BEGO45-H	3D-BEGO45		-													
1804	1804-BEGO55	Single Bridge	LL2-BEGO55-H			CMFIT10-BEGO55-H	3D-BEGO55		-													
1901	Mytis	Arrow	3.7/4.0			1901-AR3740	Single Bridge		OS-AR4760		LL2-BEGO55-H				CMFIT10-AR3740-H	-		AC-AR5	-			
1902						1902-AR4760	Single Bridge				LL2-AR4760-H				CMFIT10-AR4760-H	-			-			
2001	KYOCERA	POI	3.7	2001-POI37	Single Bridge	SC-POI37		LL2-POI37-H		CMFIT10-POI37-H	-	AC-POI375	-	-								
2002				2002-POI42	Single Bridge			LL2-POI42-H		CMFIT10-POI42-H	-	AC-POI425	-									
2101	Biohorizon	Laser Lok	3.0	2101-BH30	Single Bridge	OS-BH30		LL2-BH30-H2		CMFIT10-BH30-H2	-	3D-BH30	AC-BH305	-	-							
2201	Anthorlyr	AXIOM Bone Level	3.4/4.0/4.6/5.2 (Universal)	2201-ANG	Single Bridge	SC-ANG		LL2-ANG-E		CMFIT10-ANG-E	-	3D-ANG	AC-ANG5	ANG-E	-							
2301	WARANTEC	Oneplant	4.3/4.8/5.3 (Standard)	2301-WOP	Single Bridge	SC-WOP		LL2-ANG-N		CMFIT10-WOP-O	-	-	AC-WOP5	-	-							
2302				2302-SUMW	Single Bridge			LL2-SUMW-H2	CMFIT10-SUMW-H2	-	AC-SUMW5	-										
2401	Tekka	In-Kone	3.5/4.0/4.5/5.5 (Standard)	2401-TEK	Single Bridge	SC-TEK		-		CMFIT10-TEK-H	-	3D-TEK	AC-TEK5	TEK-H	-							
2501	Keystone	PrimaComex / Max-Trilobe	3.5	2501-KEP35	Single Bridge	SC-KEP35		LL2-KEP35-T2	LS2-KEP35-T2	CMFIT10-KEP35-T2	CMFIT14-KEP35-T2	3D-KEP35	AC-KEP5	-	-	-						
2502				2502-KEP41	Single Bridge			LL2-KEP41-T2		CMFIT10-KEP41-T2	CMFIT14-KEP41-T2	3D-KEP41		-								
2503				2503-KEP50	Single Bridge			LL2-KEP50-T2		CMFIT10-KEP50-T2	CMFIT14-KEP50-T2	3D-KEP50		-								
2504				2504-KEP50	Single Bridge			LL2-KEP50-N2		CMFIT10-KEP50-N2	CMFIT14-KEP50-N2	3D-KEP50		-								
2505				2505-KEP65	Single Bridge			LL2-KEP65-T2		CMFIT10-KEP65-T2	CMFIT14-KEP65-T2	3D-KEP65		-								
2505				2505-KEP65	Single Bridge			LL2-KEP65-N2		CMFIT10-KEP65-N2	CMFIT14-KEP65-N2	3D-KEP65		-								
2601	Neodent	GM	3.5/4.3/5.0 (Standard)	2601-NCM	Single Bridge	SC-NCM		-		CMFIT10-NCM-H	-	3D-NCM	-	-	-							
2602				2602-NCG	Single Bridge			LL2-NCG-H2	CMFIT10-NCG-H2	3D-NCG	AC-NCG5	-										

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMFit 10mm	CMFit 14mm	3D Analog	Screw (5ea Pack)	CCM Abutment (Co-Cr base)	Sirona Tibase	
2701	BICON	BICON Implant	2.0	2701-BC20	Bridge	SC-BC20	-	LL2-BC20-N	-	CMFIT10-BC20-N	-	-	-	-	-	
2702			2.5	2702-BC25	Bridge	SC-BC25	-	LL2-BC25-N	-	CMFIT10-BC25-N	-	-	-	-	-	-
2703			3.0	2703-BC30	Bridge	SC-BC30	-	LL2-BC30-N	-	CMFIT10-BC30-N	-	-	-	-	-	-
2801	Medentika	IPS	3.5/4.0/4.5/5.0	2801-MTR	Single	SC-MTR	-	-	-	CMFIT10-MTR-R	-	-	AC-MTR5	-	-	
2901	GC	Aadva Implant	4.0/5.0 (Regular /Wide)	2901-GAD	Single	SC-GAD	-	-	-	CMFIT10-GAD-H	-	-	AC-GAD5	-	-	
3001	BIOTECH	KONFAC	3.0	3001-BK30	Single	SC-BK30	-	LL2-BK30-T	-	CMFIT10-BK30-T	-	3D-BK30	AC-BK305	-	-	
3002			3.6/4.2/4.8/5.4	3002-BK40	Single	SC-BK40	-	LL2-BK40-T	-	CMFIT10-BK40-T	-	3D-BK40	AC-BK405	-	-	-
			3.75/4.1/4.8	3101-ICX	Single	SC-ICX	-	LL2-ICX-H	-	CMFIT10-ICX-H	-	3D-ICX	AC-ICX5	-	-	-
3201	B&B Denial	Dura-Vit	4.0/3.5/5.0	3201-BNB	Single	-	LL2-BNB-H	-	CMFIT10-BNB-H	-	3D-BNB	AC-BNB5	-	-		
3301	Bredent Medical	Sky	3.5 / 4.0	3301-SKR	Single	-	LL2-SKR-T	-	CMFIT10-SKR-T	-	3D-SKR	AC-SKR5	-	-		
3401	Necess	ProActive	3.5/4.0/4.5/0.5/5.5	3401-NEP	Single	-	LL2-NEP-T	-	CMFIT10-NEP-T	-	3D-NEP	AC-NEP5	-	-		
3501	Champions	(R)evolution	3.5/4.0/4.5/5.5	3501-CHR	Single	-	GMS-CHRL	-	-	CMFIT10-CHR-H	-	3D-CHR	AC-CHR5	-	-	
3601	KYOCERA	FINESIA Bone Level	3.2 / 3.4	3601-FBLN	Single	-	GMS-FBLNL	-	CMFIT10-FBLNH	-	3D-FBLN	AC-FBLN5	-	-	-	
3602			3.7 / 4.2	3602-FBLR	Single	-	GMS-FBLRL	-	LL2-FBLR-H	-	CMFIT10-FBLRH	-	3D-FBLR	AC-FBLR5	-	-

Geo MultiUnit List (Scanbody / Cylinder / 3D Analog)

NO	Compatible to Implant Systems	Platform	Library Code	Type	SC(OS) Scanbody	LB2-short	LB2-long	Temporary Cylinder	3D Analog	Screw(5ea)	CCM Abutment (Co-Cr base)
203	Dentium Screw Abutment	4.5	0203-MUDE45	Bridge	SC-MUDE	LB2-MUDE45	-	-	-	AC-MUDE5	-
204		5.5	0204-MUDE55			LB2-MUDE55	-	-	-	-	-
706	Straumann TissueLevel SynOcta coping	4.8	0706-MUIT	Bridge	SC-MUIT	LB2-MUIT	-	-	-	AC-MUIT5	-
707		3.4	0707-MUSR35			LB2-MUSR35	-	3D-MUSR35	-	-	AC-MUSR5
708	Straumann Bonelevel SRA	4.6	0708-MUSR46	Bridge	SC-MUSR46	LB2-MUSR46	-	-	3D-MUSR46	-	-
812	Nobel Multi Unit	4.8	0812-MUTB			SC-MUTB	LB2-MUTB2	LB2-MUTB	TA-MUTB	3D-MUTB	AC-MUTB5
907	Dentsply Ankylos Multi-Unit	mini	0907-AKLM	Bridge	SCB-16	LB2-AKLM	-	-	-	AC-AKL5	-
908		regular	0908-AKLR			LB2-AKLR	-	-	-	-	AC-MUTSV5
1206	Zimmer Tapered Screw-Vent Tapered Coping	4.5	1206-MUTSV	Bridge	SC-MUTSV	LB2-MUTSV	-	-	3D-MUTSV	AC-MUTSV5	-
108	Osstem Multi-Abutment	4.8	0108-MUTO			SC2-MUTO	LB2-MUTO	LB2-MUTOL	-	-	AC-MUT05

Geo Transfer Scanbody List

NO	Compatible to Implant Systems	Platform	Type	Section				
				Diameter GH / Height	4.5mm	5.0mm	6.0mm	7.0mm
101	Osstem Implant TS Mini	3.0/3.5	Single	2.0 / 5.0	GTS-SUM4020	GTS-SUM4520		
				4.0 / 7.0	GTS-SUM4040	GTS-SUM4540		
				6.0 / 9.0	GTS-SUM4060	GTS-SUM4560		
102	Osstem Implant TS Standard	4.0/4.5/5.0-	Single	2.0 / 5.0	GTS-SURO4520	GTS-SURO5020	GTS-SURO6020	GTS-SURO7020
				4.0 / 7.0	GTS-SURO4540	GTS-SURO5040	GTS-SURO6040	GTS-SURO7040
				6.0 / 9.0	GTS-SURO4560	GTS-SURO5060	GTS-SURO6060	GTS-SURO7060
201	Dentium Super-Line Implantium	3.7/4.0/4.5-	Single	2.0 / 5.0	GTS-DER4020	GTS-DER4520	GTS-DER5520	GTS-DER6520
				4.0 / 7.0	GTS-DER4040	GTS-DER4540	GTS-DER5540	GTS-DER6540
				6.0 / 9.0	GTS-DER4060	GTS-DER4560	GTS-DER5560	GTS-DER6560

Guide Tool	Compatible to Implant Systems	Diameter						
		4.0mm	4.5mm	5.0mm	5.5mm	6.0mm	6.5mm	7.0mm
Guide Tool	Osstem Implant TS Mini	GTST-S40-120H	GTST-S45-120H	-	-	-	-	-
	Osstem Implant TS Standard	-	GTST-S45-120H	GTST-S50-120H	-	GTST-S60-120H	-	GTST-S70-120H
	Dentium Super-Line Implantium	GTST-D40-127H	GTST-D45-127H	-	GTST-D55-127H	-	GTST-D65-127H	-

GeOMEDI



Digital dentistry has revolutionized the dental industry forever. To provide a better CAD-CAM experience to our customers, we founded GeoMedi in 2009.

GeoMedi's core value is diversity, expertise, and customer satisfaction.

We have one of the most diverse compatible implant systems which meet most of your requirements. We strive to develop more implant systems with superior quality.

(주) 지오메디

- 경기도 의왕시 이미로 40 (포일동, 인덕원IT밸리) D동 B112, 702~705, 709-1, 810-2, 811호
- Tel : 02-588-2871
- Fax : 02-6280-2849
- E-mail : sales@geomedi.com
- WebSite : www.geomedi.com

GeoMedi Co., Ltd.

- #B112, 702~705, 709-1, 810-2, 811, Unit D, Indeokwon IT Valley, 40, Imi-ro, Uiwang-si, Gyeonggi-do, Korea
- Tel : 82-2-588-2871
- Fax : 82-2-6280-2849
- E-mail : sales@geomedi.com
- WebSite : www.geomedi.com



GeOMEDI

PREMIUM CAD/CAM PRODUCTS